Package 'affxparser'

April 9, 2025

```
Version 1.79.1
```

Depends R (>= 2.14.0)

Suggests R.oo (>= 1.22.0), R.utils (>= 2.7.0), AffymetrixDataTestFiles

Title Affymetrix File Parsing SDK

Author Henrik Bengtsson [aut], James Bullard [aut], Robert Gentleman [ctb], Kasper Daniel Hansen [aut, cre], Jim Hester [ctb], Martin Morgan [ctb]

Maintainer Kasper Daniel Hansen <kasperdanielhansen@gmail.com>

Description Package for parsing Affymetrix files (CDF, CEL, CHP, BPMAP, BAR). It provides methods for fast and memory efficient parsing of Affymetrix files using the Affymetrix' Fusion SDK. Both ASCII- and binary-based files are supported. Currently, there are methods for reading chip definition file (CDF) and a cell intensity file (CEL). These files can be read either in full or in part. For example, probe signals from a few probesets can be extracted very quickly from a set of CEL files into a convenient list structure.

Note Fusion SDK v1.1.2

License LGPL (>= 2)

LazyLoad yes

URL https://github.com/HenrikBengtsson/affxparser

BugReports https://github.com/HenrikBengtsson/affxparser/issues

biocViews Infrastructure, DataImport, Microarray, ProprietaryPlatforms, OneChannel

git_url https://git.bioconductor.org/packages/affxparser

git_branch devel

git_last_commit d2f1dc2

git_last_commit_date 2025-01-06

Repository Bioconductor 3.21

Date/Publication 2025-04-09

2 Contents

Contents

affxparser-package	
1. Dictionary	
2. Cell coordinates and cell indices	6
9. Advanced - Cell-index maps for reading and writing	8
applyCdfGroupFields	10
applyCdfGroups	11
arrangeCelFilesByChipType	14
cdfAddBaseMmCounts	15
cdfAddPlasqTypes	16
cdfAddProbeOffsets	17
cdfGetFields	18
cdfGetGroups	19
cdfGtypeCelToPQ	19
edfHeaderToCelHeader	
edfMergeAlleles	
edfMergeStrands	
edfMergeToQuartets	
edfOrderBy	
edfOrderColumnsBy	
edfSetDimension	
compareCdfs	
compareCels	
convertCdf	
convertCel	
copyCel	
createCel	
findCdf	
findFiles	
nvertMap	
sCelFile	
parseDatHeaderString	
readBpmap	
readCcg	
readCcgHeader	
readCdf	
readCdfCellIndices	
readCdfDataFrame	
readCdfGroupNames	
readCdfHeader	
readCdfIsPm	
readCdfNbrOfCellsPerUnitGroup	
readCdfQc	
readCdfUnitNames	
readCdfUnits	
readCdfUnitsWriteMap	
eadCatOmiswitteMap	50

affxparser-package 3

	parser-package Package affxparser	
Index		86
	writeTpmap	84
	writeCelHeader	83
	writeCdfUnits	
	writeCdfQcUnits	
	writeCdfHeader	
	writeCdf	
	updateCelUnits	
	updateCel	
	readPgfHeader	
	readPgfEnv	
	readPgf	
	readClfHeader	
	readClfEnv	
	readChp	
	readCelUnits	
	readCelRectangle	
	readCelIntensities	
	readCelHeader	

Description

The **affxparser** package provides methods for fast and memory efficient parsing of Affymetrix files [1] using the Affymetrix' Fusion SDK [2,3]. Both traditional ASCII- and binary (XDA)-based files are supported, as well as Affymetrix future binary format "Calvin". The efficiency of the parsing is dependent on whether a specific file is binary or ASCII.

Currently, there are methods for reading chip definition file (CDF) and a cell intensity file (CEL). These files can be read either in full or in part. For example, probe signals from a few probesets can be extracted very quickly from a set of CEL files into a convenient list structure.

To get started

To get started, see:

- 1. readCelUnits() reads one or several Affymetrix CEL file probeset by probeset.
- 2. readCel() reads an Affymetrix CEL file. by probe.
- 3. readCdf() reads an Affymetrix CDF file. by probe.
- 4. readCdfUnits() reads an Affymetrix CDF file unit by unit.
- 5. readCdfCellIndices() Like readCdfUnits(), but returns cell indices only, which is often enough to read CEL files unit by unit.
- 6. applyCdfGroups() Re-arranges a CDF structure.
- 7. findCdf() Locates an Affymetrix CDF file by chip type. This page also describes how to setup default search path for CDF files.

4 affxparser-package

Setting up the CDF search path

Some of the functions in this package search for CDF files automatically by scanning certain directories. To add directories to the default search path, see instructions in findCdf().

Future Work

Other Affymetrix files can be parsed using the Fusion SDK. Given sufficient interest we will implement this, e.g. DAT files (image files).

Running examples

In order to run the examples, data files must exists in the current directory. Otherwise, the example scripts will do nothing. Most of the examples requires a CDF file or a CEL file, or both. Make sure the CDF file is of the same chip type as the CEL file.

Affymetrix provides data sets of different types at http://www.affymetrix.com/support/datasets. affx that can be used. There are both small are very large data sets available.

Technical details

This package implements an interface to the Fusion SDK from Affymetrix.com. This SDK (software development kit) is an open source library used for parsing the various files formats used by the Affymetrix platform.

The intention is to provide interfaces to most if not all file formats which may be parsed using Fusion.

The SDK supports parsing of all the different versions of a specific file format. This means that ASCII, binary as well as the new binary format (codename Calvin) used by Affymetrix is supported through a single API. We also expect any future changes to the file formats to be reflected in the SDK, and subsequently in this package.

However, as the current Fusion SDK does not support compressed files, neither does **affxparser**. This is in contrast to some of the existing code in **affy** and relatives (see below for links).

In general we aim to provide functions returning all information in the respective files. Currently it seems that future Affymetrix chip designs may consists of so many features that returning all information will lead to an unnecessary overhead in the case a user only wants access to a subset. We have tried to make this possible.

For older file, certain entries in the files have been removed from newer specifications, and the SDK does not provide utilities for reading these entries. This includes for instance the FEAT column of CDF files.

Currently the package as well as the Fusion SDK is in beta stage. Bugs may be related to either codebase. We are very interested in users being unable to compile/parse files using this library - this includes users with custom chip designs.

In addition, since we aim to return all information stored in the file (and accessible using the Fusion SDK) we would like reports from users being unable to do that.

The efficiency of the underlying code may vary with the version of the file being parsed. For example, we currently report the number of outliers present in a CEL file when reading the header of the file using readCelHeader. In order to obtain this information from text based CEL files

1. Dictionary 5

(version 2), the entire file needs to be read into memory. With version 3 of the file format, this information is stored in the header.

With the introduction of the Fusion SDK (and the next version of their file formats) Affymetrix has made it possible to use multibyte character sets. This implies that character information may be inaccessible if the compiler used to compile the C++ code does not support multibyte character sets (specifically we require that the R installation has defined the macro SUPPORT_MCBS in the Rconfig.h header file). For example GCC needs to be version 3.4 or greater on Solaris.

In the info subdirectory of the package installation, information regarding changes to the Fusion SDK is stored, e.g.

```
pathname <- system.file("info", "changes2fusion.txt", package="affxparser")
file.show(pathname)</pre>
```

Acknowledgments

We would like to thanks Ken Simpson (WEHI, Melbourne) and Seth Falcon (FHCRC, Seattle) for feedback and code contributions.

License

The releases of this package is licensed under LGPL version 2.1 or newer. This applies also to the Fusion SDK.

Author(s)

Henrik Bengtsson [aut], James Bullard [aut], Robert Gentleman [ctb], Kasper Daniel Hansen [aut, cre], Martin Morgan [ctb]

References

- [1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, April, 2006. http://www.affymetrix.com/support/developer/
- [2] Affymetrix Inc, Fusion Software Developers Kit (SDK), 2006. http://www.affymetrix.com/support/developer/fusion/
- [3] Henrik Bengtsson, unofficial archive of Affymetrix Fusion Software Developers Kit (SDK), https://github.com/HenrikBengtsson/Affx-Fusion-SDK
- 1. Dictionary

Description

This part describes non-obvious terms used in this package.

affxparser The name of this package.

API Application program interface, which describes the functional interface of underlying methods.

block (aka group).

BPMAP A file format containing information related to the design of the tiling arrays.

Calvin A special binary file format.

CDF A file format: chip definition file.

CEL A file format: cell intensity file.

cell (aka feature) A probe.

cell index An integer that identifies a probe uniquely.

chip An array.

chip type An identifier specifying a chip design uniquely, e.g. "Mapping50K_Xba240".

DAT A file format: contains pixel intensity values collected from an Affymetrix GeneArray scanner.

feature A probe.

Fusion SDK Open-source software development kit (SDK) provided by Affymetrix to access their data files.

group (aka block) Defines a unique subset of the cells in a unit. Expression arrays typically only have one group per unit, whereas SNP arrays have either two or four groups per unit, one for each of the two allele times possibly repeated for both strands.

MM Mismatch-match, e.g. MM probe.

PGF A file format: probe group file.

TPMAP A file format storing the relationship between (PM,MM) pairs (or PM probes) and positions on a set of sequences.

QC Quality control, e.g. QC probes and QC probe sets.

unit A probeset.

XDA A file format, aka as the binary file format.

2. Cell coordinates and cell indices

2. Cell coordinates and cell indices

Description

This part describes how Affymetrix cells, also known as probes or features, are addressed.

Cell coordinates

In Affymetrix data files, cells are uniquely identified by there *cell coordinates*, i.e. (x,y). For an array with N * K cells in N rows and K columns, the x coordinate is an integer in [0,K-1], and the y coordinate is an integer in [0,N-1]. The cell in the upper-left corner has coordinate (x,y)=(0,0) and the one in the lower-right corner (x,y)=(K-1,N-1).

Cell indices and cell-index offsets

To simplify addressing of cells, a coordinate-to-index function is used so that each cell can be addressed using a single integer instead (of two). Affymetrix defines the *cell index*, i, of cell (x, y) as

$$i = K * y + x + 1,$$

where one is added to give indices in [1, N*K]. Continuing, the above definition means that cells are ordered row by row, that is from left to right and from top to bottom, starting at the upper-left corner. For example, with a chip layout (N, K) = (1600, 1600) the cell at (x, y) = (0, 0) has index i=1, and the cell at (x, y) = (1599, 1599) has index i = 2560000. A cell at (x, y) = (1498, 3) has index i = 6299.

Given the cell index i, the coordinate (x, y) can be calculated as

$$y = floor((i-1)/K)$$

$$x = (i-1) - K * y.$$

Continuing the above example, the coordinate for cell i = 1 is be found to be (x, y) = (0, 0), for cell i = 2560000 it is (x, y) = (1599, 1599), for cell i = 6299 is it (x, y) = (1498, 3).

Converting between cell indices and (x,y) coordinates in R

Although not needed to use the methods in this package, to get the cell indices for the cell coordinates or vice versa, see xy2indices() and indices2xy() in the affy package.

Note on the zero-based "index" field of Affymetrix CDF files

An Affymetrix CDF file provides information on which cells should be grouped together. To identify these groups of cells, the cells are specified by their (x,y) coordinates, which are stored as zero-based coordinates in the CDF file.

All methods of the **affxparser** package make use of these (x,y) coordinates, and some methods make it possible to read them as well. However, it is much more common that the methods return cell indices *calculated* from the (x,y) coordinates as explained above.

In order to conveniently work with cell indices in R, the convention in *affxparser* is to use *one-based* indices. Hence the addition (and subtraction) of 1:s in the above equations. This is all taken care of by **affxparser**.

Note that, in addition to (x,y) coordinates, a CDF file also contains a one-based "index" for each cell. This "index" is redundant to the (x,y) coordinate and can be calculated analogously to the above *cell index* while leaving out the addition (subtraction) of 1:s. Importantly, since this "index" is redundant (and exists only in CDF files), we have decided to treat this field as an internal field. Methods of **affxparser** do neither provide access to nor make use of this internal field.

Author(s)

Henrik Bengtsson

9. Advanced - Cell-index maps for reading and writing

9. Advanced - Cell-index maps for reading and writing

Description

This part defines read and write maps that can be used to remap cell indices before reading and writing data from and to file, respectively.

This package provides methods to create read and write (cell-index) maps from Affymetrix CDF files. These can be used to store the cell data in an optimal order so that when data is read it is read in contiguous blocks, which is faster.

In addition to this, read maps may also be used to read CEL files that have been "reshuffled" by other software. For instance, the dChip software (http://www.dchip.org/) rotates Affymetrix Exon, Tiling and Mapping 500K data. See example below how to read such data "unrotated".

For more details how cell indices are defined, see 2. Cell coordinates and cell indices.

Motivation

When reading data from file, it is faster to read the data in the order that it is stored compared with, say, in a random order. The main reason for this is that the read arm of the hard drive has to move more if data is not read consecutively. Same applies when writing data to file. The read and write cache of the file system may compensate a bit for this, but not completely.

In Affymetrix CEL files, cell data is stored in order of cell indices. Moreover, (except for a few early chip types) Affymetrix randomizes the locations of the cells such that cells in the same unit (probeset) are scattered across the array. Thus, when reading CEL data arranged by units using for instance readCelUnits(), the order of the cells requested is both random and scattered.

Since CEL data is often queried unit by unit (except for some probe-level normalization methods), one can improve the speed of reading data by saving data such that cells in the same unit are stored together. A *write map* is used to remap cell indices to file indices. When later reading that data back, a *read map* is used to remap file indices to cell indices. Read and write maps are described next.

Definition of read and write maps

Consider cell indices i = 1, 2, ..., N * K and file indices j = 1, 2, ..., N * K. A read map is then a bijective (one-to-one) function h() such that

$$i = h(j),$$

and the corresponding write map is the inverse function $h^{-1}()$ such that

$$j = h^{-1}(i).$$

Since the mapping is required to be bijective, it holds that $i=h(h^{-1}(i))$ and that $j=h^{-1}(h(j))$. For example, consider the "reversing" read map function h(j)=N*K-j+1. The write map function is $h^{-1}(i)=N*K-i+1$. To verify the bijective property of this map, we see that $h(h^{-1}(i))=h(N*K-i+1)=N*K-(N*K-i+1)+1=i$ as well as $h^{-1}(h(j))=h^{-1}(N*K-j+1)=N*K-(N*K-j+1)+1=j$.

Read and write maps in R

In this package, read and write maps are represented as integer vectors of length N * K with *unique* elements in $\{1, 2, ..., N * K\}$. Consider cell and file indices as in previous section.

For example, the "reversing" read map in previous section can be represented as

```
readMap <- (N*K):1
```

Given a vector j of file indices, the cell indices are the obtained as i = readMap[j]. The corresponding write map is

```
writeMap <- (N*K):1
```

and given a vector i of cell indices, the file indices are the obtained as j = writeMap[i].

Note also that the bijective property holds for this mapping, that is i == readMap[writeMap[i]] and i == writeMap[readMap[i]] are both TRUE.

Because the mapping is bijective, the write map can be calculated from the read map by:

```
writeMap <- order(readMap)</pre>
```

and vice versa:

```
readMap <- order(writeMap)</pre>
```

Note, the invertMap() method is much faster than order().

Since most algorithms for Affymetrix data are based on probeset (unit) models, it is natural to read data unit by unit. Thus, to optimize the speed, cells should be stored in contiguous blocks of units. The methods readCdfUnitsWriteMap() can be used to generate a *write map* from a CDF file such that if the units are read in order, readCelUnits() will read the cells data in order. Example:

```
Find any CDF file
cdfFile <- findCdf()

# Get the order of cell indices
indices <- readCdfCellIndices(cdfFile)
indices <- unlist(indices, use.names=FALSE)

# Get an optimal write map for the CDF file</pre>
```

```
writeMap <- readCdfUnitsWriteMap(cdfFile)

# Get the read map
readMap <- invertMap(writeMap)

# Validate correctness
indices2 <- readMap[indices] # == 1, 2, 3, ..., N*K</pre>
```

Warning, do not misunderstand this example. It can not be used improve the reading speed of default CEL files. For this, the data in the CEL files has to be rearranged (by the corresponding write map).

Reading rotated CEL files

It might be that a CEL file was rotated by another software, e.g. the dChip software rotates Affymetrix Exon, Tiling and Mapping 500K arrays 90 degrees clockwise, which remains rotated when exported as CEL files. To read such data in a non-rotated way, a read map can be used to "unrotate" the data. The 90-degree clockwise rotation that dChip effectively uses to store such data is explained by:

```
h <- readCdfHeader(cdfFile)
# (x,y) chip layout rotated 90 degrees clockwise
nrow <- h$cols
ncol <- h$rows
y <- (nrow-1):0
x <- rep(1:ncol, each=nrow)
writeMap <- as.vector(y*ncol + x)</pre>
```

Thus, to read this data "unrotated", use the following read map:

```
readMap <- invertMap(writeMap)
data <- readCel(celFile, indices=1:10, readMap=readMap)</pre>
```

Author(s)

Henrik Bengtsson

applyCdfGroupFields

Applies a function to a list of fields of each group in a CDF structure

Description

Applies a function to a list of fields of each group in a CDF structure.

applyCdfGroups 11

Usage

```
applyCdfGroupFields(cdf, fcn, ...)
```

Arguments

cdf A CDF list structure.

fcn A function that takes a list structure of fields and returns an updated list of

fields.

... Arguments passed to the fcn function.

Value

Returns an updated CDF list structure.

Author(s)

Henrik Bengtsson

See Also

```
applyCdfGroups().
```

applyCdfGroups

Applies a function over the groups in a CDF structure

Description

Applies a function over the groups in a CDF structure.

Usage

```
applyCdfGroups(cdf, fcn, ...)
```

Arguments

cdf A CDF list structure.

fcn A function that takes a list structure of group elements and returns an updated

list of groups.

... Arguments passed to the fcn function.

Value

Returns an updated CDF list structure.

12 applyCdfGroups

Pre-defined restructuring functions

Generic: • cdfGetFields() - Gets a subset of groups fields in a CDF structure.

- cdfGetGroups() Gets a subset of groups in a CDF structure.
- cdfOrderBy() Orders the fields according to the value of another field in the same CDF group.
- cdfOrderColumnsBy() Orders the columns of fields according to the values in a certain row of another field in the same CDF group.

Designed for SNP arrays: • cdfAddBaseMmCounts() - Adds the number of allele A and allele B mismatching nucleotides of the probes in a CDF structure.

- cdfAddProbeOffsets() Adds probe offsets to the groups in a CDF structure.
- cdfGtypeCelToPQ() Function to imitate Affymetrix' gtype_cel_to_pq software.
- cdfMergeAlleles() Function to join CDF allele A and allele B groups strand by strand.
- cdfMergeStrands() Function to join CDF groups with the same names.

We appreciate contributions.

Author(s)

Henrik Bengtsson

Examples

```
if (require("AffymetrixDataTestFiles")) {
cdfFile <- findCdf("Mapping10K_Xba131")</pre>
# Identify the unit index from the unit name
unitName <- "SNP_A-1509436"
unit <- which(readCdfUnitNames(cdfFile) == unitName)</pre>
# Read the CDF file
cdf0 <- readCdfUnits(cdfFile, units=unit, stratifyBy="pmmm", readType=FALSE, readDirection=FALSE)</pre>
cat("Default CDF structure:\n")
print(cdf0)
# Tabulate the information in each group
cdf <- readCdfUnits(cdfFile, units=unit)</pre>
cdf <- applyCdfGroups(cdf, lapply, as.data.frame)</pre>
print(cdf)
# Infer the (true or the relative) offset for probe quartets.
cdf <- applyCdfGroups(cdf0, cdfAddProbeOffsets)</pre>
cat("Probe offsets:\n")
print(cdf)
```

applyCdfGroups 13

```
# Identify the number of nucleotides that mismatch the
# allele A and the allele B sequences, respectively.
cdf <- applyCdfGroups(cdf, cdfAddBaseMmCounts)</pre>
cat("Allele A & B target sequence mismatch counts:\n")
print(cdf)
# Combine the signals from the sense and the anti-sense
# strands in a SNP CEL files.
# First, join the strands in the CDF structure.
cdf <- applyCdfGroups(cdf, cdfMergeStrands)</pre>
cat("Joined CDF structure:\n")
print(cdf)
# Rearrange values of group fields into quartets. This
# requires that the values are already arranged as PMs and MMs.
cdf <- applyCdfGroups(cdf0, cdfMergeAlleles)</pre>
cat("Probe quartets:\n")
print(cdf)
# Get the x and y cell locations (note, zero-based)
x <- unlist(applyCdfGroups(cdf, cdfGetFields, "x"), use.names=FALSE)</pre>
y <- unlist(applyCdfGroups(cdf, cdfGetFields, "y"), use.names=FALSE)</pre>
# Validate
ncol <- readCdfHeader(cdfFile)$cols</pre>
cells <- as.integer(y*ncol+x+1)</pre>
cells <- sort(cells)</pre>
cells0 <- readCdfCellIndices(cdfFile, units=unit)</pre>
cells0 <- unlist(cells0, use.names=FALSE)</pre>
cells0 <- sort(cells0)</pre>
stopifnot(identical(cells0, cells))
```

arrangeCelFilesByChipType

Moves CEL files to subdirectories with names corresponding to the chip types

Description

Moves CEL files to subdirectories with names corresponding to the chip types according to the CEL file headers. For instance, a HG_U95Av2 CEL file with pathname "data/foo.CEL" will be moved to subdirectory celFiles/HG_U95Av2/.

Usage

```
arrangeCelFilesByChipType(pathnames=list.files(pattern = "[.](cel|CEL)$"),
  path="celFiles/", aliases=NULL, ...)
```

Arguments

pathnames	A character vector of CEL pathnames to be moved.
path	A character string specifying the root output directory, which in turn will contain chip-type subdirectories. All directories will be created, if missing.
aliases	A named character string with chip type aliases. For instance, aliases=c("Focus"="HG-Focus") will treat a CEL file with chiptype label 'Focus' (early-access name) as if it was 'HG-Focus' (official name).
	Not used.

Value

Returns (invisibly) a named character vector of the new pathnames with the chip types as the names. Files that could not be moved or where not valid CEL files are set to missing values.

Author(s)

Henrik Bengtsson

See Also

The chip type is inferred from the CEL file header, cf. readCelHeader().

cdfAddBaseMmCounts 15

cdfAddBaseMmCounts	Adds the number of allele A and allele B mismatching nucleotides of
	the probes in a CDF structure

Description

Adds the number of allele A and allele B mismatching nucleotides of the probes in a CDF structure.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Identifies the number of nucleotides (bases) in probe sequences that mismatch the the target sequence for allele A and the allele B, as used by [1].

Usage

```
cdfAddBaseMmCounts(groups, ...)
```

Arguments

groups A list structure with groups. Each group must contain the fields thase, phase,

and offset (from cdfAddProbeOffsets()).

... Not used.

Details

Note that the above counts can be inferred from the CDF structure alone, i.e. no sequence information is required. Consider a probe group interrogating allele A. First, all PM probes matches the allele A target sequence perfectly regardless of shift. Moreover, all these PM probes mismatch the allele B target sequence at exactly one position. Second, all MM probes mismatches the allele A sequence at exactly one position. This is also true for the allele B sequence, *except* for an MM probe with zero offset, which only mismatch at one (the middle) position. For a probe group interrogating allele B, the same rules applies with labels A and B swapped. In summary, the mismatch counts for PM probes can take values 0 and 1, and for MM probes they can take values 0, 1, and 2.

Value

Returns a list structure with the same number of groups as the groups argument. To each group, two fields is added:

mmACount The number of nucleotides in the probe sequence that mismatches the target

sequence of allele A.

mmBCount The number of nucleotides in the probe sequence that mismatches the target

sequence of allele B.

Author(s)

Henrik Bengtsson

16 cdfAddPlasqTypes

References

[1] LaFramboise T, Weir BA, Zhao X, Beroukhim R, Li C, Harrington D, Sellers WR, and Meyerson M. *Allele-specific amplification in cancer revealed by SNP array analysis*, PLoS Computational Biology, Nov 2005, Volume 1, Issue 6, e65.

[2] Affymetrix, *Understanding Genotyping Probe Set Structure*, 2005. http://www.affymetrix.com/support/developer/whitepapers/genotyping_probe_set_structure.affx

See Also

To add required probe offsets, cdfAddProbeOffsets(). applyCdfGroups().

cdfAddPlasqTypes

Adds the PLASQ types for the probes in a CDF structure

Description

Adds the PLASQ types for the probes in a CDF structure.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Usage

```
cdfAddPlasqTypes(groups, ...)
```

Arguments

groups A list structure with groups. Each group must contain the fields thase, phase,

and expos.

... Not used.

Details

This function identifies the number of nucleotides (bases) in probe sequences that mismatch the the target sequence for allele A and the allele B, as used by PLASQ [1], and adds an integer [0,15] interpreted as one of 16 probe types. In PLASQ these probe types are referred to as: 0=MMoBR, 1=MMoBF, 2=MMcBR, 3=MMcBF, 4=MMoAR, 5=MMoAF, 6=MMcAR, 7=MM-cAF, 8=PMoBR, 9=PMoBF, 10=PMcBR, 11=PMcBF, 12=PMoAR, 13=PMoAF, 14=PMcAR, 15=PM-cAF.

Pseudo rule for finding out the probe-type value:

- PM/MM: For MMs add 0, for PMs add 8.
- A/B: For Bs add 0, for As add 4.
- o/c: For shifted (o) add 0, for centered (c) add 2.

cdfAddProbeOffsets 17

• R/F: For antisense (R) add 0, for sense (F) add 1.

```
Example: (PM,A,c,R) = 8 + 4 + 2 + 0 = 14 (=PMcAR)
```

Value

Returns a list structure with the same number of groups as the groups argument. To each group, one fields is added:

```
plasqType A vector of integers in [0,15].
```

Author(s)

Henrik Bengtsson

References

[1] LaFramboise T, Weir BA, Zhao X, Beroukhim R, Li C, Harrington D, Sellers WR, and Meyerson M. *Allele-specific amplification in cancer revealed by SNP array analysis*, PLoS Computational Biology, Nov 2005, Volume 1, Issue 6, e65.

cdfAddProbeOffsets

Adds probe offsets to the groups in a CDF structure

Description

Adds probe offsets to the groups in a CDF structure.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Usage

```
cdfAddProbeOffsets(groups, ...)
```

Arguments

groups A list structure with groups. Each group must contain the fields tbase, and

expos.

... Not used.

Value

Returns a list structure with half the number of groups as the groups argument (since allele A and allele B groups have been joined).

Author(s)

Henrik Bengtsson

18 cdfGetFields

References

[1] Affymetrix, $Understanding\ Genotyping\ Probe\ Set\ Structure$, 2005. http://www.affymetrix.com/support/developer/whitepapers/genotyping_probe_set_structure.affx

See Also

```
applyCdfGroups().
```

cdfGetFields

Gets a subset of groups fields in a CDF structure

Description

Gets a subset of groups fields in a CDF structure.

This function is designed to be used with applyCdfGroups().

Usage

```
cdfGetFields(groups, fields, ...)
```

Arguments

groups A list of groups.

fields A character vector of names of fields to be returned.

... Not used.

Details

Note that an error is *not* generated for missing fields. Instead the field is returned with value NA. The reason for this is that it is much faster.

Value

Returns a list structure of groups.

Author(s)

Henrik Bengtsson

See Also

applyCdfGroups().

cdfGetGroups 19

cdfGetGroups

Gets a subset of groups in a CDF structure

Description

Gets a subset of groups in a CDF structure.

This function is designed to be used with applyCdfGroups().

Usage

```
cdfGetGroups(groups, which, ...)
```

Arguments

groups A list of groups.

which An integer or character vector of groups be returned.

... Not used.

Value

Returns a list structure of groups.

Author(s)

Henrik Bengtsson

See Also

```
applyCdfGroups().
```

cdfGtypeCelToPQ

Function to imitate Affymetrix' gtype_cel_to_pq software

Description

Function to imitate Affymetrix' gtype_cel_to_pq software.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Usage

```
cdfGtypeCelToPQ(groups, ...)
```

20 cdfHeaderToCelHeader

Arguments

groups A list structure with groups.

... Not used.

Value

Returns a list structure with a single group. The fields in this groups are in turn vectors (all of equal length) where the elements are stored as subsequent quartets (PMA, MMA, PMB, MMB) with all forward-strand quartets first followed by all reverse-strand quartets.

Author(s)

Henrik Bengtsson

References

[1] Affymetrix, *Understanding Genotyping Probe Set Structure*, 2005. http://www.affymetrix.com/support/developer/whitepapers/genotyping_probe_set_structure.affx

See Also

```
applyCdfGroups().
```

cdfHeaderToCelHeader

Creates a valid CEL header from a CDF header

Description

Creates a valid CEL header from a CDF header.

Usage

```
cdfHeaderToCelHeader(cdfHeader, sampleName="noname", date=Sys.time(), ..., version="4")
```

Arguments

cdfHeader A CDF list structure.

sampleName The name of the sample to be added to the CEL header.

date The (scan) date to be added to the CEL header.

... Not used

version The file-format version of the generated CEL file. Currently only version 4 is

supported.

Value

Returns a CDF list structure.

cdfMergeAlleles 21

Author(s)

Henrik Bengtsson

cdfMergeAlleles

Function to join CDF allele A and allele B groups strand by strand

Description

Function to join CDF allele A and allele B groups strand by strand.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Usage

```
cdfMergeAlleles(groups, compReverseBases=FALSE, collapse="", ...)
```

Arguments

groups A list structure with groups.

compReverseBases

If $\ensuremath{\mathsf{TRUE}}$, the group names, which typically are names for bases, are turned into

their complementary bases for the reverse strand.

collapse The character string used to collapse the allele A and the allele B group names.

... Not used.

Details

Allele A and allele B are merged into a matrix where first row hold the elements for allele A and the second elements for allele B.

Value

Returns a list structure with the two groups forward and reverse, if the latter exists.

Author(s)

Henrik Bengtsson

References

[1] Affymetrix, *Understanding Genotyping Probe Set Structure*, 2005. http://www.affymetrix.com/support/developer/whitepapers/genotyping_probe_set_structure.affx

```
applyCdfGroups().
```

22 cdfMergeStrands

cdfMergeStrands

Function to join CDF groups with the same names

Description

Function to join CDF groups with the same names.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

This can be used to join the sense and anti-sense groups of the same allele in SNP arrays.

Usage

```
cdfMergeStrands(groups, ...)
```

Arguments

```
groups A list structure with groups.
... Not used.
```

Details

If a unit has two strands, they are merged such that the elements for the second strand are concatenated to the end of the elements of first strand (This is done separately for the two alleles).

Value

Returns a list structure with only two groups.

Author(s)

Henrik Bengtsson

References

```
[1] Affymetrix, Understanding Genotyping Probe Set Structure, 2005. http://www.affymetrix.com/support/developer/whitepapers/genotyping_probe_set_structure.affx
```

```
applyCdfGroups().
```

cdfMergeToQuartets 23

 ${\tt cdfMergeToQuartets}$

Function to re-arrange CDF groups values in quartets

Description

Function to re-arrange CDF groups values in quartets.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Note, this requires that the group values have already been arranged in PMs and MMs.

Usage

```
cdfMergeToQuartets(groups, ...)
```

Arguments

groups A list structure with groups.

... Not used.

Value

Returns a list structure with the two groups forward and reverse, if the latter exists.

Author(s)

Henrik Bengtsson

References

[1] Affymetrix, *Understanding Genotyping Probe Set Structure*, 2005. http://www.affymetrix.com/support/developer/whitepapers/genotyping_probe_set_structure.affx

```
applyCdfGroups().
```

24 cdfOrderColumnsBy

cdfOrderBy Orders the fields according to the value of another field in the same CDF group
--

Description

Orders the fields according to the value of another field in the same CDF group.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

Usage

```
cdfOrderBy(groups, field, ...)
```

Arguments

groups A list of groups.

field The field whose values are used to order the other fields.

... Optional arguments passed order().

Value

Returns a list structure of groups.

Author(s)

Henrik Bengtsson

See Also

```
cdfOrderColumnsBy(). applyCdfGroups().
```

cdf0rderColumnsBy

Orders the columns of fields according to the values in a certain row of another field in the same CDF group

Description

Orders the columns of fields according to the values in a certain row of another field in the same CDF group. Note that this method requires that the group fields are matrices.

This function is design to be used with applyCdfGroups() on an Affymetrix Mapping (SNP) CDF list structure.

cdfSetDimension 25

Usage

```
cdfOrderColumnsBy(groups, field, row=1, ...)
```

Arguments

groups A list of groups.

field The field whose values in row row are used to order the other fields.

row The row of the above field to be used to find the order.

... Optional arguments passed order().

Value

Returns a list structure of groups.

Author(s)

Henrik Bengtsson

See Also

```
cdfOrderBy(). applyCdfGroups().
```

 ${\tt cdfSetDimension}$

Sets the dimension of an object

Description

Sets the dimension of an object.

This function is designed to be used with applyCdfGroupFields().

Usage

```
cdfSetDimension(field, dim, ...)
```

Arguments

field An R object.

dim An integer vector.

... Not used.

Value

Returns a list structure of groups.

Author(s)

Henrik Bengtsson

26 compareCdfs

See Also

```
applyCdfGroupFields().
```

compareCdfs

Compares the contents of two CDF files

Description

Compares the contents of two CDF files.

Usage

```
compareCdfs(pathname, other, quick=FALSE, verbose=0, ...)
```

Arguments

pathname The pathname of the first CDF file. other The pathname of the seconds CDF file. quick

If TRUE, only a subset of the units are compared, otherwise all units are com-

verbose An integer. The larger the more details are printed.

Not used. . . .

Details

The comparison is done with an upper-limit memory usage, regardless of the size of the CDFs.

Value

Returns TRUE if the two CDF are equal, otherwise FALSE. If FALSE, the attribute reason contains a string explaining what difference was detected, and the attributes value1 and value2 contain the two objects/values that differs.

Author(s)

Henrik Bengtsson

See Also

convertCdf().

compareCels 27

|--|

Description

Compares the contents of two CEL files.

Usage

```
compareCels(pathname, other, readMap=NULL, otherReadMap=NULL, verbose=0, ...)
```

Arguments

pathname The pathname of the first CEL file.

other The pathname of the seconds CEL file.

readMap An optional read map for the first CEL file.

otherReadMap An optional read map for the second CEL file.

verbose An integer. The larger the more details are printed.

... Not used.

Value

Returns TRUE if the two CELs are equal, otherwise FALSE. If FALSE, the attribute reason contains a string explaining what difference was detected, and the attributes value1 and value2 contain the two objects/values that differs.

Author(s)

Henrik Bengtsson

See Also

convertCel().

28 convertCdf

convertCdf	Converts a CDF into the same CDF but with another format

Description

Converts a CDF into the same CDF but with another format. Currently only CDF files in version 4 (binary/XDA) can be written. However, any input format is recognized.

Usage

```
convertCdf(filename, outFilename, version="4", force=FALSE, ..., .validate=TRUE,
   verbose=FALSE)
```

Arguments

filename	The pathname of the original CDF file.
outFilename	The pathname of the destination CDF file. If the same as the source file, an exception is thrown.
version	The version of the output file format.
force	If FALSE, and the version of the original CDF is the same as the output version, the new CDF will not be generated, otherwise it will.
• • •	Not used.
.validate	If TRUE, a consistency test between the generated and the original CDF is performed. Note that the memory overhead for this can be quite large, because two complete CDF structures are kept in memory at the same time.
verbose	If TRUE, extra details are written while processing.

Value

Returns (invisibly) TRUE if a new CDF was generated, otherwise FALSE.

Benchmarking of ASCII and binary CDFs

Binary CDFs are much faster to read than ASCII CDFs. Here are some example for reading complete CDFs (the difference is even larger when reading CDFs in subsets):

- HG-U133A (22283 units): ASCII 11.7s (9.3x), binary 1.20s (1x).
- Hu6800 (7129 units): ASCII 3.5s (6.1x), binary 0.57s (1x).

Confirmed conversions to binary (XDA) CDFs

The following chip types have been converted using convertCdf() and then verified for correctness using compareCdfs(): ASCII-to-binary: HG-U133A, Hu6800. Binary-to-binary: Test3.

Author(s)

Henrik Bengtsson

convertCel 29

See Also

See compareCdfs() to compare two CDF files. writeCdf().

Examples

```
if (require("AffymetrixDataTestFiles")) {
                                 # START #
chipType <- "Test3"</pre>
cdfFiles <- findCdf(chipType, firstOnly=FALSE)</pre>
cdfFiles <- list(</pre>
 ASCII=grep("ASCII", cdfFiles, value=TRUE),
 XDA=grep("XDA", cdfFiles, value=TRUE)
)
outFile <- file.path(tempdir(), sprintf("%s.cdf", chipType))</pre>
convertCdf(cdfFiles$ASCII, outFile, verbose=TRUE)
}
                                 # STOP #
```

convertCel

Converts a CEL into the same CEL but with another format

Description

Converts a CEL into the same CEL but with another format. Currently only CEL files in version 4 (binary/XDA) can be written. However, any input format is recognized.

Usage

```
convertCel(filename, outFilename, readMap=NULL, writeMap=NULL, version="4",
   newChipType=NULL, ..., .validate=FALSE, verbose=FALSE)
```

Arguments

filename The pathname of the original CEL file.

outFilename The pathname of the destination CEL file. If the same as the source file, an

exception is thrown.

readMap An optional read map for the input CEL file.
writeMap An optional write map for the output CEL file.

version The version of the output file format.

30 convertCel

newChipType (Only for advanced users who fully understands the Affymetrix CEL file format!) An optional string for overriding the chip type (label) in the CEL file header.

... Not used.

.validate If TRUE, a consistency test between the generated and the original CEL is performed.

verbose If TRUE, extra details are written while processing.

Value

Returns (invisibly) TRUE if a new CEL was generated, otherwise FALSE.

Benchmarking of ASCII and binary CELs

Binary CELs are much faster to read than ASCII CELs. Here are some example for reading complete CELs (the difference is even larger when reading CELs in subsets):

• To do

WARNING: Changing the chip type label

The newChipType argument changes the label in the part of DAT header that specifies the chip type of the CEL file. Note that it does not change anything else in the CEL file. This type of relabeling is valid for updating the chip type *label* of CEL files that where generated during, say, an "Early Access" period leading to a different chip type label than what more recent CEL files of the same physical chip type have.

Author(s)

Henrik Bengtsson

See Also

```
createCel().
```

Examples

copyCel 31

copyCel

Copies a CEL file

Description

Copies a CEL file.

The file must be a valid CEL file, if not an exception is thrown.

Usage

```
copyCel(from, to, overwrite=FALSE, ...)
```

Arguments

from The filename of the CEL file to be copied.

to The filename of destination file.

overwrite If FALSE and the destination file already exists, an exception is thrown, otherwise

not.

... Not used.

Value

Return TRUE if file was successfully copied, otherwise FALSE.

Author(s)

Henrik Bengtsson

```
isCelFile().
```

32 createCel

createCel Creates an empty CEL file	createCel	Creates an empty CEL file	
-------------------------------------	-----------	---------------------------	--

Description

Creates an empty CEL file.

Usage

createCel(filename, header, nsubgrids=0, overwrite=FALSE, ..., cdf=NULL, verbose=FALSE)

Arguments

filename The filename of the CEL file to be created. header A list structure describing the CEL header, similar to the structure returned by readCelHeader(). This header can be of any CEL header version. overwrite If FALSE and the file already exists, an exception is thrown, otherwise the file is created. The number of subgrids. nsubgrids Not used. . . . cdf (optional) The pathname of a CDF file for the CEL file to be created. If given, the CEL header (argument header) is validated against the CDF header, otherwise not. If TRUE, a CDF file is located automatically based using findCdf(header\$chiptype). verbose An integer specifying how much verbose details are outputted.

Details

Currently only binary (v4) CEL files are supported. The current version of the method does not make use of the Fusion SDK, but its own code to create the CEL file.

Value

Returns (invisibly) the pathname of the file created.

Redundant fields in the CEL header

There are a few redundant fields in the CEL header. To make sure the CEL header is consistent, redundant fields are cleared and regenerated. For instance, the field for the total number of cells is calculated from the number of cell rows and columns.

Author(s)

Henrik Bengtsson

findCdf 33

Examples

```
if (require("AffymetrixDataTestFiles")) {
                                       # START #
# Search for first available ASCII CEL file
path <- system.file("rawData", package="AffymetrixDataTestFiles")</pre>
files <- findFiles(pattern="[.](cel|CEL)$", path=path, recursive=TRUE, firstOnly=FALSE)
files <- grep("ASCII", files, value=TRUE)</pre>
file <- files[1]
# Read the CEL header
hdr <- readCelHeader(file)</pre>
# Assert that we found an ASCII CEL file, but any will do
stopifnot(hdr$version == 3)
# Create a CEL v4 file of the same chip type
outFile <- file.path(tempdir(), "zzz.CEL")</pre>
if (file.exists(outFile))
 file.remove(outFile)
createCel(outFile, hdr, overwrite=TRUE)
str(readCelHeader(outFile))
# Verify correctness by update and re-read a few cells
intensities <- as.double(1:100)</pre>
indices <- seq(along=intensities)</pre>
updateCel(outFile, indices=indices, intensities=intensities)
value <- readCel(outFile, indices=indices)$intensities</pre>
stopifnot(identical(intensities, value))
```

findCdf

Search for CDF files in multiple directories

Description

Search for CDF files in multiple directories.

34 findCdf

Usage

```
findCdf(chipType=NULL, paths=NULL, recursive=TRUE, pattern="[.](c|C)(d|D)(f|F)$", ...)
```

Arguments

chipType A character string of the chip type to search for.

paths A character vector of paths to be searched. The current directory is always searched at the beginning. If NULL, default paths are searched. For more details, see below.

recursive If TRUE, directories are searched recursively.

recursive If TRUE, directories are searched recursively.

pattern A regular expression file name pattern to match.

Additional arguments passed to findFiles().

Details

Note, the current directory is always searched first, but never recursively (unless it is added to the search path explicitly). This provides an easy way to override other files in the search path.

If paths is NULL, then a set of default paths are searched. The default search path constitutes:

```
    getOption("AFFX_CDF_PATH")
    Sys.getenv("AFFX_CDF_PATH")
```

One of the easiest ways to set system variables for R is to set them in an .Renviron file, e.g.

```
# affxparser: Set default CDF path
AFFX_CDF_PATH=${AFFX_CDF_PATH};M:/Affymetrix_2004-100k_trios/cdf
AFFX_CDF_PATH=${AFFX_CDF_PATH};M:/Affymetrix_2005-500k_data/cdf
```

See Startup for more details.

Value

Returns a vector of the full pathnames of the files found.

Author(s)

Henrik Bengtsson

See Also

This method is used internally by readCelUnits() if the CDF file is not specified.

findFiles 35

Examples

```
if (require("AffymetrixDataTestFiles")) {
                                 # START #
# Find a specific CDF file
cdfFile <- findCdf("Mapping10K_Xba131")</pre>
print(cdfFile)
# Find the first CDF file (no matter what it is)
cdfFile <- findCdf()</pre>
print(cdfFile)
# Find all CDF files in search path and display their headers
cdfFiles <- findCdf(firstOnly=FALSE)</pre>
for (cdfFile in cdfFiles) {
 cat("=======\n")
 hdr <- readCdfHeader(cdfFile)</pre>
 str(hdr)
}
```

findFiles

Finds one or several files in multiple directories

Description

Finds one or several files in multiple directories.

Usage

```
findFiles(pattern=NULL, paths=NULL, recursive=FALSE, firstOnly=TRUE, allFiles=TRUE, ...)
```

Arguments

```
pattern A regular expression file name pattern to match.

A character vector of paths to be searched.

recursive If TRUE, the directory structure is searched breath-first, in lexicographic order.

firstOnly If TRUE, the method returns as soon as a matching file is found, otherwise not.

If FALSE, files and directories starting with a period will be skipped, otherwise not.

Arguments passed to list.files().
```

36 invertMap

Value

Returns a vector of the full pathnames of the files found.

Paths

The paths argument may also contain paths specified as semi-colon (";") separated paths, e.g. "/usr/;usr/bin/;.;".

Windows Shortcut links

If package **R.utils** is available and loaded, Windows Shortcut links (*.lnk) are recognized and can be used to imitate links to directories elsewhere. For more details, see filePath.

Author(s)

Henrik Bengtsson

invertMap

Inverts a read or a write map

Description

Inverts a read or a write map.

Usage

```
invertMap(map, ...)
```

Arguments

```
map An integer vector.
... Not used.
```

Details

An map is defined to be a vector of n with unique finite values in [1, n]. Finding the inverse of a map is the same as finding the rank of each element, cf. order(). However, this method is much faster, because it utilizes the fact that all values are unique and in [1, n]. Moreover, for any map it holds that taking the inverse twice will result in the same map.

Value

Returns an integer vector.

Author(s)

Henrik Bengtsson

isCelFile 37

See Also

To generate an optimized write map for a CDF file, see readCdfUnitsWriteMap().

Examples

```
set.seed(1)
# Simulate a read map for a chip with 1.2 million cells
nbrOfCells <- 1200000
readMap <- sample(nbr0fCells)</pre>
# Get the corresponding write map
writeMap <- invertMap(readMap)</pre>
# A map inverted twice should be equal itself
stopifnot(identical(invertMap(writeMap), readMap))
# Another example illustrating that the write map is the
# inverse of the read map
idx <- sample(nbr0fCells, size=1000)</pre>
stopifnot(identical(writeMap[readMap[idx]], idx))
# invertMap() is much faster than order()
t1 <- system.time(invertMap(readMap))[3]</pre>
cat(sprintf("invertMap() : %5.2fs [ 1.00x]\n", t1))
t2 <- system.time(writeMap2 <- sort.list(readMap, na.last=NA, method="quick"))[3]
cat(sprintf("'quick sort' : \%5.2fs \ [\%5.2fx]\n", \ t2, \ t2/t1))
stopifnot(identical(writeMap, writeMap2))
t3 <- system.time(writeMap2 <- order(readMap))[3]
cat(sprintf("order()
                          : %5.2fs [%5.2fx]\n", t3, t3/t1))
stopifnot(identical(writeMap, writeMap2))
# Clean up
rm(nbrOfCells, idx, readMap, writeMap, writeMap2)
```

isCelFile

Checks if a file is a CEL file or not

Description

Checks if a file is a CEL file or not.

Usage

```
isCelFile(filename, ...)
```

Arguments

```
filename A filename.
... Not used.
```

Value

Returns TRUE if a CEL file, otherwise FALSE. ASCII (v3), binary (v4;XDA), and binary (CCG v1;Calvin) CEL files are recognized. If file does not exist, an exception is thrown.

Author(s)

Henrik Bengtsson

See Also

```
readCel(), readCelHeader(), readCelUnits().
```

```
parseDatHeaderString Parses a DAT header string
```

Description

Parses a DAT header string.

Usage

```
parseDatHeaderString(header, timeFormat="%m/%d/%y %H:%M:%S", ...)
```

Arguments

header A character string.

timeFormat The format string used to parse the timestamp. For more details, see strptime().

If NULL, no parsing is done.

. . . Not used.

Value

Returns named list structure.

Author(s)

Henrik Bengtsson

See Also

```
readCelHeader().
```

readBpmap 39

|--|--|

Description

Parses (parts of) a Bpmap (binary probe mapping) file from Affymetrix.

Usage

```
readBpmap(filename, seqIndices = NULL, readProbeSeq = TRUE, readSeqInfo
= TRUE, readPMXY = TRUE, readMMXY = TRUE, readStartPos = TRUE,
readCenterPos = FALSE, readStrand = TRUE, readMatchScore = FALSE,
readProbeLength = FALSE, verbose = 0)
readBpmapHeader(filename)
readBpmapSeqinfo(filename, seqIndices = NULL, verbose = 0)
```

Arguments

_	
filename	The filename as a character.
seqIndices	\boldsymbol{A} vector of integers, detailing the indices of the sequences being read. If NULL, the entire file is being read.
readProbeSeq	Do we read the probe sequences.
readSeqInfo	Do we read the sequence information (a list containing information such as sequence name, number of hits etc.)
readPMXY	Do we read the (x,y) coordinates of the PM-probes.
readMMXY	Do we read the (x,y) coordinates of the MM-probes (only relevant if the file has MM information)
readStartPos	Do we read the start position of the probes.
readCenterPos	Do we return the start position of the probes.
readStrand	Do we return the strand of the hits.
readMatchScore readProbeLength	Do we return the matchscore.
	Doe we return the probelength.
verbose	How verbose do we want to be.

Details

readBpmap reads a BPMAP file, which is a binary file containing information about a given probe's location in a sequence. Here sequence means some kind of reference sequence, typically a chromosome or a scaffold. readBpmapHeader reads the header of the BPMAP file, and readBpmapSeqinfo reads the sequence info of the sequences (so this function is merely a convenience function).

40 readCcg

Value

For readBpmap: A list of lists, one list for every sequence read. The components of the sequence lists, depends on the argument of the function call. For readBpmapheader a list with two components version and numSequences. For readBpmapSeqinfo a list of lists containing the sequence info.

Author(s)

Kasper Daniel Hansen

See Also

tpmap2bpmap for information on how to write Bpmap files.

readCcg

Reads an Affymetrix Command Console Generic (CCG) Data file

Description

Reads an Affymetrix Command Console Generic (CCG) Data file. The CCG data file format is also known as the Calvin file format.

Usage

```
readCcg(pathname, verbose=0, .filter=NULL, ...)
```

Arguments

pathname The pathname of the CCG file.

verbose An integer specifying the verbose level. If 0, the file is parsed quietly. The

higher numbers, the more details.

.filter A list.... Not used.

Details

Note, the current implementation of this methods does not utilize the Affymetrix Fusion SDK library. Instead, it is implemented in R from the file format definition [1].

Value

A named list structure consisting of ...

readCcg 41

About the CCG file format

A CCG file, consists of a "file header", a "generic data header", and "data" section, as outlined here:

- File Header
- Generic Data Header (for the file)
 - 1. Generic Data Header (for the files 1st parent)
 - (a) Generic Data Header (for the files 1st parents 1st parent)
 - (b) Generic Data Header (for the files 1st parents 2nd parent)
 - (c) ..
 - (d) Generic Data Header (for the files 1st parents Mth parent)
 - 2. Generic Data Header (for the files 2nd parent)
 - 3. ...
 - 4. Generic Data Header (for the files Nth parent)
- Data
 - 1. Data Group #1
 - (a) Data Set #1
 - Parameters
 - Column definitions
 - Matrix of data
 - (b) Data Set #2
 - (c) ...
 - (d) Data Set #L
 - 2. Data Group #2
 - 3. ...
 - 4. Data Group #K

Author(s)

Henrik Bengtsson

References

[1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, April, 2006. http://www.affymetrix.com/support/developer/

See Also

```
readCcgHeader(). readCdfUnits().
```

42 readCcgHeader

(CCG) file

Description

Reads an the header of an Affymetrix Command Console Generic (CCG) file.

Usage

```
readCcgHeader(pathname, verbose=0, .filter=list(fileHeader = TRUE, dataHeader = TRUE),
...)
```

Arguments

pathname	The pathname of the CCG file.
verbose	An integer specifying the verbose level. If 0, the file is parsed quietly. The higher numbers, the more details.
.filter	A list.
	Not used.

Details

Note, the current implementation of this methods does not utilize the Affymetrix Fusion SDK library. Instead, it is implemented in R from the file format definition [1].

Value

A named list structure consisting of ...

Author(s)

Henrik Bengtsson

References

```
[1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, April, 2006. http://www.affymetrix.com/support/developer/
```

See Also

```
readCcg().
```

readCdf 43

readCdf Parsing a CDF file using Affymetrix Fusion SDK
Tarsing a CDF file using Affymenta Fusion SDK

Description

Parsing a CDF file using Affymetrix Fusion SDK. This function parses a CDF file using the Affymetrix Fusion SDK. *This function will most likely be replaced by the more general* readCdfUnits() *function*.

Usage

```
readCdf(filename, units=NULL,
    readXY=TRUE, readBases=TRUE,
    readIndexpos=TRUE, readAtoms=TRUE,
    readUnitType=TRUE, readUnitDirection=TRUE,
    readUnitNumber=TRUE, readUnitAtomNumbers=TRUE,
    readGroupAtomNumbers=TRUE, readGroupDirection=TRUE,
    readIndices=FALSE, readIsPm=FALSE,
    stratifyBy=c("nothing", "pmmm", "pm", "mm"),
    verbose=0)
```

not are retrieved, otherwise not.

Arguments

readIsPm

filename The filename of the CDF file. An integer vector of unit indices specifying which units to be read. If NULL, units all units are read. readXY If TRUE, cell row and column (x,y) coordinates are retrieved, otherwise not. readBases If TRUE, cell P and T bases are retrieved, otherwise not. readIndexpos If TRUE, cell indexpos are retrieved, otherwise not. readUnitType If TRUE, unit types are retrieved, otherwise not. readUnitDirection If TRUE, unit directions are retrieved, otherwise not. readUnitNumber If TRUE, unit numbers are retrieved, otherwise not. readUnitAtomNumbers If TRUE, unit atom numbers are retrieved, otherwise not. readGroupAtomNumbersIf TRUE, group atom numbers are retrieved, otherwise not. readGroupDirection If TRUE, group directions are retrieved, otherwise not. readIndices If TRUE, cell indices calculated from the row and column (x,y) coordinates are retrieved, otherwise not. Note that these indices are *one-based*.

If TRUE, cell flags indicating whether the cell is a perfect-match (PM) probe or

44 readCdf

stratifyBy A character string specifying which and how elements in group fields are re-

turned. If "nothing", elements are returned as is, i.e. as vectors. If "pm"/"mm", only elements corresponding to perfect-match (PM) / mismatch (MM) probes are returned (as vectors). If "pmmm", elements are returned as a matrix where the first row holds elements corresponding to PM probes and the second corresponding to MM probes. Note that in this case, it is assumed that there are equal number of PMs and MMs; if not, an error is generated. Moreover, the PMs and MMs may not even be paired, i.e. there is no guarantee that the two elements in

a column corresponds to a PM-MM pair.

verbose An integer specifying the verbose level. If 0, the file is parsed quietly. The

higher numbers, the more details.

Value

A list with one component for each unit. Every component is again a list with three components

groups This is again a list with one component for each group (also called block). The

information on each group is a list with 5 components, x, y, pbase, tbase,

expos.

type type of the unit.
direction direction of the unit.

Cell indices are one-based

Note that in **affxparser** all *cell indices* are by convention *one-based*, which is more convenient to work with in R. For more details on one-based indices, see 2. Cell coordinates and cell indices.

Note

This version of the function does not return information on the QC probes. This will be added in a (near) future release. In addition we expect the header to be part of the returned object.

So expect changes to the structure of the value of the function in next release. Please contact the developers for details.

Author(s)

James Bullard and Kasper Daniel Hansen.

References

[1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, June 14, 2005. http://www.affymetrix.com/support/developer/

See Also

It is recommended to use readCdfUnits() instead of this method. readCdfHeader() for getting the header of a CDF file.

readCdfCellIndices 45

readCdfCellIndices	Reads (one-based) cell indices of units (probesets) in an Affymetrix CDF file
--------------------	---

Description

Reads (one-based) cell indices of units (probesets) in an Affymetrix CDF file.

Usage

```
readCdfCellIndices(filename, units=NULL, stratifyBy=c("nothing", "pmmm", "pm", "mm"),
   verbose=0)
```

Arguments

O .	
filename	The filename of the CDF file.
units	An integer vector of unit indices specifying which units to be read. If NULL, all units are read.
stratifyBy	A character string specifying which and how elements in group fields are returned. If "nothing", elements are returned as is, i.e. as vectors. If "pm"/"mm", only elements corresponding to perfect-match (PM) / mismatch (MM) probes are returned (as vectors). If "pmmm", elements are returned as a matrix where the first row holds elements corresponding to PM probes and the second corresponding to MM probes. Note that in this case, it is assumed that there are equal number of PMs and MMs; if not, an error is generated. Moreover, the PMs and MMs may not even be paired, i.e. there is no guarantee that the two elements in a column corresponds to a PM-MM pair.
verbose	An integer specifying the verbose level. If 0, the file is parsed quietly. The higher numbers, the more details.

Value

A named list where the names corresponds to the names of the units read. Each unit element of the list is in turn a list structure with one element groups which in turn is a list. Each group element in groups is a list with a single field named indices. Thus, the structure is

46 readCdfDataFrame

```
+- unit #2
.
+- unit #J
```

This is structure is compatible with what readCdfUnits() returns.

Note that these indices are one-based.

Cell indices are one-based

Note that in **affxparser** all *cell indices* are by convention *one-based*, which is more convenient to work with in R. For more details on one-based indices, see 2. Cell coordinates and cell indices.

Author(s)

Henrik Bengtsson

See Also

readCdfUnits().

readCdfDataFrame	Reads units (probesets) from an Affymetrix CDF file	
------------------	---	--

Description

Reads units (probesets) from an Affymetrix CDF file. Gets all or a subset of units (probesets).

Usage

readCdfDataFrame(filename, units=NULL, groups=NULL, cells=NULL, fields=NULL, drop=TRUE, verbose=0)

Arguments

filename	The filename of the CDF file.
units	An integer vector of unit indices specifying which units to be read. If NULL, all are read.
groups	An integer vector of group indices specifying which groups to be read. If NULL, all are read.
cells	An integer vector of cell indices specifying which cells to be read. If NULL, all are read.
fields	A character vector specifying what fields to read. If NULL, all unit, group and cell fields are returned.
drop	If TRUE and only one field is read, then a vector (rather than a single-column data.frame) is returned.
verbose	An integer specifying the verbose level. If 0, the file is parsed quietly. The higher numbers, the more details.

readCdfGroupNames 47

Value

An NxK data. frame or a vector of length N.

Author(s)

Henrik Bengtsson

References

[1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, June 14, 2005. http://www.affymetrix.com/support/developer/

See Also

For retrieving the CDF as a list structure, see readCdfUnits.

Examples

```
if (require("AffymetrixDataTestFiles")) {
# Find any CDF file
cdfFile <- findCdf()</pre>
units <- 101:120
fields <- c("unit", "unitName", "group", "groupName", "cell")</pre>
df <- readCdfDataFrame(cdfFile, units=units, fields=fields)</pre>
stopifnot(identical(sort(unique(df$unit)), units))
fields <- c("unit", "unitName", "unitType")</pre>
fields <- c(fields, "group", "groupName")
fields <- c(fields, "x", "y", "cell", "pbase", "tbase")</pre>
df <- readCdfDataFrame(cdfFile, units=units, fields=fields)</pre>
stopifnot(identical(sort(unique(df$unit)), units))
# STOP #
```

readCdfGroupNames

Reads group names for a set of units (probesets) in an Affymetrix CDF file

48 readCdfHeader

Description

Reads group names for a set of units (probesets) in an Affymetrix CDF file.

This is for instance useful for SNP arrays where the nucleotides used for the A and B alleles are the same as the group names.

Usage

readCdfGroupNames(filename, units=NULL, truncateGroupNames=TRUE, verbose=0)

Arguments

filename The filename of the CDF file.

units An integer vector of unit indices specifying which units to be read. If NULL,

all units are read.

truncateGroupNames

A logical variable indicating whether unit names should be stripped from the

beginning of group names.

verbose An integer specifying the verbose level. If 0, the file is parsed quietly. The

higher numbers, the more details.

Value

A named list structure where the names of the elements are the names of the units read. Each element is a character vector with group names for the corresponding unit.

Author(s)

Henrik Bengtsson

See Also

readCdfUnits().

readCdfHeader

Reads the header associated with an Affymetrix CDF file

Description

Reads the header of an Affymetrix CDF file using the Fusion SDK.

Usage

readCdfHeader(filename)

Arguments

filename name of the CDF file.

readCdfIsPm 49

Value

A named list with the following components:

rows the number of rows on the chip.

cols the number of columns on the chip.

probesets the number of probesets on the chip.

qcprobesets the number of QC probesets on the chip.

reference the reference sequence (this component only exists for resequencing chips).

chiptype the type of the chip.
filename the name of the cdf file.

Author(s)

James Bullard and Kasper Daniel Hansen

See Also

```
readCdfUnits().
```

Examples

```
for (zzz in 0) {
# Find any CDF file
cdfFile <- findCdf()
if (is.null(cdfFile))
    break

header <- readCdfHeader(cdfFile)
print(header)
} # for (zzz in 0)</pre>
```

readCdfIsPm

Checks if cells in a CDF file are perfect-match probes or not

Description

Checks if cells in a CDF file are perfect-match probes or not.

Usage

```
readCdfIsPm(filename, units=NULL, verbose=0)
```

Arguments

filename of the CDF file.

units An integer vector of unit indices specifying which units to be read. If NULL,

all units are read.

verbose An integer specifying the verbose level. If 0, the file is parsed quietly. The

higher numbers, the more details.

Value

A named list of named logical vectors. The name of the list elements are unit names and the names of the logical vector are group names.

Author(s)

Henrik Bengtsson

readCdfNbrOfCellsPerUnitGroup

Gets the number of cells (probes) that each group of each unit in a

CDF file

Description

Gets the number of cells (probes) that each group of each unit in a CDF file.

Usage

readCdfNbrOfCellsPerUnitGroup(filename, units=NULL, verbose=0)

Arguments

filename of the CDF file.

units An integer vector of unit indices specifying which units to be read. If NULL,

all units are read.

verbose An integer specifying the verbose level. If 0, the file is parsed quietly. The

higher numbers, the more details.

Value

A named list of named integer vectors. The name of the list elements are unit names and the names of the integer vector are group names.

Author(s)

Henrik Bengtsson

```
if (require("AffymetrixDataTestFiles")) {
cdfFile <- findCdf("Mapping10K_Xba131")</pre>
groups <- readCdfNbrOfCellsPerUnitGroup(cdfFile)</pre>
# Number of units read
print(length(groups))
## 11564
# Details on two units
print(groups[56:57])
## $`SNP_A-1516438`
## SNP_A-1516438C SNP_A-1516438T SNP_A-1516438C SNP_A-1516438T
            10
                                                      10
                          10
                                        10
##
## $`SNP_A-1508602`
## SNP_A-1508602A SNP_A-1508602G SNP_A-1508602A SNP_A-1508602G
            10
                          10
# Number of groups with different number of cells
print(table(unlist(groups)))
##
    10
          60
## 46240
# Number of cells per unit
nbr0fCellsPerUnit <- unlist(lapply(groups, FUN=sum))</pre>
print(table(nbrOfCellsPerUnit))
nbrOfCellsPerUnit
##
     40
## 11560
# Number of groups per unit
nbrOfGroupsPerUnit <- unlist(lapply(groups, FUN=length))</pre>
# Details on a few units
print(nbr0fGroupsPerUnit[20:30])
## SNP_A-1512666 SNP_A-1512740 SNP_A-1512132 SNP_A-1516082 SNP_A-1511962
            4
                        4
                                     4
                                                   4
## SNP_A-1515637 SNP_A-1515878 SNP_A-1518789 SNP_A-1518296 SNP_A-1519701
##
            4
                          4
                                      4
                                                   4
## SNP_A-1511743
##
# Number of units for each unique number of groups
print(table(nbrOfGroupsPerUnit))
```

52 readCdfQc

```
## nbrOfGroupsPerUnit
     1
##
     4 11560
x <- list()
for (size in unique(nbr0fGroupsPerUnit)) {
 subset <- groups[nbr0fGroupsPerUnit==size]</pre>
 t <- matrix(unlist(subset), nrow=size)</pre>
 colnames(t) <- names(subset)</pre>
 x[[as.character(size)]] <- t
 rm(subset, t)
# Check if there are any quartet units where the number
# of cells in Group 1 & 2 or Group 3 & 4 does not have
# the same number of cells.
# Group 1 & 2
print(sum(x[["4"]][1,]-x[["4"]][2,] != 0))
# Group 3 & 4
print(sum(x[["4"]][3,]-x[["4"]][4,] != 0))
# STOP #
```

readCdfQc

Reads the QC units of CDF file

Description

Reads the QC units of CDF file.

Usage

```
readCdfQc(filename, units = NULL, verbose = 0)
```

Arguments

filename name of the CDF file.

units The QC unit indices as a vector of integers. NULL indicates that all units should

be read.

verbose how verbose should the output be. 0 means no output, with higher numbers

being more verbose.

Value

A list with one component for each QC unit.

readCdfUnitNames 53

Author(s)

Kasper Daniel Hansen

See Also

readCdf().

readCdfUnitNames

Reads unit (probeset) names from an Affymetrix CDF file

Description

Gets the names of all or a subset of units (probesets) in an Affymetrix CDF file. This can be used to get a map between unit names an the internal unit indices used by the CDF file.

Usage

```
readCdfUnitNames(filename, units=NULL, verbose=0)
```

Arguments

filename of the CDF file.

units An integer vector of unit indices specifying which units to be read. If NULL,

all units are read.

verbose An integer specifying the verbose level. If 0, the file is parsed quietly. The

higher numbers, the more details.

Value

A character vector of unit names.

Author(s)

```
Henrik Bengtsson (http://www.braju.com/R/)
```

See Also

```
readCdfUnits().
```

```
## Not run: See help(readCdfUnits) for an example
```

54 readCdfUnits

readCdfUnits	Reads units (probesets) from an Affymetrix CDF file	

Description

Reads units (probesets) from an Affymetrix CDF file. Gets all or a subset of units (probesets).

Usage

```
readCdfUnits(filename, units=NULL, readXY=TRUE, readBases=TRUE, readExpos=TRUE,
readType=TRUE, readDirection=TRUE, stratifyBy=c("nothing", "pmmm", "pm", "mm"),
readIndices=FALSE, verbose=0)
```

Arguments

Bu		
filename	The filename of the CDF file.	
units	An integer vector of unit indices specifying which units to be read. If NULL, all units are read.	
readXY	If TRUE, cell row and column (x,y) coordinates are retrieved, otherwise not.	
readBases	If TRUE, cell P and T bases are retrieved, otherwise not.	
readExpos	If TRUE, cell "expos" values are retrieved, otherwise not.	
readType	If TRUE, unit types are retrieved, otherwise not.	
readDirection	If TRUE, unit and group directions are retrieved, otherwise not.	
stratifyBy	A character string specifying which and how elements in group fields are returned. If "nothing", elements are returned as is, i.e. as vectors. If "pm"/"mm", only elements corresponding to perfect-match (PM) / mismatch (MM) probes are returned (as vectors). If "pmmm", elements are returned as a matrix where the first row holds elements corresponding to PM probes and the second corresponding to MM probes. Note that in this case, it is assumed that there are equal number of PMs and MMs; if not, an error is generated. Moreover, the PMs and MMs may not even be paired, i.e. there is no guarantee that the two elements in a column corresponds to a PM-MM pair.	
readIndices	If TRUE, cell indices <i>calculated</i> from the row and column (x,y) coordinates are retrieved, otherwise not. Note that these indices are <i>one-based</i> .	

Value

verbose

A named list where the names corresponds to the names of the units read. Each element of the list is in turn a list structure with three components:

higher numbers, the more details.

An integer specifying the verbose level. If 0, the file is parsed quietly. The

readCdfUnits 55

groups	A list with one component for each group (also called block). The information on each group is a list of up to seven components: x, y, pbase, tbase, expos, indices, and direction. All fields but the latter have the same number of values as there are cells in the group. The latter field has only one value indicating the direction for the whole group.
type	An integer specifying the type of the unit, where 1 is "expression", 2 is "genotyping", 3 is "CustomSeq", and 4 "tag".
direction	An integer specifying the direction of the unit, which defines if the probes are interrogating the sense or the anti-sense target, where 0 is "no direction", 1 is "sense", and 2 is "anti-sense".

Cell indices are one-based

Note that in **affxparser** all *cell indices* are by convention *one-based*, which is more convenient to work with in R. For more details on one-based indices, see 2. Cell coordinates and cell indices.

Author(s)

James Bullard and Kasper Daniel Hansen. Modified by Henrik Bengtsson to read any subset of units and/or subset of parameters, to stratify by PM/MM, and to return cell indices.

References

[1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, June 14, 2005. http://www.affymetrix.com/support/developer/

See Also

```
readCdfCellIndices().
```

readCdfUnitsWriteMap Generates an Affymetrix cell-index write map from a CDF file

Description

Generates an Affymetrix cell-index write map from a CDF file.

The purpose of this method is to provide a re-ordering of cell elements such that cells in units (probesets) can be stored in contiguous blocks. When reading cell elements unit by unit, minimal file re-position is required resulting in a faster reading.

Note: At the moment does this package not provide methods to write/reorder CEL files. In the meanwhile, you have to write and re-read using your own file format. That's not too hard using writeBin() and readBin().

Usage

```
readCdfUnitsWriteMap(filename, units=NULL, ..., verbose=FALSE)
```

Arguments

filename The pathname of the CDF file.

units An integer vector of unit indices specifying which units to listed first. All

other units are added in order at the end. If NULL, units are in order.

... Additional arguments passed to readCdfUnits().

verbose Either a logical, a numeric, or a Verbose object specifying how much ver-

bose/debug information is written to standard output. If a Verbose object, how detailed the information is is specified by the threshold level of the object. If a numeric, the value is used to set the threshold of a new Verbose object. If TRUE, the threshold is set to -1 (minimal). If FALSE, no output is written (and neither

is the **R.utils** package required).

Value

A integer vector which is a write map.

Author(s)

Henrik Bengtsson

See Also

To invert maps, see invertMap(). readCel() and readCelUnits().

```
if (require("AffymetrixDataTestFiles")) {
# Find any CDF file
cdfFile <- findCdf()</pre>
# Create a cell-index map (for writing)
writeMap <- readCdfUnitsWriteMap(cdfFile)</pre>
# Inverse map to be used to read cell elements such that, when read
# read unit by unit, they are read much faster.
readMap <- invertMap(writeMap)</pre>
# Validate the two maps
stopifnot(identical(readMap[writeMap], 1:length(readMap)))
cat("Summary of the \"randomness\" of the cell indices:\n")
moves <- diff(readMap) - 1
cat(sprintf("Number of unnecessary file re-positioning: %d (%.1f%%)\n",
              sum(moves != 0), 100*sum(moves != 0)/length(moves)))
cat(sprintf("Extra positioning: %.1fGb\n", sum(abs(moves))/1024^3))
smallMoves <- moves[abs(moves) <= 25];</pre>
largeMoves <- moves[abs(moves) > 25];
layout(matrix(1:2))
main <- "Non-signed file moves required in unorded file"</pre>
hist(smallMoves, nclass=51, main=main, xlab="moves <=25 bytes")</pre>
hist(largeMoves, nclass=101, main="", xlab="moves >25 bytes")
# Clean up
layout(1)
rm(cdfFile, readMap, writeMap, moves, smallMoves, largeMoves, main)
# STOP #
```

```
# Function to read Affymetrix probeset annotations
readAffymetrixProbesetAnnotation <- function(pathname, ...) {</pre>
 # Get headers
 header <- scan(pathname, what="character", sep=",", quote="\"",
                                                    quiet=TRUE, nlines=1);
 # Read only a subset of columns (unique to this example)
 cols <- c("Probe Set ID"="probeSet",</pre>
           "Chromosome"="chromosome",
           "Physical Position"="physicalPosition",
           "dbSNP RS ID"="dbSnpId");
 colClasses <- rep("NULL", length(header));</pre>
 colClasses[header %in% names(cols)] <- "character";</pre>
 # Read the data (this is what takes time)
 df <- read.table(pathname, colClasses=colClasses, header=TRUE, sep=",",</pre>
        quote="\"", na.strings="---", strip.white=TRUE, check.names=FALSE,
                blank.lines.skip=FALSE, fill=FALSE, comment.char="", ...);
 # Re-order columns
 df <- df[,match(names(cols),colnames(df))];</pre>
 colnames(df) <- cols;</pre>
 # Use "Probe Set ID" as rownames. Note that if we use 'row.names=1'
 # or similar something goes wrong. /HB 2006-03-06
 rownames(df) <- df[[1]];</pre>
 df <- df[,-1];</pre>
 # Change types of columns
 df[[1]] <- factor(df[[1]], levels=c(1:22,"X","Y",NA), ordered=TRUE);</pre>
 df[[2]] <- as.integer(df[[2]]);</pre>
 df;
} # readAffymetrixProbesetAnnotation()
# Main
for (zz in 1) {
# Chip to be remapped
chipType <- "Mapping50K_Xba240"</pre>
annoFile <- paste(chipType, "_annot.csv", sep="")</pre>
cdfFile <- findCdf(chipType)</pre>
if (is.null(cdfFile) || !file.exists(annoFile))
 break;
```

readCel 59

```
# Read SNP location details
snpInfo <- readAffymetrixProbesetAnnotation(annoFile)</pre>
# Order by chromsome and then physical position
o <- order(snpInfo[[1]], snpInfo[[2]])</pre>
snpInfo <- snpInfo[o,]</pre>
rm(o)
# Read unit names in CDF file
unitNames <- readCdfUnitNames(cdfFile)</pre>
# The CDF unit indices sorted by chromsomal position
units <- match(rownames(snpInfo), unitNames)</pre>
# ...and cell indices in the same order
writeMap <- readCdfUnitsWriteMap(cdfFile, units=units)</pre>
# Inverse map to be used to write cell elements such that, if they
# later are read unit by unit, they are read in contiguous blocks.
readMap <- invertMap(writeMap)</pre>
# Clean up
rm(chipType, annoFile, cdfFile, snpInfo, unitNames, units, readMap, writeMap)
} # for (zz in 1)
```

readCel

Reads an Affymetrix CEL file

Description

This function reads all or a subset of the data in an Affymetrix CEL file.

Usage

```
readCel(filename,
    indices = NULL,
    readHeader = TRUE,
    readXY = FALSE, readIntensities = TRUE,
    readStdvs = FALSE, readPixels = FALSE,
    readOutliers = TRUE, readMasked = TRUE,
    readMap = NULL,
    verbose = 0,
    .checkArgs = TRUE)
```

60 readCel

Arguments

filename the name of the CEL file.

indices a vector of indices indicating which features to read. If the argument is NULL all

features will be returned.

readXY a logical: will the (x,y) coordinates be returned.

readIntensities

a logical: will the intensities be returned.

readStdvs a logical: will the standard deviations be returned.
readPixels a logical: will the number of pixels be returned.

readOutliers a logical: will the outliers be return.

readMasked a logical: will the masked features be returned.
readHeader a logical: will the header of the file be returned.

readMap A vector remapping cell indices to file indices. If NULL, no mapping is used.

verbose how verbose do we want to be. 0 is no verbosity, higher numbers mean more

verbose output. At the moment the values 0, 1 and 2 are supported.

. checkArgs If TRUE, the arguments will be validated, otherwise not. Warning: This should

only be used if the arguments have been validated elsewhere!

Value

A CEL files consists of a *header*, a set of *cell values*, and information about *outliers* and masked cells.

The cell values, which are values extract for each cell (aka feature or probe), are the (x,y) coordinate, intensity and standard deviation estimates, and the number of pixels in the cell. If readIndices=NULL, cell values for all cells are returned, Only cell values specified by argument readIndices are returned.

This value returns a named list with components described below:

header The header of the CEL file. Equivalent to the output from readCelHeader, see

the documentation for that function.

x, y (cell values) Two integer vectors containing the x and y coordinates associated

with each feature.

intensities (cell value) A numeric vector containing the intensity associated with each fea-

ture.

stdvs (cell value) A numeric vector containing the standard deviation associated with

each feature.

pixels (cell value) An integer vector containing the number of pixels associated with

each feature.

outliers An integer vector of indices specifying which of the queried cells that are

flagged as outliers. Note that there is a difference between outliers=NULL and outliers=integer(0); the last case happens when readOutliers=TRUE but

there are no outliers.

readCel 61

masked

An integer vector of indices specifying which of the queried cells that are flagged as masked. Note that there is a difference between masked=NULL and masked=integer(0); the last case happens when readMasked=TRUE but there are no masked features.

The elements of the cell values are ordered according to argument indices. The lengths of the cell-value elements equals the number of cells read.

Which of the above elements that are returned are controlled by the readNnn arguments. If FALSE, the corresponding element above is NULL, e.g. if readStdvs=FALSE then stdvs is NULL.

Outliers and masked cells

The Affymetrix image analysis software flags cells as outliers and masked. This method does not return these flags, but instead vectors of cell indices listing which cells of the queried cells are outliers and masked, respectively. The current community view seems to be that this should be done based on statistical modeling of the actual probe intensities and should be based on the choice of preprocessing algorithm. Most algorithms are only using the intensities from the CEL file.

Memory usage

The Fusion SDK allocates memory for the entire CEL file, when the file is accessed (but does not actually read the file into memory). Using the indices argument will therefore only affect the memory use of the final object (as well as speed), not the memory allocated in the C function used to parse the file. This should be a minor problem however.

Troubleshooting

It is considered a bug if the file contains information not accessible by this function, please report it.

Author(s)

James Bullard and Kasper Daniel Hansen

See Also

readCelHeader() for a description of the header output. Often a user only wants to read the intensities, look at readCelIntensities() for a function specialized for that use.

```
for (zzz in 0) { # Only so that 'break' can be used
# Scan current directory for CEL files
celFiles <- list.files(pattern="[.](c|C)(e|E)(1|L)$")
if (length(celFiles) == 0)
    break;
celFile <- celFiles[1]
# Read a subset of cells</pre>
```

62 readCelHeader

```
idxs <- c(1:5, 1250:1500, 450:440)
cel <- readCel(celFile, indices=idxs, readOutliers=TRUE)
str(cel)
# Clean up
rm(celFiles, celFile, cel)
} # for (zzz in 0)</pre>
```

readCelHeader

Parsing the header of an Affymetrix CEL file

Description

Reads in the header of an Affymetrix CEL file using the Fusion SDK.

Usage

```
readCelHeader(filename)
```

Arguments

filename the name of the CEL file.

Details

This function returns the header of a CEL file. Affymetrix operates with different versions of this file format. Depending on what version is being read, different information is accessible.

Value

A named list with components described below. The entries are obtained from the Fusion SDK interface functions. We try to obtain all relevant information from the file.

filename the name of the cel file.
version the version of the cel file.

cols the number of columns on the chip.
rows the number of rows on the chip.

total the total number of features on the chip. Usually equal to rows times cols, but

since it is a separate attribute in the SDK we decided to include it anyway.

algorithm the algorithm used to create the CEL file.

parameters the parameters used in the algorithm. Seems to be semi-colon separated.

chiptype the type of the chip.

header the entire header of the CEL file. Only available for non-calvin format files.

datheader the entire dat header of the CEL file. This contains for example a date.

readCelIntensities 63

librarypackage the library package name of the file. Empty for older versions.

cellmargin a parameter used to generate the CEL file. According to Affymetrix, it desig-

nates the number of pixels to ignore around the feature border when calculating

the intensity value (the number of pixels ignored are cellmargin divided by 2).

noutliers the number of features reported as outliers.

nmasked the number of features reported as masked.

Note

Memory usage: the Fusion SDK allocates memory for the entire CEL file, when the file is accessed. The memory footprint of this function will therefore seem to be (rather) large.

Speed: CEL files of version 2 (standard text files) needs to be completely read in order to report the number of outliers and masked features.

Author(s)

James Bullard and Kasper Daniel Hansen

See Also

readCel() for reading in the entire CEL file. That function also returns the header. See affxparserInfo for general comments on the package and the Fusion SDK.

Examples

```
# Scan current directory for CEL files
files <- list.files(pattern="[.](c|C)(e|E)(1|L)$")
if (length(files) > 0) {
   header <- readCelHeader(files[1])
   print(header)
   rm(header)
}
# Clean up
rm(files)</pre>
```

readCelIntensities

Reads the intensities contained in several Affymetrix CEL files

Description

Reads the intensities of several Affymetrix CEL files (as opposed to readCe1() which only reads a single file).

Usage

```
readCelIntensities(filenames, indices = NULL, ..., verbose = 0)
```

64 readCelIntensities

Arguments

indices a vector of which indices should be read. If the argument is NULL all features

will be returned.

... Additional arguments passed to readCel().

verbose an integer: how verbose do we want to be, higher means more verbose.

Details

The function will initially allocate a matrix with the same memory footprint as the final object.

Value

A matrix with a number of rows equal to the length of the indices argument (or the number of features on the entire chip), and a number of columns equal to the number of files. The columns are ordered according to the filenames argument.

Note

Currently this function builds on readCel(), and simply calls this function multiple times. If testing yields sufficient reasons for doing so, it may be re-implemented in C++.

Author(s)

James Bullard and Kasper Daniel Hansen

See Also

readCel() for a discussion of a more versatile function, particular with details of the indices argument.

```
# Scan current directory for CEL files
files <- list.files(pattern="[.](c|C)(e|E)(1|L)$")
if (length(files) >= 2) {
  cel <- readCelIntensities(files[1:2])
  str(cel)
  rm(cel)
}
# Clean up
rm(files)</pre>
```

readCelRectangle 65

readCelRectangle	Reads a spatial subset of probe-level data from Affymetrix CEL files

Description

Reads a spatial subset of probe-level data from Affymetrix CEL files.

Usage

```
readCelRectangle(filename, xrange=c(0, Inf), yrange=c(0, Inf), ..., asMatrix=TRUE)
```

Arguments

filename	The pathname of the CEL file.
xrange	A numeric vector of length two giving the left and right coordinates of the cells to be returned.
yrange	A numeric vector of length two giving the top and bottom coordinates of the cells to be returned.
	Additional arguments passed to readCel().
asMatrix	If TRUE, the CEL data fields are returned as matrices with element (1,1) corresponding to cell (xrange[1],yrange[1]).

Value

A named list CEL structure similar to what readCel(). In addition, if asMatrix is TRUE, the CEL data fields are returned as matrices, otherwise not.

Author(s)

Henrik Bengtsson

See Also

The readCel() method is used internally.

66 readCelUnits

readCelUnits

Reads probe-level data ordered as units (probesets) from one or several Affymetrix CEL files

Description

Reads probe-level data ordered as units (probesets) from one or several Affymetrix CEL files by using the unit and group definitions in the corresponding Affymetrix CDF file.

Usage

```
readCelUnits(filenames, units=NULL, stratifyBy=c("nothing", "pmmm", "pm", "mm"),
  cdf=NULL, ..., addDimnames=FALSE, dropArrayDim=TRUE, transforms=NULL, readMap=NULL,
  verbose=FALSE)
```

Arguments

filenames	The filenames of the CEL files.
units	An integer vector of unit indices specifying which units to be read. If NULL, all units are read.
stratifyBy	Argument passed to low-level method readCdfCellIndices.
cdf	A character filename of a CDF file, or a CDF list structure. If NULL, the

CDF file is searched for by findCdf() first starting from the current directory

and then from the directory where the first CEL file is.

readCelUnits 67

... Arguments passed to low-level method readCel, e.g. readXY and readStdvs.

returned CEL structure in bytes increases by 30-40% with dimension names.

dropArrayDim If TRUE and only one array is read, the elements of the group field do not have

an array dimension.

transforms A list of exactly length(filenames) functions. If NULL, no transformation

is performed. Intensities read are passed through the corresponding transform

function before being returned.

readMap A vector remapping cell indices to file indices. If NULL, no mapping is used.

verbose Either a logical, a numeric, or a Verbose object specifying how much ver-

bose/debug information is written to standard output. If a Verbose object, how detailed the information is is specified by the threshold level of the object. If a numeric, the value is used to set the threshold of a new Verbose object. If TRUE, the threshold is set to -1 (minimal). If FALSE, no output is written (and neither

is the **R.utils** package required).

Value

A named list with one element for each unit read. The names corresponds to the names of the units read. Each unit element is in turn a list structure with groups (aka blocks). Each group contains requested fields, e.g. intensities, stdvs, and pixels. If more than one CEL file is read, an extra dimension is added to each of the fields corresponding, which can be used to subset by CEL file.

Note that neither CEL headers nor information about outliers and masked cells are returned. To access these, use readCelHeader() and readCel().

Author(s)

Henrik Bengtsson

References

[1] Affymetrix Inc, Affymetrix GCOS 1.x compatible file formats, June 14, 2005. http://www.affymetrix.com/support/developer/

See Also

Internally, readCelHeader(), readCdfUnits() and readCel() are used.

68 readChp

readChp

A function to read Affymetrix CHP files

Description

This function will parse any type of CHP file and return the results in a list. The contents of the list will depend on the type of CHP file that is parsed and readers are referred to Affymetrix documentation of what should be there, and how to interpret it.

Usage

```
readChp(filename, withQuant = TRUE)
```

Arguments

filename The name of the CHP file to read.

withQuant A boolean value, currently largely unused.

Details

This is an interface to the Affymetrix Fusion SDK. The Affymetrix documentation should be consulted for explicit details.

Value

A list is returned. The contents of the list depend on the type of CHP file that was read. Users may want to translate the different outputs into specific containers.

Troubleshooting

It is considered a bug if the file contains information not accessible by this function, please report it.

Author(s)

R. Gentleman

readClf 69

See Also

```
readCel
```

Examples

readClf

Parsing a CLF file using Affymetrix Fusion SDK

Description

This function parses a CLF file using the Affymetrix Fusion SDK. CLF (chip layout) files contain information associating probe ids with chip x- and y- coordinates.

Usage

```
readClf(file)
```

Arguments

file

character(1) providing a path to the CLF file to be input.

Value

An list. The header element is always present.

header	A list with information about the CLF file. The list contains elements described in the CLF file format document referenced below.
dims	A length-two integer vector of chip x- and y-coordinates.
id	An integer vector of length prod(dims) containing probe identifiers.
x	An integer vector of length prod(dims) containing x-coordinates corresponding to the entries in id.
У	An integer vector of length prod(dims) containing y-coordinates corresponding to the entries in id.

Author(s)

Martin Morgan

70 readClfEnv

See Also

https://www.affymetrix.com/support/developer/fusion/File_Format_CLF_aptv161.pdf describes CLF file content.

readClfEnv

Parsing a CLF file using Affymetrix Fusion SDK

Description

This function parses a CLF file using the Affymetrix Fusion SDK. CLF (chip layout) files contain information associating probe ids with chip x- and y- coordinates.

Usage

```
readClfEnv(file, readBody = TRUE)
```

Arguments

file character(1) providing a path to the CLF file to be input.

readBody logical(1) indicating whether the entire file should be parsed (TRUE) or only

the file header information describing the chips to which the file is relevant.

Value

An environment. The header element is always present; the remainder are present when readBody=TRUE.

header A list with information about the CLF file. The list contains elements described

in the CLF file format document referenced below.

dims A length-two integer vector of chip x- and y-coordinates.

id An integer vector of length prod(dims) containing probe identifiers.

x An integer vector of length prod(dims) containing x-coordinates corresponding

to the entries in id.

y An integer vector of length prod(dims) containing y-coordinates corresponding

to the entries in id.

Author(s)

Martin Morgan

See Also

https://www.affymetrix.com/support/developer/fusion/File_Format_CLF_aptv161.pdf describes CLF file content.

readClfHeader 71

readClfHeader	Read the header of a CLF file.

Description

Reads the header of a CLF file. The exact information stored in this file can be viewed in the readClfEnv() documentation which reads the header in addition to the body.

Usage

```
readClfHeader(file)
```

Arguments

file file a CLF file

Value

A list of header elements.

readPgf Parsing a PGF file using Affymetrix Fusion SDK

Description

This function parses a PGF file using the Affymetrix Fusion SDK. PGF (probe group) files describe probes present within probe sets, including the type (e.g., pm, mm) of the probe and probeset.

Usage

```
readPgf(file, indices = NULL)
```

Arguments

file character(1) providing a path to the PGF file to be input. indices integer(n) a vector of indices of the probesets to be read.

Value

An list. The header element is always present; the remainder are present when readBody=TRUE.

The elements present when readBody=TRUE describe probe sets, atoms, and probes. Elements within probe sets, for instance, are coordinated such that the ith index of one vector (e.g., probesetId) corresponds to the ith index of a second vector (e.g., probesetType). The atoms contained within probeset i are in positions probesetStartAtom[i]:(probesetStartAtom[i+1]-1) of the atom vectors. A similar map applies to probes within atoms, using atomStartProbe as the index.

The PGF file format includes optional elements; these elements are always present in the list, but with appropriate default values.

72 readPgfEnv

header A list with information about the PGF file. The list contains elements described

in the PGF file format document referenced below.

probesetId integer vector of probeset identifiers.

probesetType character vector of probeset types. Types are described in the PGF file format

document.

probesetName character vector of probeset names.

probesetStartAtom

integer vector of the start index (e.g., in the element atomId of atoms belonging

to this probeset).

atomId integer vector of atom identifiers.

atomExonPosition

integer vector of probe interrogation position relative to the target sequence.

atomStartProbe integer vector of the start index (e.g., in the element probeId of probes belong-

ing to this atom).

probeId integer vector of probe identifiers.

probeType character vector of probe types. Types are described in the PGF file format

document.

probeGcCount integer vector of probe GC content.

probeLength integer vector of probe lengths.

probeInterrogationPosition

integer vector of the position, within the probe, at which interrogation occurs.

probeSequence character vector of the probe sequence.

Author(s)

Martin Morgan

See Also

https://www.affymetrix.com/support/developer/fusion/File_Format_PGF_aptv161.pdf describes PGF file content.

The internal function .pgfProbeIndexFromProbesetIndex provides a map between the indices of probe set entries and the indices of the probes contained in the probe set.

readPgfEnv Parsing a PGF file using Affymetrix Fusion SDK

Description

This function parses a PGF file using the Affymetrix Fusion SDK. PGF (probe group) files describe probes present within probe sets, including the type (e.g., pm, mm) of the probe and probeset.

readPgfEnv 73

Usage

```
readPgfEnv(file, readBody = TRUE, indices = NULL)
```

Arguments

file character(1) providing a path to the PGF file to be input.

readBody logical(1) indicating whether the entire file should be parsed (TRUE) or only

the file header information describing the chips to which the file is relevant.

indices integer(n) vector of positive integers indicating which probesets to read. These

integers must be sorted (increasing) and unique.

Value

An environment. The header element is always present; the remainder are present when readBody=TRUE.

The elements present when readBody=TRUE describe probe sets, atoms, and probes. Elements within probe sets, for instance, are coordinated such that the ith index of one vector (e.g., probesetId) corresponds to the ith index of a second vector (e.g., probesetType). The atoms contained within probeset i are in positions probesetStartAtom[i]:(probesetStartAtom[i+1]-1) of the atom vectors. A similar map applies to probes within atoms, using atomStartProbe as the index.

The PGF file format includes optional elements; these elements are always present in the environment, but with appropriate default values.

header A list with information about the PGF file. The list contains elements described

in the PGF file format document referenced below.

probesetId integer vector of probeset identifiers.

probesetType character vector of probeset types. Types are described in the PGF file format

document.

probesetName character vector of probeset names.

probesetStartAtom

integer vector of the start index (e.g., in the element atomId of atoms belonging

to this probeset).

atomId integer vector of atom identifiers.

atomExonPosition

integer vector of probe interrogation position relative to the target sequence.

atomStartProbe integer vector of the start index (e.g., in the element probeId of probes belong-

ing to this atom).

probeId integer vector of probe identifiers.

probeType character vector of probe types. Types are described in the PGF file format

document.

probeGcCount integer vector of probe GC content.

probeLength integer vector of probe lengths.

probeInterrogationPosition

integer vector of the position, within the probe, at which interrogation occurs.

probeSequence character vector of the probe sequence.

74 readPgfHeader

Author(s)

Martin Morgan

See Also

https://www.affymetrix.com/support/developer/fusion/File_Format_PGF_aptv161.pdf describes PGF file content.

The internal function .pgfProbeIndexFromProbesetIndex provides a map between the indices of probe set entries and the indices of the probes contained in the probe set.

readPgfHeader

Read the header of a PGF file into a list.

Description

This function reads the header of a PGF file into a list more details on what the exact fields are can be found in the details section.

Usage

```
readPgfHeader(file)
```

Arguments

file

file: A file in PGF format

Details

https://www.affymetrix.com/support/developer/fusion/File_Format_PGF_aptv161.pdf

Value

A list corresponding to the elements in the header.

updateCel 75

|--|

Description

Updates a CEL file.

Usage

```
updateCel(filename, indices=NULL, intensities=NULL, stdvs=NULL, pixels=NULL,
writeMap=NULL, ..., verbose=0)
```

Arguments

filename	The filename of the CEL file.
indices	A numeric vector of cell (probe) indices specifying which cells to updated. If NULL, all indices are considered.
intensities	A numeric vector of intensity values to be stored. Alternatively, it can also be a named data. frame or matrix (or list) where the named columns (elements) are the fields to be updated.
stdvs	A optional numeric vector.
pixels	A optional numeric vector.
writeMap	An optional write map.
	Not used.
verbose	An integer specifying how much verbose details are outputted.

Details

Currently only binary (v4) CEL files are supported. The current version of the method does not make use of the Fusion SDK, but its own code to navigate and update the CEL file.

Value

Returns (invisibly) the pathname of the file updated.

Author(s)

Henrik Bengtsson

76 updateCel

Examples

```
if (require("AffymetrixDataTestFiles")) {
# Search for some available Calvin CEL files
path <- system.file("rawData", package="AffymetrixDataTestFiles")</pre>
files <- findFiles(pattern="[.](cel|CEL)$", path=path, recursive=TRUE, firstOnly=FALSE)\\
files <- grep("FusionSDK_HG-U133A", files, value=TRUE)</pre>
files <- grep("Calvin", files, value=TRUE)</pre>
file <- files[1]
# Convert to an XDA CEL file
filename <- file.path(tempdir(), basename(file))</pre>
if (file.exists(filename))
 file.remove(filename)
convertCel(file, filename)
fields <- c("intensities", "stdvs", "pixels")</pre>
# Cells to be updated
idxs <- 1:2
# Get CEL header
hdr <- readCelHeader(filename)</pre>
# Get the original data
cel <- readCel(filename, indices=idxs, readStdvs=TRUE, readPixels=TRUE)</pre>
print(cel[fields])
cel0 <- cel
# Square-root the intensities
updateCel(filename, indices=idxs, intensities=sqrt(cel$intensities))
cel <- readCel(filename, indices=idxs, readStdvs=TRUE, readPixels=TRUE)</pre>
print(cel[fields])
# Update a few cell values by a data frame
data <- data.frame(</pre>
 intensities=cel0$intensities,
 stdvs=c(201.1, 3086.1)+0.5,
 pixels=c(9,9+1)
updateCel(filename, indices=idxs, data)
# Assert correctness of update
cel <- readCel(filename, indices=idxs, readStdvs=TRUE, readPixels=TRUE)</pre>
```

updateCel 77

```
print(cel[fields])
for (ff in fields) {
 stopifnot(all.equal(cel[[ff]], data[[ff]], .Machine$double.eps^0.25))
# Update a region of the CEL file
# Load pre-defined data
side <- 306
pathname <- system.file("extras/easternEgg.gz", package="affxparser")</pre>
con <- gzfile(pathname, open="rb")</pre>
z <- readBin(con=con, what="integer", size=1, signed=FALSE, n=side^2)</pre>
close(con)
z <- matrix(z, nrow=side)</pre>
side <- min(hdr$cols - 2*22, side)</pre>
z <- as.double(z[1:side,1:side])</pre>
x <- matrix(22+0:(side-1), nrow=side, ncol=side, byrow=TRUE)</pre>
idxs \leftarrow as.vector((1 + x) + hdr$cols*t(x))
# Load current data in the same region
z0 <- readCel(filename, indices=idxs)$intensities</pre>
# Mix the two data sets
z < -(0.3*z^2 + 0.7*z0)
# Update the CEL file
updateCel(filename, indices=idxs, intensities=z)
# Make some spatial changes
rotate270 \leftarrow function(x, ...) {
 x \leftarrow t(x)
 nc <- ncol(x)
 if (nc < 2) return(x)
 x[,nc:1,drop=FALSE]
# Display a spatial image of the updated CEL file
cel <- readCelRectangle(filename, xrange=c(0,350), yrange=c(0,350))</pre>
z <- rotate270(cel$intensities)</pre>
sub <- paste("Chip type:", cel$header$chiptype)</pre>
image(z, col=gray.colors(256), axes=FALSE, main=basename(filename), sub=sub)
text(x=0, y=1, labels="(0,0)", adj=c(0,-0.7), cex=0.8, xpd=TRUE)
text(x=1, y=0, labels="(350,350)", adj=c(1,1.2), cex=0.8, xpd=TRUE)
# Clean up
file.remove(filename)
rm(files, cel, cel0, idxs, data, ff, fields, rotate270)
# STOP #
```

78 updateCelUnits

updateCelUnits Updates a CEL file unit by unit
--

Description

Updates a CEL file unit by unit.

Please note that, contrary to readCelUnits(), this method can only update a single CEL file at the time.

Usage

```
updateCelUnits(filename, cdf=NULL, data, ..., verbose=0)
```

Arguments

filename	The filename of the CEL file.
cdf	A (optional) CDF list structure either with field indices or fields x and y. If NULL, the unit names (and from there the cell indices) are inferred from the names of the elements in data.
data	A list structure in a format similar to what is returned by readCelUnits() for a single CEL file only.
•••	Optional arguments passed to readCdfCellIndices(), which is called if cdf is not given.
verbose	An integer specifying how much verbose details are outputted.

Value

Returns what updateCel() returns.

Working with re-arranged CDF structures

Note that if the cdf structure is specified the CDF file is *not* queried, but all information about cell x and y locations, that is, cell indices is expected to be in this structure. This can be very useful when one work with a cdf structure that originates from the underlying CDF file, but has been restructured for instance through the applyCdfGroups() method, and data correspondingly. This update method knows how to update such structures too.

Author(s)

Henrik Bengtsson

See Also

Internally, updateCel() is used.

updateCelUnits 79

Examples

```
if (require("AffymetrixDataTestFiles")) {
# Search for some available Calvin CEL files
path <- system.file("rawData", package="AffymetrixDataTestFiles")</pre>
files <- findFiles(pattern="[.](cel|CEL)$", path=path, recursive=TRUE, firstOnly=FALSE)
files <- grep("FusionSDK_Test3", files, value=TRUE)</pre>
files <- grep("Calvin", files, value=TRUE)</pre>
file <- files[1]
# Convert to an XDA CEL file
pathname <- file.path(tempdir(), basename(file))</pre>
if (file.exists(pathname))
 file.remove(pathname)
convertCel(file, pathname)
# Check for the CDF file
hdr <- readCelHeader(pathname)</pre>
cdfFile <- findCdf(hdr$chiptype)</pre>
hdr <- readCdfHeader(cdfFile)</pre>
nbrOfUnits <- hdr$nunits</pre>
print(nbr0fUnits);
# Example: Read and re-write the same data
units <- c(101, 51)
data1 <- readCelUnits(pathname, units=units, readStdvs=TRUE)</pre>
cat("Original data:\n")
str(data1)
updateCelUnits(pathname, data=data1)
data2 <- readCelUnits(pathname, units=units, readStdvs=TRUE)</pre>
cat("Updated data:\n")
str(data2)
stopifnot(identical(data1, data2))
# Example: Random read and re-write "stress test"
for (kk in 1:10) {
 nunits <- sample(min(1000,nbrOfUnits), size=1)</pre>
 units <- sample(nbr0fUnits, size=nunits)</pre>
 cat(sprintf("%02d. Selected %d random units: reading", kk, nunits));
 t <- system.time({</pre>
   data1 <- readCelUnits(pathname, units=units, readStdvs=TRUE)</pre>
```

80 writeCdf

writeCdf

}, gcFirst=TRUE)[3]

Creates a binary CDF file

Description

This function creates a binary CDF file given a valid CDF structure containing all necessary elements.

Warning: The API for this function is likely to be changed in future versions.

Usage

```
writeCdf(fname, cdfheader, cdf, cdfqc, overwrite=FALSE, verbose=0)
```

Arguments

fname name of the CDF file.

cdfheader A list with a structure equal to the output of readCdfHeader.

cdf A list with a structure equal to the output of readCdf.
cdfqc A list with a structure equal to the output of readCdfQc.

overwrite Overwrite existing file?

verbose how verbose should the output be. 0 means no output, with higher numbers

being more verbose.

Details

This function has been validated mainly by reading in various ASCII or binary CDF files which are written back as new CDF files, and compared element by element with the original files.

Value

This function is used for its byproduct: creating a CDF file.

writeCdfHeader 81

Author(s)

Kasper Daniel Hansen

See Also

To read the CDF "regular" and QC units with all necessary fields and values for writing a CDF file, see readCdf, readCdfQc() and readCdfHeader. To compare two CDF files, see compareCdfs.

writeCdfHeader Writes a CDF header

Description

Writes a CDF header. This method is not intended to be used explicitly. To write a CDF, use writeCdf() instead.

Usage

writeCdfHeader(con, cdfHeader, unitNames, qcUnitLengths, unitLengths, verbose=0)

Arguments

con An open connection to which nothing has been written.

cdfHeader A CDF header list structure.

unitNames A character vector of all unit names.

qcUnitLengths An integer vector of all the number of bytes in each of the QC units.

unitLengths An integer vector of all the number of bytes in each of the (ordinary) units.

verbose An integer specifying how much verbose details are outputted.

Value

Returns nothing.

Author(s)

Henrik Bengtsson

See Also

This method is called by writeCdf(). See also writeCdfQcUnits() and writeCdfUnits().

82 writeCdfUnits

writeCdfQcUnits Writes CDF QC units

Description

Writes CDF QC units. This method is not intended to be used explicitly. To write a CDF, use writeCdf() instead.

Usage

```
writeCdfQcUnits(con, cdfQcUnits, verbose=0)
```

Arguments

con An open connection to which a CDF header already has been written by writeCdfHeader().

cdfQcUnits A list structure of CDF QC units as returned by readCdf() (not readCdfUnits()).

verbose An integer specifying how much verbose details are outputted.

Value

Returns nothing.

Author(s)

Henrik Bengtsson

See Also

This method is called by writeCdf(). See also writeCdfHeader() and writeCdfUnits().

writeCdfUnits Writes CDF units

Description

Writes CDF units. *This method is not intended to be used explicitly. To write a CDF, use* writeCdf() *instead.*

Usage

```
writeCdfUnits(con, cdfUnits, verbose=0)
```

writeCelHeader 83

Arguments

con An open connection to which a CDF header and QC units already have been

written by writeCdfHeader() and writeCdfQcUnits(), respectively.

cdfUnits A list structure of CDF units as returned by readCdf() (not readCdfUnits()).

verbose An integer specifying how much verbose details are outputted.

Value

Returns nothing.

Author(s)

Henrik Bengtsson

See Also

This method is called by writeCdf(). See also writeCdfHeader() and writeCdfQcUnits().

writeCelHeader

Writes a CEL header to a connection

Description

Writes a CEL header to a connection.

Usage

```
writeCelHeader(con, header, outputVersion=c("4"), ...)
```

Arguments

con A connection.

header A list structure describing the CEL header, similar to the structure returned by

readCelHeader().

outputVersion A character string specifying the output format. Currently only CEL version

4 (binary;XDA) are supported.

... Not used.

Details

Currently only CEL version 4 (binary; XDA) headers can be written.

Value

Returns (invisibly) the pathname of the file created.

84 writeTpmap

Redundant fields

The CEL v4 header contains redundant information. To avoid inconsistency this method generates such redundant values from the original values. This is consistent to how the CEL reader in Fusion SDK does it, cf. readCelHeader(). The redundant information is in the (CEL v3) header field, which contains the CEL header information as it would appear in the CEL v3 format. This in turn contains a DAT header field reproducing the DAT header from the image analysis. It is from this DAT header that the chip type is extracted.

Author(s)

Henrik Bengtsson

writeTpmap

Writes BPMAP and TPMAP files

Description

Writes BPMAP and TPMAP files.

Usage

```
writeTpmap(filename, bpmaplist, verbose = 0)
tpmap2bpmap(tpmapname, bpmapname, verbose = 0)
```

Arguments

filename The filename.

bpmaplist A list structure similar to the result of readBpmap.

tpmapname Filename of the TPMAP file.

bpmapname Filename of the BPMAP file.

verbose How verbose do we want to be.

Details

writeTpmap writes a text probe map file, while tpmap2bpmap converts such a file to a binary probe mapping file. Somehow Affymetrix has different names for the same structure, depending on whether the file is binary or text. I have seen many TPMAP files referred to as BPMAP files.

Value

These functions are called for their side effects (creating files).

Author(s)

Kasper Daniel Hansen

writeTpmap 85

See Also

readBpmap

Index

* IO	updateCelUnits, 78
compareCdfs, 26	writeCdf, 80
compareCels, 27	writeCdfHeader, 81
convertCdf, 28	writeCdfQcUnits,82
convertCel, 29	writeCdfUnits,82
copyCel, 31	writeCelHeader,83
createCel, 32	writeTpmap,84
findCdf, 33	* documentation
findFiles, 35	1. Dictionary, 5
invertMap, 36	2. Cell coordinates and cell
isCelFile, 37	indices, 6
parseDatHeaderString, 38	9. Advanced - Cell-index maps for
readBpmap, 39	reading and writing, 8
readCcg, 40	* file
readCcgHeader, 42	compareCdfs, 26
readCdf, 43	compareCels, 27
readCdfCellIndices, 45	convertCdf, 28
readCdfDataFrame, 46	convertCel, 29
readCdfGroupNames, 47	copyCel, 31
readCdfHeader, 48	createCel, 32
readCdfIsPm,49	findCdf, 33
<pre>readCdfNbrOfCellsPerUnitGroup, 50</pre>	findFiles, 35
readCdfQc, 52	invertMap, 36
readCdfUnitNames, 53	isCelFile, 37
readCdfUnits, 54	parseDatHeaderString,38
readCdfUnitsWriteMap, 56	readBpmap, 39
readCel, 59	readCcg, 40
readCelHeader, 62	readCcgHeader, 42
readCelIntensities, 63	readCdf, 43
readCelRectangle, 65	readCdfCellIndices, 45
readCelUnits,66	readCdfDataFrame,46
readChp, 68	readCdfGroupNames, 47
readClf, 69	readCdfHeader,48
readClfEnv, 70	readCdfIsPm, 49
readClfHeader,71	readCdfNbrOfCellsPerUnitGroup, 50
readPgf, 71	readCdfQc, 52
readPgfEnv, 72	readCdfUnitNames, 53
readPgfHeader, 74	readCdfUnits,54
updateCel, 75	readCdfUnitsWriteMap, 56

INDEX 87

readCel,59	writeCdfHeader, 81
readCelHeader, 62	writeCdfQcUnits,82
readCelIntensities, 63	writeCdfUnits,82
readCelRectangle, 65	writeCelHeader, 83
readCelUnits,66	* package
readChp, 68	affxparser-package, 3
readClf, 69	* programming
readClfEnv, 70	applyCdfGroupFields, 10
readPgf, 71	applyCdfGroups, 11
readPgfEnv, 72	arrangeCelFilesByChipType, 14
readPgfHeader, 74	cdfAddBaseMmCounts, 15
updateCel, 75	cdfAddPlasqTypes, 16
updateCelUnits, 78	cdfAddProbeOffsets, 17
writeCdf, 80	cdfGetFields, 18
writeCdfHeader, 81	cdfGetGroups, 19
writeCdfQcUnits, 82	cdfGtypeCelToPQ, 19
writeCdfUnits,82	cdfHeaderToCelHeader, 20
writeCelHeader, 83	cdfMergeAlleles, 21
writeTpmap, 84	cdfMergeStrands, 22
internal	cdfMergeToQuartets, 23
9. Advanced - Cell-index maps for	cdf0rderBy, 24
reading and writing, 8	cdfOrderColumnsBy, 24
arrangeCelFilesByChipType, 14	cdfSetDimension, 25
cdfAddBaseMmCounts, 15	copyCel, 31
cdfAddPlasqTypes, 16	isCelFile, 37
cdfAddProbeOffsets, 17	parseDatHeaderString,38
cdfGetFields, 18	1. Dictionary, 5
cdfGetGroups, 19	2. Cell coordinates and cell indices, 6
cdfGtypeCelToPQ, 19	9. Advanced - Cell-index maps for
cdfHeaderToCelHeader, 20	reading and writing, 8
cdfMergeAlleles, 21	
cdfMergeStrands, 22	affxparser (affxparser-package), 3
cdfMergeToQuartets, 23	affxparser-package, 3
_	<pre>applyCdfBlocks (applyCdfGroups), 11</pre>
cdfOrderBy, 24	applyCdfGroupFields, 10, 25, 26
cdfOrderColumnsBy, 24	applyCdfGroups, <i>3</i> , <i>11</i> , 11, <i>15</i> – <i>25</i> , <i>78</i>
cdfSetDimension, 25	arrangeCelFilesByChipType, 14
copyCel, 31	
findFiles, 35	cdfAddBaseMmCounts, 12, 15
invertMap, 36	cdfAddPlasqTypes, 16
isCelFile, 37	cdfAddProbeOffsets, <i>12</i> , <i>15</i> , <i>16</i> , 17
parseDatHeaderString, 38	cdfGetFields, <i>12</i> , 18
readCdf, 43	cdfGetGroups, 12, 19
readCdfDataFrame, 46	cdfGtypeCelToPQ, 12, 19
readCdfIsPm, 49	cdfHeaderToCelHeader, 20
readCdfNbrOfCellsPerUnitGroup, 50	cdfMergeAlleles, 12, 21
readCdfQc, 52	cdfMergeStrands, 12, 22
readCdfUnitsWriteMap, 56	cdfMergeToQuartets, 23
writeCdf, 80	cdf0rderBy, <i>12</i> , <i>24</i> , <i>25</i>

88 INDEX

cdfOrderColumnsBy, 12, 24, 24	readCdfGroupNames, 47
cdfSetDimension, 25	readCdfHeader, <i>44</i> , 48, <i>81</i>
character, 14, 18, 19, 21, 34, 35, 38, 44–46,	readCdfIsPm,49
48, 53, 54, 66, 81, 83	<pre>readCdfNbrOfCellsPerUnitGroup, 50</pre>
compareCdfs, 26, 29, 81	readCdfQc, 52, 81
compareCels, 27	readCdfUnitNames, 53
connection, $81-83$	readCdfUnits, 3, 41, 43, 44, 46–49, 53, 54,
convertCdf, 26, 28	56, 67, 82, 83
	readCdfUnitsWriteMap, 9, 37, 56
convertCel, 27, 29	
copyCel, 31	readCe1, 3, 38, 57, 59, 63–65, 67, 69
createCel, 30, 32	readCelHeader, 14, 32, 38, 61, 62, 67, 83, 84
1	readCelIntensities, 61, 63
data.frame, 46, 47, 75	readCelRectangle, 65
EN CE 26 20 20 22 25 20 56 67	readCelUnits, 3, 8, 9, 34, 38, 57, 66, 78
FALSE, 26–28, 30–32, 35, 38, 56, 67	readChp, 68
filePath, 36	readClf, 69
findCdf, 3, 4, 33, 66	readClfEnv, 70, 71
findFiles, <i>34</i> , <i>35</i>	readClfHeader, 71
function, <i>11</i> , <i>15–19</i> , <i>21–25</i> , <i>67</i>	readPgf, 71
	readPgfEnv, 72
integer, 9, 17, 19, 25–27, 32, 36, 40, 42–46,	readPgfHeader, 74
48, 50, 53–56, 66, 75, 78, 81–83	
invertMap, 9, 36, 57	Startup, <i>34</i>
isCelFile, <i>31</i> , 37	strptime, 38
	00. p 010, 0 0
list, 11, 15–25, 32, 38, 40, 42, 45, 47, 48, 50,	tpmap2bpmap, 40
54, 55, 65–67, 75, 78, 81–83	tpmap2bpmap(writeTpmap), 84
list.files, 35	TRUE, 9, 21, 26–28, 30–32, 34, 35, 38, 43, 46,
logical, 48, 50, 56, 67	54, 56, 65, 67
	31, 30, 03, 07
matrix, 21, 75	updateCel, 75, 78
	updateCelUnits, 78
NA, <i>18</i>	apadececionics, 76
NULL, 34, 38, 43, 45, 46, 48, 50, 53, 54, 56, 60,	vector, 9, 14, 17–19, 34–36, 43–48, 50, 53,
66, 67, 75, 78	54, 56, 60, 65–67, 75, 81
numeric, 56, 65, 67, 75	Verbose, 56, 67
1.4	ver bose, 50, 07
order, 24, 25, 36	writeCdf, 29, 80, 81-83
2. 22., 2., 2., 2.	writeCdfHeader, 81, 82, 83
parseDatHeaderString, 38	writeCdfQcUnits, <i>81</i> , 82, <i>83</i>
3) - 1	
readBin, 56	writeCdfUnits, 81, 82, 82
readBpmap, 39, 85	writeCelHeader, 83
readBpmapHeader (readBpmap), 39	writeTpmap,84
readBpmapSeqinfo(readBpmap), 39	
readCcg, 40, 42	xy2indices, 7
readCcgHeader, 41, 42	
readCdf, 3, 43, 53, 81–83	
readCdfCellIndices, 3, 45, 55, 66, 78	
readCdfDataFrame, 46	