

Package ‘msdata’

September 23, 2012

Version 0.1.10

Date 2012-04-25

Title Various Mass Spectrometry raw data example files

Author Steffen Neumann <sneumann@ipb-halle.de>

Maintainer Steffen Neumann <sneumann@ipb-halle.de>

Suggests xcms

ZipData no

Description Ion Trap positive ionization mode data in mzData file format.

Subset from 500-850 m/z and 1190-1310 seconds, incl. MS2 and MS3,intensity threshold 100.000. Extracts from FTICR Apex III, m/z 400-450.

Subset of UPLC - Bruker micrOTOFq data, both mzData and mzML.

biocViews ExperimentData

License GPL (>= 2)

R topics documented:

msdata *Sample FTICR, LC/MS and MS\$^n\$ data*

Description

x object containing a subset of LC/MS raw data from a Thermo Finnigan LCQ Deca XP. The data is a subset from 500-850 m/z and 1190-1310 seconds, incl. MS2 and MS3, intensity threshold 100.000. It was collected in positive ionization mode.

`xs` object containing a subset of FTICR data from a Bruker APEx III FTICR. The data is a subset from 400-450 m/z, collected in positive ionization mode.

Usage

```
data(xs)
```

Format

The format is:

```
xs
```

Details

The corresponding raw mzdata files are located in the `fticr` and `iontrap` subdirectory of this package.

See Also

[xcmsSet](#), [xcmsRaw](#)

Examples

```
## The directory with the mzData LC/MS files
data(xs)
mzdatapath <- file.path(.find.package("msdata"), "iontrap")
mzdatapath
files <- list.files(mzdatapath, recursive = TRUE, full.names = TRUE)
files

if (require(xcms)) {

## xcmsSet Summary
show(xs)

## Access raw data file
x <- xcmsRaw(files[1])
x

}
```

Index

*Topic **datasets**

 msdata, [1](#)

 msdata, [1](#)

 x (msdata), [1](#)

 xcmsRaw, [2](#)

 xcmsSet, [2](#)

 xs (msdata), [1](#)