

Package ‘AffyCompatible’

March 25, 2013

Type Package

Title Affymetrix GeneChip software compatibility

Version 1.18.2

Author Martin Morgan, Robert Gentleman

Maintainer Martin Morgan <mtmorgan@fhcrc.org>

Depends R (>= 2.7.0), XML (>= 2.8-1), RCurl (>= 0.8-1), methods

Imports Biostrings

Description This package provides an interface to Affymetrix chip annotation and sample attribute files. The package allows an easy way for users to download and manage local data bases of Affymetrix NetAffx annotation files. The package also provides access to GeneChip Operating System (GCOS) and GeneChip Command Console (AGCC)-compatible sample annotation files.

License Artistic-2.0

biocViews Infrastructure, Microarray, OneChannel

R topics documented:

AffyCompatible-package	2
accessors	2
AffyCompatible-class	3
Internal-class	3
NetAffxResource	4
NetAffxResource-class	5
readArr	7
readMage	8
readNetAffx	9
readXml	9
xclass	10
xmlValue.XMLAttributeValue	11
zzz-ARRCompatibility-class	12
zzz-DTTCompatibility-class	13
zzz-NetAffxCompatibility-class	14

Index	16
--------------	-----------

AffyCompatible-package *Affymetrix chip annotation and sample attribute management*

Description

This package manages Affymetrix NetAffx annotation files, and parses MAGE and ARR sample attribute files into R data structures.

NetAffx provides a way to manage the diverse annotation files provided by Affymetrix for their arrays.

ARR files provide attribute information about samples used in an experiment. These files are created, for instance, by the GeneChip Command Console or compliant applications, or

MAGE files also provide information about samples. Support in this package is for MAGE produced by Affymetrix GCOS compliant applications, and emphasizes sample description and processing.

Details

Package details are available in the DESCRIPTION file of the package, or with the R command `packageDescription("AffyCompatible")`

Overviews of how to use the functionality of this package are available in the package vignettes. See `vignette(package="AffyCompatible")`.

Author(s)

Written by Martin Morgan.

Maintainer: Martin Morgan <martin.t.morgan@gmail.com>

accessors

Access NetAffxResource-class content

Description

These functions exist to access the information represented in objects of `linkS4class{NetAffxResource}` and related classes.

Usage

```

annotationsFile(x)
directory(x)
password(x)
user(x)

```

Arguments

x An object for which the accessor is defined.

Value

See, e.g., `help('directory,NetAffxResource-method')` for values returned by specific methods.

Author(s)

Martin Morgan

See Also

[NetAffxResource](#)

AffyCompatible-class *Virtual classes defining the "AffyCompatibility" hierarchy*

Description

The classes `AffyCompatible`, `DTTCompatibility`, `ARRCompatibility`, and `NetAffxCompatibility` define the overall class hierarchy in this package. Information about these classes is not usually relevant to the user.

Objects from the Class

These are virtual classes: no objects may be created from them.

Methods

show signature(object = "AffyCompatible"): for each slot, provide the slot name, and then invoke show for the slot content.

Author(s)

Martin Morgan

Examples

```
showClass(getClass("NetAffxCompatibility",
  where=getNamespace("AffyCompatible")))
```

Internal-class *(Internal) Represent 'typed associations' from XML documents*

Description

This internal class represents 'typed associates' (vectors of identically-typed objects) that appear as vectors to the R user.

Objects from the Class

Objects are created internally, during XML document parsing.

Slots

.Data: Object of class "list" containing the association objects

.extends: Object of class "character" and length 1 describing the class each element in .Data must extend

Extends

Class "[list](#)", from data part. Class "[vector](#)", by class "[list](#)", distance 2.

Methods

[signature(x = ".TypedAssociation", i = "numeric", j = "missing"): Subset the typed association, returning a typed association with elements i.

[signature(x = ".TypedAssociation", i = "ANY", j = "ANY"): This method is present to stop access to subsetting methods for class list; it has no user-visible functionality.

show signature(object = ".TypedAssociation"): Compactly display the typed associate, e.g., for `.extneds=Foo` and length 5, display `Foo(5)`.

Author(s)

Martin Morgan

NetAffxResource

Create and access a NetAffxResource database

Description

Use this function to create and access a NetAffxResource database.

Usage

```
NetAffxResource(user = character(0),
  password = character(0),
  affxLicence = "FHCRC0607",
  directory = tempdir(), ...)
```

Arguments

user	Length 1 character vector containing the Affymetrix user name authorized to access the NetAffx resource.
password	Length 1 character vector containing the Affymetrix password of the user.
affxLicence	Length 1 character vector containing the Affymetrix licence (typically application specific, i.e., the specified licence covers AffyCompatible).
directory	(Optional) length 1 character vector providing the full path to the local directory where the NetAffx database and files will be managed. Defaults to a session-specific temporary directory.
...	Additional arguments, not currently used.

Details

NetAffxResource creates a [NetAffxResource](#) object that manages a data base of NetAffx resources. The data base is located on the user file system as indicated by the directory argument. It is populated when the user requests specific files.

Users should provide a directory argument, so that the data base can be reused across sessions. This means that large files are only downloaded once.

Additional details are provided in the vignette accompanying this package.

Value

[NetAffxResource](#).

Author(s)

Martin Morgan

See Also

[NetAffxResource](#). The vignette 'Annotations with NetAffx'.

NetAffxResource-class *Manage NetAffx annotations*

Description

This class is for accessing and retrieving NetAffx annotation files for Affymetrix chips. It requires that the user has established an account with Affymetrix. Detailed usage is described in the 'Annotations with NetAffx' vignette.

Objects from the Class

Objects are created by calls of the form `NetAffxResource(user="user@research.org", password="terces")`, as illustrated below and in the 'Annotations with NetAffx' vignette.

Slots

user: Object of class "character" containing the full user name (typically an email address) obtained when the user registers with Affymetrix.

password: Object of class "character" containing the password associated with the Affymetrix user name. The password is entered and stored in plain text, and should not be considered secure within Bioconductor.

affxLicence: Object of class "character" containing the licence, obtained from Affymetrix, for retrieving NetAffx annotations. The licence is usually software specific (e.g., for the Affy-Compatible package of Bioconductor). Users will not usually change the default value.

directory: Object of class "character", containing the local directory path where the NetAffx data base and annotation files will reside.

annotationsFile: Object of class "character" containing the file name of the NetAffx data base. Users will not usually change the default value.

affxUrl: Object of class "character" containing the URL used to download the NetAffx annotations file. Users will not usually change the default value.

Extends

Class "[NetAffxCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "NetAffxCompatibility", distance 2.

Methods

[[signature(x = "NetAffxResource", i = "numeric", j = "missing")
[[signature(x = "NetAffxResource", i = "character", j = "missing"): Retrieve an object representing all NetAffx annotations available for an array. The argument *i* may be a numeric index or character string corresponding to an array name returned by `names(x)`.

[[signature(x = "NetAffxResource", i = "numeric", j = "numeric")

[[signature(x = "NetAffxResource", i = "character", j = "numeric")

[[signature(x = "NetAffxResource", i = "numeric", j = "character")

[[signature(x = "NetAffxResource", i = "character", j = "character"): Retrieve an object representing a specific annotation for an array. The argument *i* may be a numeric index or character string corresponding to an array name returned by `names(x)`. The argument *j* may be a numeric index or character string corresponding to an annotation returned by `description(x[[i]])`.

affxLicence signature(x = "NetAffxResource"): the character string representation of the licence used to access the NetAffx site.

affxUrl signature(x = "NetAffxResource"): the url used to access the NetAffx site.

annotationsFile signature(x = "NetAffxResource"): the file name of NetAffx data base, as stored on the local disk of the user.

directory signature(x = "NetAffxResource"): the directory in which the NetAffx data base and annotation files are stored on the local disk of the user.

names signature(x = "NetAffxResource"): the names of the arrays for which annotations are available.

password signature(x = "NetAffxResource"): The password associated with the user name obtained when the user registered with Affymetrix. The password is entered and stored as plain text, and is not secure.

readAnnotation signature(netAffxResource = "NetAffxResource", array = "numeric", annotation = "numeric")

readAnnotation signature(netAffxResource = "NetAffxResource", array = "numeric", annotation = "character")

readAnnotation signature(netAffxResource = "NetAffxResource", array = "character", annotation = "numeric")

readAnnotation signature(netAffxResource = "NetAffxResource", array = "character", annotation = "character")

readAnnotation signature(netAffxResource = "NetAffxResource", array = "missing", annotation = "AffxAnnotation")

Read a specific annotation file as an R object. For the first four methods, *array*, if present, indexes or names the array while *annotation* indexes or names the annotation. All methods also accept *content = "logical"* to indicate whether the content of the file should be returned if possible (*content=TRUE*, default) or only the path to the local file (*content=FALSE*) and *update = "logical"* indicating whether the file should be read from disk if available (*update=FALSE*, default) or retrieved from the NetAffx web site (*update=TRUE*). `readAnnotation` parses CSV, Tabular, and PSI files to data frames, and reads FASTA files to a list of description / sequence pairs (see `Biostrings::readFASTA`).

The `...` argument is passed to corresponding read functions. CSV files are parsed using `read.csv(conn, ...)`, Tabular files with `read.delim(conn, ...)`, PSI files with `read.delim(conn, header=FALSE, skip=1, ...)` and FASTA with `readFASTA(conn, ...)`.

show signature(object = "NetAffxResource"): compactly show the NetAffxResource object.

user signature(x = "NetAffxResource"): The user name (typically email address) obtained when the user registers with the Affymetrix web site.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette ‘Annotations with NetAffx’. Affymetrix user login site <https://www.affymetrix.com/site/login/login.affx>.

[.Rprofile](#) for ways to automatically set variables (such as password, and directory) for each session.

Examples

```
## password must be obtained by the user
rsrc <- NetAffxResource(user="mtmorgan@fhcrc.org", password=password)
head(names(rsrc))
affxDescription(rsrc[["Bovine"]])
## paradigm to display all elements of 'annotation'
sapply(affxAnnotation(rsrc[["Bovine"]]), force)
## obtain and read a single annotation
anno <- rsrc[["Porcine", "Annotations, CSV format"]]
df <- readAnnotation(rsrc, annotation=anno)
```

readArr

Read ARR sample attribute files

Description

This functions reads the content of Affymetrix ARR sample attribute files created by GeneChip(R) Command Console. Objects are represented as R classes. Navigation and extraction of data are described in the vignette ‘Retrieving MAGE and ARR sample attributes’.

Usage

```
readArr(flz, ...)
```

Arguments

`flz` A character vector containing the paths to files to be read.
`...` Additional arguments, currently unused.

Value

An object or list of objects of class `ArrayFileSet`, depending on whether one or several file paths are provided.

Author(s)

Martin Morgan

See Also

The vignette ‘Retrieving MAGE and ARR sample attributes’ for detailed illustration of use. [xclass](#) and [readXml](#) for advanced access methods.

Examples

```
dir <- system.file("extdata", "ARR", package="AffyCompatible")
fls <- list.files(dir, full=TRUE)
readArr(fls[[1]])
readArr(fls)
```

readMage

Read DTT (MAGE) sample attribute files

Description

This functions reads the content of Affymetrix DTT (MAGE) sample attribute files created by GCOS. Objects are represented as R classes. Navigation and extraction of data are described in the vignette ‘Retrieving MAGE and ARR sample attributes’.

Usage

```
readMage(fls, ...)
```

Arguments

`fls` A character vector containing the paths to files to be read.
`...` Additional arguments, currently unused.

Value

An object or list of objects of class `ArrayFileSet`, depending on whether one or several file paths are provided.

Author(s)

Martin Morgan

See Also

The vignette ‘Retrieving MAGE and ARR sample attributes’ for detailed illustration of use.
[xclass](#) and [readXml](#) for advanced access methods.

Examples

```
dir <- system.file("extdata", "DTT", package="AffyCompatible")
fls <- list.files(dir, full=TRUE)
basename(fls)
readArr(fls[[1]])
```

readNetAffx	<i>Update or read a NetAffx data base</i>
-------------	-------------------------------------------

Description

This function provides a way to read or update a NetAffx data base. The usual use for this is to force the data base description to be updated (by accessing the Affymetrix web site) or to retrieve the entire list of NetAffx resource data base. A more normal use case is to

Usage

```
readNetAffx(netAffxResource, update = FALSE, ...)
```

Arguments

netAffxResource	An object of class NetAffxResource containing the path to the local data base, and permissions for access the Affymetrix NetAffx resource.
update	An object of logical(1) indicating whether the data base is to be updated (TRUE) or not.
...	Additional arguments, not currently used.

Value

An object representing the annotations found; see [NetAffxResource](#) for additional details.

Author(s)

Martin Morgan

Examples

```
## retrieve all annotations
## Not run:
example("NetAffxResource-class")
readNetAffx(rsrc)

## End(Not run)
```

readXml	<i>(Advanced) Read XML files of Affymetrix sample attributes or annotation</i>
---------	--------------------------------------------------------------------------------

Description

This functions reads the content of Affymetrix DTT (MAGE) or ARR sample files, or AffxNetAffx-AnnotFile array annotation files, return an internal representation of the underlying XML. The return value can be used with [xclass](#) to rapidly create R objects representing deeply nested nodes of the document. The return value can also be used with `xpathApply` in the XML package for access to character-based representation of the data.

Usage

```
readXml(fl, ...)
```

Arguments

`fl` A length-1 character vector giving the full path to the file to be parsed.
`...` Additional arguments, passed to `xmlTreeParse`.

Value

An object of class `XMLInternalDocument`.

Author(s)

Martin Morgan

See Also

The vignettes ‘Retrieving MAGE and ARR sample attributes’ and ‘Annotation retrieval with NetAffxResource’ for detailed illustration of use.

[xclass](#) for instantiating objects from xpath queries.

Examples

```
mageDir <- system.file("extdata", "DTT", package="AffyCompatible")
mageFiles <- list.files(mageDir, full=TRUE)
xml <- readXml(mageFiles[[1]])
xclass(xml, "//MeasuredBioAssay")[[1]]
```

xclass

(Advanced) Construct R objects from XML document queries

Description

This function takes an `xml` object, created from [readXml](#), and an xpath query referencing an element of the document, and returns an R object representing the element.

The xpath query must retrieve an element, not an attribute.

Usage

```
xclass(doc, xpathq, prefix = character(0))
```

Arguments

`doc` An `XMLInternalDocument` retrieved using [readXml](#).
`xpathq` A length-1 character string representing an xpath query into the document.
`prefix` A prefix used to map, when necessary, between the XML class name and the corresponding R class. For NetAffx-based classes, the prefix is ‘Affx’; no prefix is required for DTT or ARR documents.

Value

An R object corresponding to the class(es) referenced by the xpath query.

Author(s)

Martin Morgan

References

<http://www.w3.org/TR/xpath> provides a very useful description of xpath. Additional references are in the vignette ‘Retrieving MAGE and ARR sample attributes’.

See Also

The vignettes ‘Retrieving MAGE and ARR sample attributes’ and ‘Annotation retrieval with NetAffxResource’ for detailed illustration of use.

[readXml](#) for appropriately reading XML documents.

Examples

```
example(readXml)
```

xmlValue.XMLAttributeValue

Additional functions defined in AffyCompatible

Description

This page documents additional, miscellaneous, functions defined by AffyCompatible

Usage

```
## S3 method for class 'XMLAttributeValue'
xmlValue(x, ignoreComments = FALSE, recursive = TRUE,
         encoding=XML:::CE_NATIVE)
```

Arguments

x	An object of class XMLAttributeValue.
ignoreComments	Logical; ignored by the xmlValue.XMLAttributeValue method
recursive	ignored
encoding	File encoding

Value

A character vector (usually length 1) containing the character representation of the XML attribute.

zzz-ARRCompatibility-class

(Auto-generated) classes and methods for navigating ARR sample attributes.

Description

These classes and methods facilitate navigation of the Affymetrix ARR sample attribute files. ARR files are generated by the Affymetrix AGCC (GeneChip(R) Command Console) application. Specifically, the classes are R representations of the Affymetrix classes defined in ArraySetAndTemplate-File.dtd. The 'show' methods for these classes indicate how objects are to be accessed; details are in the 'Retrieving MAGE and ARR sample attributes' vignette.

Objects from the Class

Objects are created internally, by parsing XML retrieved from ARR, and are not created directly by the user.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "[ARRCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "ARRCompatibility", distance 2.

Methods

'Accessor' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke `accessor(object)`.

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using `[, []`. A useful paradigm for displaying all instances is illustrated in the vignette.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Retrieving MAGE and ARR sample attributes'.

Examples

```
## defined and auto-generated classes
names(slot(getClass("ARRCompatibility",
  where=getNamespace("AffyCompatible")),
  "subclasses"))
```

zzz-DTTCompatibility-class

(Auto-generated) classes and methods for navigating DTT sample attributes.

Description

These classes and methods facilitate navigation of the Affymetrix DTT sample attribute files. DTT files are generated by the Affymetrix GCOS application. Specifically, the classes are R representations of the Affymetrix classes defined in MAGE-ML.dtd. The 'show' methods for these classes indicate how objects are to be accessed; details are in the 'Retrieving MAGE and ARR sample attributes' vignette.

Objects from the Class

Objects are created internally, by parsing XML retrieved from ARR, and are not created directly by the user.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "[DTTCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "DTTCompatibility", distance 2.

Methods

'Accessor' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke `accessor(object)`.

Some objects contains several instances of another object. These 'typed associations' are represented by the [.TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using `[, [[`. A useful paradigm for displaying all instances is illustrated in the vignette.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Retrieving MAGE and ARR sample attributes'.

Examples

```
## defined and auto-generated classes
names(slot(getClass("DTTCompatibility",
  where=getNamespace("AffyCompatible")),
  "subclasses"))
```

zzz-NetAffxCompatibility-class

(Auto-generated) classes and methods for navigating NetAffx resources

Description

These classes and methods facilitate navigation of the NetAffx resource data base. They are R representations of the Affymetrix classes defined in NetAffxAnnotFileList.dtd. The 'show' methods for these class indicate how objects are to be accessed; details are in the 'Annotations with NetAffx' vignette.

The auto-generated classes are discoverable as indicated in the example below.

Objects from the Class

Objects are created internally, by parsing XML retrieved from the local NetAffx data base or from the Affymetrix web site.

Slots

Objects have a variety of slots. To discover slot names and content type, see the slotNames portion of the example, below.

Extends

Class "[NetAffxCompatibility](#)", directly. Class "[AffyCompatible](#)", by class "[NetAffxCompatibility](#)", distance 2.

Methods

'Read' methods are defined on auto-generated classes to access slot content.

Accessors are apparent from displaying the object; each line of the display is of the form accessor: content. To retrieve the content of an instance object, invoke accessor(object).

Some objects contains several instances of another object. These 'typed associations' are represented by the [TypedAssociation](#) class, and are navigated or subset as one would a standard R vector, e.g., using [, [[]]. A useful paradigm for displaying all instances is illustrated on the [NetAffxResource](#) page.

Author(s)

Martin Morgan <mtmorgan@fhcrc.org>

See Also

Vignette 'Annotations with NetAffx'. Affymetrix user login site <https://www.affymetrix.com/site/login/login.affx>.

[NetAffxResource](#) for creating and managing the NetAffx data base.

[NetAffxResource](#) for navigating the data base.

Examples

```
## defined and auto-generated classes
getClass("NetAffxCompatibility", where=getNamespace("AffyCompatible"))

## slots of an auto-generated class
getSlots(getClass("AffxNetAffxAnnotFileList",
  where=getNamespace("AffyCompatible")))

## Not run:
example("NetAffxrResource-class")

## End(Not run)
```

Index

*Topic **classes**

- AffyCompatible-class, 3
- Internal-class, 3
- NetAffxResource-class, 5
- readArr, 7
- readMage, 8
- zzz-ARRCompatibility-class, 12
- zzz-DTTCompatibility-class, 13
- zzz-NetAffxCompatibility-class, 14

*Topic **database**

- NetAffxResource, 4
- readNetAffx, 9

*Topic **manip**

- readXml, 9
- xclass, 10
- xmlValue.XMLAttributeValue, 11

*Topic **methods**

- accessors, 2

*Topic **package**

- AffyCompatible-package, 2

.Rprofile, 7

.TypedAssociation, 12–14

.TypedAssociation-class (Internal-class), 3

[,TypedAssociation,ANY,ANY-method
(Internal-class), 3

[,TypedAssociation,numeric,missing,ANY-method
(Internal-class), 3

[,TypedAssociation,numeric,missing-method
(Internal-class), 3

[[,NetAffxResource,character,character-method
(NetAffxResource-class), 5

[[,NetAffxResource,character,missing-method
(NetAffxResource-class), 5

[[,NetAffxResource,character,numeric-method
(NetAffxResource-class), 5

[[,NetAffxResource,numeric,character-method
(NetAffxResource-class), 5

[[,NetAffxResource,numeric,missing-method
(NetAffxResource-class), 5

[[,NetAffxResource,numeric,numeric-method
(NetAffxResource-class), 5

accession (zzz-DTTCompatibility-class), 13

accession,DatabaseEntry-method

(zzz-DTTCompatibility-class), 13

accessions_ assnlist

(zzz-DTTCompatibility-class), 13

accessions_ assnlist,BibliographicReference-method

(zzz-DTTCompatibility-class), 13

Accessions_ assnlist-class

(zzz-DTTCompatibility-class), 13

accessionVersion

(zzz-DTTCompatibility-class), 13

accessionVersion,DatabaseEntry-method

(zzz-DTTCompatibility-class), 13

accessors, 2

action (zzz-DTTCompatibility-class), 13

action,Audit-method

(zzz-DTTCompatibility-class), 13

action_ assn (zzz-DTTCompatibility-class),
13

action_ assn,Treatment-method

(zzz-DTTCompatibility-class), 13

Action_ assn-class

(zzz-DTTCompatibility-class), 13

actionMeasurement_ assn

(zzz-DTTCompatibility-class), 13

actionMeasurement_ assn,Treatment-method

(zzz-DTTCompatibility-class), 13

ActionMeasurement_ assn-class

(zzz-DTTCompatibility-class), 13

activityDate (zzz-DTTCompatibility-class),
13

activityDate,ProtocolApplication-method

(zzz-DTTCompatibility-class), 13

address (zzz-DTTCompatibility-class), 13

address,Organization-method

(zzz-DTTCompatibility-class), 13

address,Person-method

(zzz-DTTCompatibility-class), 13

adjustments_ assnlist

(zzz-DTTCompatibility-class), 13

adjustments_ assnlist,ArrayManufactureDeviation-method

(zzz-DTTCompatibility-class), 13

Adjustments_ assnlist-class

(zzz-DTTCompatibility-class), 13

- affiliation_ assnref
(zzz-DTTCompatibility-class), 13
- affiliation_ assnref,Person-method
(zzz-DTTCompatibility-class), 13
- Affiliation_ assnref-class
(zzz-DTTCompatibility-class), 13
- affxAnnotation
(zzz-NetAffxCompatibility-class), 14
- affxAnnotation,AffxArray-method
(zzz-NetAffxCompatibility-class), 14
- AffxAnnotation-class
(zzz-NetAffxCompatibility-class), 14
- affxArray (zzz-NetAffxCompatibility-class), 14
- affxArray,AffxNetAffxAnnotFileList-method
(zzz-NetAffxCompatibility-class), 14
- affxArray,Array_ assnlist-method
(zzz-DTTCompatibility-class), 13
- AffxArray-class
(zzz-NetAffxCompatibility-class), 14
- affxCategory
(zzz-DTTCompatibility-class), 13
- affxCategory,OntologyEntry-method
(zzz-DTTCompatibility-class), 13
- affxComment
(zzz-NetAffxCompatibility-class), 14
- affxComment,AffxAnnotation-method
(zzz-NetAffxCompatibility-class), 14
- affxComment,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- AffxComment-class
(zzz-NetAffxCompatibility-class), 14
- affxCompression
(zzz-NetAffxCompatibility-class), 14
- affxCompression,AffxUrl-method
(zzz-NetAffxCompatibility-class), 14
- affxCrc32 (zzz-NetAffxCompatibility-class), 14
- affxCrc32,AffxFile-method
(zzz-NetAffxCompatibility-class), 14
- affxCrc32,AffxUrl-method
(zzz-NetAffxCompatibility-class), 14
- affxDate (zzz-NetAffxCompatibility-class), 14
- affxDate,AffxFile-method
(zzz-NetAffxCompatibility-class), 14
- affxDate,Audit-method
(zzz-DTTCompatibility-class), 13
- affxDescription
(zzz-NetAffxCompatibility-class), 14
- affxDescription,AffxAnnotation-method
(zzz-NetAffxCompatibility-class), 14
- affxDescription,AffxArray-method
(zzz-NetAffxCompatibility-class), 14
- affxEnd (zzz-DTTCompatibility-class), 13
- affxEnd,CompositePosition-method
(zzz-DTTCompatibility-class), 13
- affxEnd,ReporterPosition-method
(zzz-DTTCompatibility-class), 13
- affxEnd,SequencePosition-method
(zzz-DTTCompatibility-class), 13
- affxFile (zzz-NetAffxCompatibility-class), 14
- affxFile,AffxAnnotation-method
(zzz-NetAffxCompatibility-class), 14
- AffxFile-class
(zzz-NetAffxCompatibility-class), 14
- affxImage (zzz-DTTCompatibility-class), 13
- affxImage,PhysicalBioAssayData_ assnlist-method
(zzz-DTTCompatibility-class), 13
- affxLicence
(zzz-NetAffxCompatibility-class), 14
- affxLicence,NetAffxResource-method
(NetAffxResource-class), 5
- affxName (zzz-NetAffxCompatibility-class), 14
- affxName,AffxArray-method
(zzz-NetAffxCompatibility-class), 14
- affxName,AffxFile-method
(zzz-NetAffxCompatibility-class), 14
- AffxNetAffxAnnotFileList-class
(zzz-NetAffxCompatibility-class), 14

- affxOrder (*zzz-DTTCompatibility-class*), 13
- affxOrder,BioDataCube-method
(*zzz-DTTCompatibility-class*), 13
- affxOrder,Treatment-method
(*zzz-DTTCompatibility-class*), 13
- affxPerson (*zzz-DTTCompatibility-class*),
13
- affxPerson,Contact_ assnlist-method
(*zzz-DTTCompatibility-class*), 13
- affxPerson,Person_ assnlist-method
(*zzz-DTTCompatibility-class*), 13
- affxRow (*zzz-DTTCompatibility-class*), 13
- affxRow,FeatureLocation-method
(*zzz-DTTCompatibility-class*), 13
- affxRow,Zone-method
(*zzz-DTTCompatibility-class*), 13
- affxSequence
(*zzz-DTTCompatibility-class*), 13
- affxSequence,BioSequence-method
(*zzz-DTTCompatibility-class*), 13
- affxSize (*zzz-NetAffxCompatibility-class*),
14
- affxSize,AffxFile-method
(*zzz-NetAffxCompatibility-class*),
14
- affxSize,AffxUrl-method
(*zzz-NetAffxCompatibility-class*),
14
- affxStart (*zzz-DTTCompatibility-class*), 13
- affxStart,CompositePosition-method
(*zzz-DTTCompatibility-class*), 13
- affxStart,ReporterPosition-method
(*zzz-DTTCompatibility-class*), 13
- affxStart,SequencePosition-method
(*zzz-DTTCompatibility-class*), 13
- affxStatus
(*zzz-NetAffxCompatibility-class*),
14
- affxStatus,AffxNetAffxAnnotFileList-method
(*zzz-NetAffxCompatibility-class*),
14
- AffxStatus-class
(*zzz-NetAffxCompatibility-class*),
14
- affxText (*zzz-DTTCompatibility-class*), 13
- affxText,Description-method
(*zzz-DTTCompatibility-class*), 13
- affxText,Protocol-method
(*zzz-DTTCompatibility-class*), 13
- affxTitle (*zzz-DTTCompatibility-class*), 13
- affxTitle,BibliographicReference-method
(*zzz-DTTCompatibility-class*), 13
- affxTitle,Protocol-method
(*zzz-DTTCompatibility-class*), 13
- affxType (*zzz-NetAffxCompatibility-class*),
14
- affxType,AffxAnnotation-method
(*zzz-NetAffxCompatibility-class*),
14
- affxType,AffxFile-method
(*zzz-NetAffxCompatibility-class*),
14
- affxUrl (*zzz-NetAffxCompatibility-class*),
14
- affxUrl,AffxFile-method
(*zzz-NetAffxCompatibility-class*),
14
- affxUrl,NetAffxResource-method
(*NetAffxResource-class*), 5
- AffxUrl-class
(*zzz-NetAffxCompatibility-class*),
14
- affyBarcode (*zzz-ARRCompatibility-class*),
12
- affyBarcode,PhysicalArray-method
(*zzz-ARRCompatibility-class*), 12
- AffyCompatible, 5, 12–14
- AffyCompatible (*AffyCompatible-package*),
2
- AffyCompatible-class, 3
- AffyCompatible-package, 2
- analysisResults_ assnreflist
(*zzz-DTTCompatibility-class*), 13
- analysisResults_ assnreflist,Experiment-method
(*zzz-DTTCompatibility-class*), 13
- AnalysisResults_ assnreflist-class
(*zzz-DTTCompatibility-class*), 13
- annotations_ assnlist
(*zzz-DTTCompatibility-class*), 13
- annotations_ assnlist,Description-method
(*zzz-DTTCompatibility-class*), 13
- annotations_ assnlist,ExperimentalFactor-method
(*zzz-DTTCompatibility-class*), 13
- Annotations_ assnlist-class
(*zzz-DTTCompatibility-class*), 13
- annotationsFile (accessors), 2
- annotationsFile,NetAffxResource-method
(*NetAffxResource-class*), 5
- Array-class (*zzz-DTTCompatibility-class*),
13
- array_ assnlist
(*zzz-DTTCompatibility-class*), 13
- array_ assnlist,Array_ package-method
(*zzz-DTTCompatibility-class*), 13

- Array_ assnlist-class
(zzz-DTTCompatibility-class), 13
- array_ assnref
(zzz-DTTCompatibility-class), 13
- array_ assnref,BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- array_ assnref,Hybridization-method
(zzz-DTTCompatibility-class), 13
- Array_ assnref-class
(zzz-DTTCompatibility-class), 13
- array_ package
(zzz-DTTCompatibility-class), 13
- array_ package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- Array_ package-class
(zzz-DTTCompatibility-class), 13
- array_ ref (zzz-DTTCompatibility-class),
13
- array_ ref,Array_ assnref-method
(zzz-DTTCompatibility-class), 13
- array_ ref,Arrays_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- Array_ ref-class
(zzz-DTTCompatibility-class), 13
- arrayAttribute
(zzz-ARRCompatibility-class), 12
- arrayAttribute,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- ArrayAttribute-class
(zzz-ARRCompatibility-class), 12
- arrayDesign (zzz-DTTCompatibility-class),
13
- arrayDesign,ArrayDesign_ assnlist-method
(zzz-DTTCompatibility-class), 13
- ArrayDesign-class
(zzz-DTTCompatibility-class), 13
- arrayDesign_ assnlist
(zzz-DTTCompatibility-class), 13
- arrayDesign_ assnlist,ArrayDesign_ package-method
(zzz-DTTCompatibility-class), 13
- ArrayDesign_ assnlist-class
(zzz-DTTCompatibility-class), 13
- arrayDesign_ assnref
(zzz-DTTCompatibility-class), 13
- arrayDesign_ assnref,Array-method
(zzz-DTTCompatibility-class), 13
- ArrayDesign_ assnref-class
(zzz-DTTCompatibility-class), 13
- arrayDesign_ package
(zzz-DTTCompatibility-class), 13
- arrayDesign_ package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- ArrayDesign_ package-class
(zzz-DTTCompatibility-class), 13
- arrayDesign_ ref
(zzz-DTTCompatibility-class), 13
- arrayDesign_ ref,ArrayDesign_ assnref-method
(zzz-DTTCompatibility-class), 13
- ArrayDesign_ ref-class
(zzz-DTTCompatibility-class), 13
- arrayGroup,ArrayGroup_ assnlist-method
(zzz-DTTCompatibility-class), 13
- ArrayGroup-class
(zzz-DTTCompatibility-class), 13
- arrayGroup_ assnlist
(zzz-DTTCompatibility-class), 13
- arrayGroup_ assnlist,Array_ package-method
(zzz-DTTCompatibility-class), 13
- ArrayGroup_ assnlist-class
(zzz-DTTCompatibility-class), 13
- arrayGroup_ assnref
(zzz-DTTCompatibility-class), 13
- arrayGroup_ assnref,Array-method
(zzz-DTTCompatibility-class), 13
- ArrayGroup_ assnref-class
(zzz-DTTCompatibility-class), 13
- arrayGroup_ ref
(zzz-DTTCompatibility-class), 13
- arrayGroup_ ref,ArrayGroup_ assnref-method
(zzz-DTTCompatibility-class), 13
- ArrayGroup_ ref-class
(zzz-DTTCompatibility-class), 13
- arrayIdentifier
(zzz-DTTCompatibility-class), 13
- arrayIdentifier,Array-method
(zzz-DTTCompatibility-class), 13
- arrayManufacture
(zzz-DTTCompatibility-class), 13
- arrayManufacture,ArrayManufacture_ assnlist-method
(zzz-DTTCompatibility-class), 13
- ArrayManufacture-class
(zzz-DTTCompatibility-class), 13
- arrayManufacture_ assnlist
(zzz-DTTCompatibility-class), 13
- arrayManufacture_ assnlist,Array_ package-method
(zzz-DTTCompatibility-class), 13
- ArrayManufacture_ assnlist-class
(zzz-DTTCompatibility-class), 13
- arrayManufacture_ ref
(zzz-DTTCompatibility-class), 13
- arrayManufacture_ ref,Information_ assnref-method
(zzz-DTTCompatibility-class), 13

- ArrayManufacture_ref-class
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviation
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviation,ArrayManufactureDeviations_assnlist
(zzz-DTTCompatibility-class), 13
- ArrayManufactureDeviation-class
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviations_assnlist
(zzz-DTTCompatibility-class), 13
- arrayManufactureDeviations_assnlist,Array-method
(zzz-DTTCompatibility-class), 13
- ArrayManufactureDeviations_assnlist-class
(zzz-DTTCompatibility-class), 13
- arrayManufacturers_assnreflist
(zzz-DTTCompatibility-class), 13
- arrayManufacturers_assnreflist,ArrayManufactureDeviations_assnlist
(zzz-DTTCompatibility-class), 13
- ArrayManufacturers_assnreflist-class
(zzz-DTTCompatibility-class), 13
- arrayName (zzz-ARRCompatibility-class), 12
- arrayName,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- arrays_assnreflist
(zzz-DTTCompatibility-class), 13
- arrays_assnreflist,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- arrays_assnreflist,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- Arrays_assnreflist-class
(zzz-DTTCompatibility-class), 13
- ArraySetFile-class
(zzz-ARRCompatibility-class), 12
- arraySpacingX
(zzz-DTTCompatibility-class), 13
- arraySpacingX,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- arraySpacingY
(zzz-DTTCompatibility-class), 13
- arraySpacingY,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- arrayXOrigin
(zzz-DTTCompatibility-class), 13
- arrayXOrigin,Array-method
(zzz-DTTCompatibility-class), 13
- arrayYOrigin
(zzz-DTTCompatibility-class), 13
- arrayYOrigin,Array-method
(zzz-DTTCompatibility-class), 13
- ARRCompatibility, 12
- ARRCompatibility-class
(AffyCompatible-class), 3
- associations_assnlist
(zzz-DTTCompatibility-class), 13
- associations_assnlist,OntologyEntry-method
(zzz-DTTCompatibility-class), 13
- Associations_assnlist-class
(zzz-DTTCompatibility-class), 13
- audit (zzz-DTTCompatibility-class), 13
- audit,AuditTrail_assnlist-method
(zzz-DTTCompatibility-class), 13
- Audit-class (zzz-DTTCompatibility-class), 13
- auditAndSecurity_package
(zzz-DTTCompatibility-class), 13
- auditAndSecurity_package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- auditAndSecurity_package-class
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Array-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Audit-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioAssayDimension-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioSample-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioSequence-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,BioSource-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Channel-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,CompositeCompositeMap-method
(zzz-DTTCompatibility-class), 13

- auditTrail_assnlist,CompositeGroup-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Compound-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Database-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,DerivedBioAssay-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Description-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Error-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Experiment-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ExperimentDesign-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,FactorValue-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Failed-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Feature-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,FeatureDimension-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,FeatureExtraction-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Fiducial-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Hardware-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,HardwareApplication-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Hybridization-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Image-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Node-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,NodeContents-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Organization-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Parameter-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Person-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Protocol-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,PValue-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Ratio-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Reporter-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ReporterDimension-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class), 13

- auditTrail_assnlist,Security-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Software-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Transformation-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Treatment-method
(zzz-DTTCompatibility-class), 13
- auditTrail_assnlist,Zone-method
(zzz-DTTCompatibility-class), 13
- AuditTrail_assnlist-class
(zzz-DTTCompatibility-class), 13
- authors (zzz-DTTCompatibility-class), 13
- authors,BibliographicReference-method
(zzz-DTTCompatibility-class), 13

- barcode (zzz-DTTCompatibility-class), 13
- barcode,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- basis (zzz-DTTCompatibility-class), 13
- basis,SeqFeature-method
(zzz-DTTCompatibility-class), 13
- bibliographicReference
(zzz-DTTCompatibility-class), 13
- bibliographicReference,BibliographicReferences_assnlist-method
(zzz-DTTCompatibility-class), 13
- BibliographicReference-class
(zzz-DTTCompatibility-class), 13
- bibliographicReferences_assnlist
(zzz-DTTCompatibility-class), 13
- bibliographicReferences_assnlist,Description-method
(zzz-DTTCompatibility-class), 13
- BibliographicReferences_assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssay_assnlist
(zzz-DTTCompatibility-class), 13
- bioAssay_assnlist,BioAssay_package-method
(zzz-DTTCompatibility-class), 13
- BioAssay_assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssay_assnref
(zzz-DTTCompatibility-class), 13
- bioAssay_assnref,BioAssayTuple-method
(zzz-DTTCompatibility-class), 13
- BioAssay_assnref-class
(zzz-DTTCompatibility-class), 13
- bioAssay_package
(zzz-DTTCompatibility-class), 13
- bioAssay_package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- BioAssay_package-class
(zzz-DTTCompatibility-class), 13
- BioAssay_ref-class
(zzz-DTTCompatibility-class), 13
- BioAssayCreation
(zzz-DTTCompatibility-class), 13
- bioAssayCreation,BioAssayCreation_assn-method
(zzz-DTTCompatibility-class), 13
- BioAssayCreation-class
(zzz-DTTCompatibility-class), 13
- bioAssayCreation_assn
(zzz-DTTCompatibility-class), 13
- bioAssayCreation_assn,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- BioAssayCreation_assn-class
(zzz-DTTCompatibility-class), 13
- bioAssayData_assnlist
(zzz-DTTCompatibility-class), 13
- bioAssayData_assnlist,BioAssayData_package-method
(zzz-DTTCompatibility-class), 13
- BioAssayData_assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssayData_assnreflist
(zzz-DTTCompatibility-class), 13
- bioAssayData_assnreflist,Experiment-method
(zzz-DTTCompatibility-class), 13
- BioAssayData_assnreflist-class
(zzz-DTTCompatibility-class), 13
- bioAssayData_package
(zzz-DTTCompatibility-class), 13
- bioAssayData_package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- BioAssayData_package-class
(zzz-DTTCompatibility-class), 13
- BioAssayData_ref-class
(zzz-DTTCompatibility-class), 13
- bioAssayDataCluster
(zzz-DTTCompatibility-class), 13
- bioAssayDataCluster,BioAssayDataCluster_assnlist-method
(zzz-DTTCompatibility-class), 13
- BioAssayDataCluster-class
(zzz-DTTCompatibility-class), 13
- bioAssayDataCluster_assnlist
(zzz-DTTCompatibility-class), 13
- bioAssayDataCluster_assnlist,HigherLevelAnalysis_package-method
(zzz-DTTCompatibility-class), 13

- BioAssayDataCluster_ assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssayDataCluster_ref
(zzz-DTTCompatibility-class), 13
- bioAssayDataCluster_ref,AnalysisResults_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- BioAssayDataCluster_ref-class
(zzz-DTTCompatibility-class), 13
- bioAssayDataSources_ assnreflist
(zzz-DTTCompatibility-class), 13
- bioAssayDataSources_ assnreflist,Transformation-method
(zzz-DTTCompatibility-class), 13
- BioAssayDataSources_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- bioAssayDimension
(zzz-DTTCompatibility-class), 13
- bioAssayDimension,BioAssayDimension_ assnlist-method
(zzz-DTTCompatibility-class), 13
- BioAssayDimension-class
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ assnlist
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ assnlist,BioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- BioAssayDimension_ assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ assnref
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ assnref,DerivedBioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ assnref,MeasuredBioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ assnref,NodeContents-method
(zzz-DTTCompatibility-class), 13
- BioAssayDimension_ assnref-class
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ref
(zzz-DTTCompatibility-class), 13
- bioAssayDimension_ref,BioAssayDimension_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- BioAssayDimension_ref-class
(zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_ assnreflist
(zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_ assnreflist,DerivedBioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_ assnreflist,MeasuredBioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- bioAssayFactorValues_ assnreflist,PhysicalBioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- BioAssayFactorValues_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- bioAssayMap
(zzz-DTTCompatibility-class), 13
- bioAssayMap,BioAssayMap_ assnlist-method
(zzz-DTTCompatibility-class), 13
- BioAssayMap-class
(zzz-DTTCompatibility-class), 13
- bioAssayMap_ assnlist
(zzz-DTTCompatibility-class), 13
- bioAssayMap_ assnlist,BioAssayData_ package-method
(zzz-DTTCompatibility-class), 13
- BioAssayMap_ assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssayMap_ref
(zzz-DTTCompatibility-class), 13
- bioAssayMap_ref,BioAssayMaps_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- bioAssayMap_ref,DerivedBioAssayMap_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- BioAssayMap_ref-class
(zzz-DTTCompatibility-class), 13
- bioAssayMapping
(zzz-DTTCompatibility-class), 13
- bioAssayMapping,BioAssayMapping_ assnlist-method
(zzz-DTTCompatibility-class), 13
- BioAssayMapping-class
(zzz-DTTCompatibility-class), 13
- bioAssayMapping_ assn
(zzz-DTTCompatibility-class), 13
- bioAssayMapping_ assn,Transformation-method
(zzz-DTTCompatibility-class), 13
- BioAssayMapping_ assn-class
(zzz-DTTCompatibility-class), 13
- bioAssayMaps_ assnreflist
(zzz-DTTCompatibility-class), 13
- bioAssayMaps_ assnreflist,BioAssayMapping-method
(zzz-DTTCompatibility-class), 13
- BioAssayMaps_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- bioAssayMapTarget_ assnref
(zzz-DTTCompatibility-class), 13
- bioAssayMapTarget_ assnref,BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- BioAssayMapTarget_ assnref-class
(zzz-DTTCompatibility-class), 13
- BioAssays_ assnreflist
(zzz-DTTCompatibility-class), 13
- bioAssays_ assnreflist,BioAssayDimension-method
(zzz-DTTCompatibility-class), 13
- bioAssays_ assnreflist,Experiment-method
(zzz-DTTCompatibility-class), 13
- BioAssays_ assnreflist-class
(zzz-DTTCompatibility-class), 13

- bioAssayTreatment
(zzz-DTTCompatibility-class), 13
- bioAssayTreatment,BioAssayTreatments_ assnlist-method
(zzz-DTTCompatibility-class), 13
- BioAssayTreatment-class
(zzz-DTTCompatibility-class), 13
- bioAssayTreatments_ assnlist
(zzz-DTTCompatibility-class), 13
- bioAssayTreatments_ assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- BioAssayTreatments_ assnlist-class
(zzz-DTTCompatibility-class), 13
- bioAssayTuple
(zzz-DTTCompatibility-class), 13
- bioAssayTuple,BioAssayTuples_ assnlist-method
(zzz-DTTCompatibility-class), 13
- BioAssayTuple-class
(zzz-DTTCompatibility-class), 13
- bioAssayTuples_ assnlist
(zzz-DTTCompatibility-class), 13
- bioAssayTuples_ assnlist,BioDataTuples-method
(zzz-DTTCompatibility-class), 13
- BioAssayTuples_ assnlist-class
(zzz-DTTCompatibility-class), 13
- bioDataCube
(zzz-DTTCompatibility-class), 13
- bioDataCube,BioDataValues_ assn-method
(zzz-DTTCompatibility-class), 13
- BioDataCube-class
(zzz-DTTCompatibility-class), 13
- bioDataTuples
(zzz-DTTCompatibility-class), 13
- bioDataTuples,BioDataValues_ assn-method
(zzz-DTTCompatibility-class), 13
- BioDataTuples-class
(zzz-DTTCompatibility-class), 13
- bioDataValues_ assn
(zzz-DTTCompatibility-class), 13
- bioDataValues_ assn,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- bioDataValues_ assn,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- BioDataValues_ assn-class
(zzz-DTTCompatibility-class), 13
- bioEvent_ package
(zzz-DTTCompatibility-class), 13
- bioEvent_ package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- BioEvent_ package-class
(zzz-DTTCompatibility-class), 13
- biologicalCharacteristics_ assnreflist
(zzz-DTTCompatibility-class), 13
- biologicalCharacteristics_ assnreflist,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- biologicalCharacteristics_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- bioMaterial_ assnlist
(zzz-DTTCompatibility-class), 13
- bioMaterial_ assnlist,BioMaterial_ package-method
(zzz-DTTCompatibility-class), 13
- BioMaterial_ assnlist-class
(zzz-DTTCompatibility-class), 13
- bioMaterial_ assnref
(zzz-DTTCompatibility-class), 13
- bioMaterial_ assnref,BioMaterialMeasurement-method
(zzz-DTTCompatibility-class), 13
- bioMaterial_ assnref,ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- bioMaterial_ assnref,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- BioMaterial_ assnref-class
(zzz-DTTCompatibility-class), 13
- bioMaterial_ package
(zzz-DTTCompatibility-class), 13
- bioMaterial_ package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- BioMaterial_ package-class
(zzz-DTTCompatibility-class), 13
- BioMaterial_ ref-class
(zzz-DTTCompatibility-class), 13
- bioMaterialMeasurement
(zzz-DTTCompatibility-class), 13
- bioMaterialMeasurement,SourceBioMaterialMeasurements_ assn
(zzz-DTTCompatibility-class), 13
- BioMaterialMeasurement-class
(zzz-DTTCompatibility-class), 13
- bioMaterialPlateCol
(zzz-DTTCompatibility-class), 13
- bioMaterialPlateCol,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- bioMaterialPlateIdentifier
(zzz-DTTCompatibility-class), 13
- bioMaterialPlateIdentifier,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- bioMaterialPlateRow
(zzz-DTTCompatibility-class), 13
- bioMaterialPlateRow,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- bioSample (zzz-DTTCompatibility-class),
13
- bioSample,BioMaterial_ assnlist-method
(zzz-DTTCompatibility-class), 13
- BioSample-class
(zzz-DTTCompatibility-class), 13

- bioSample_ref
 - (zzz-DTTCompatibility-class), 13
- bioSample_ref,BioMaterial_assnref-method
 - (zzz-DTTCompatibility-class), 13
- BioSample_ref-class
 - (zzz-DTTCompatibility-class), 13
- bioSequence (zzz-DTTCompatibility-class), 13
- bioSequence,BioSequence_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- BioSequence-class
 - (zzz-DTTCompatibility-class), 13
- bioSequence_assnlist
 - (zzz-DTTCompatibility-class), 13
- bioSequence_assnlist,BioSequence_package-method
 - (zzz-DTTCompatibility-class), 13
- BioSequence_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- bioSequence_package
 - (zzz-DTTCompatibility-class), 13
- bioSequence_package,MAGE-ML-method
 - (zzz-DTTCompatibility-class), 13
- BioSequence_package-class
 - (zzz-DTTCompatibility-class), 13
- bioSequence_ref
 - (zzz-DTTCompatibility-class), 13
- bioSequence_ref,BiologicalCharacteristics_assnreflist-method
 - (zzz-DTTCompatibility-class), 13
- bioSequence_ref,ImmobilizedCharacteristics_assnreflist-method
 - (zzz-DTTCompatibility-class), 13
- BioSequence_ref-class
 - (zzz-DTTCompatibility-class), 13
- bioSource (zzz-DTTCompatibility-class), 13
- bioSource,BioMaterial_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- BioSource-class
 - (zzz-DTTCompatibility-class), 13
- bioSource_ref
 - (zzz-DTTCompatibility-class), 13
- bioSource_ref,BioMaterial_assnref-method
 - (zzz-DTTCompatibility-class), 13
- BioSource_ref-class
 - (zzz-DTTCompatibility-class), 13
- bqs_package
 - (zzz-DTTCompatibility-class), 13
- bqs_package,MAGE-ML-method
 - (zzz-DTTCompatibility-class), 13
- BQS_package-class
 - (zzz-DTTCompatibility-class), 13
- category_assn
 - (zzz-DTTCompatibility-class), 13
- category_assn,ExperimentalFactor-method
 - (zzz-DTTCompatibility-class), 13
- Category_assn-class
 - (zzz-DTTCompatibility-class), 13
- channel (zzz-DTTCompatibility-class), 13
- channel,Channel_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- Channel-class
 - (zzz-DTTCompatibility-class), 13
- channel_assnlist
 - (zzz-DTTCompatibility-class), 13
- channel_assnlist,BioAssay_package-method
 - (zzz-DTTCompatibility-class), 13
- Channel_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- channel_assnref
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,DerivedSignal-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,Error-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,ExpectedValue-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,Failed-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,MeasuredSignal-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,PresentAbsent-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,PValue-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,Ratio-method
 - (zzz-DTTCompatibility-class), 13
- channel_assnref,SpecializedQuantitationType-method
 - (zzz-DTTCompatibility-class), 13
- Channel_assnref-class
 - (zzz-DTTCompatibility-class), 13
- channel_ref (zzz-DTTCompatibility-class), 13
- channel_ref,Channel_assnref-method
 - (zzz-DTTCompatibility-class), 13
- channel_ref,Channels_assnreflist-method
 - (zzz-DTTCompatibility-class), 13
- Channel_ref-class
 - (zzz-DTTCompatibility-class), 13
- channels_assnreflist
 - (zzz-DTTCompatibility-class), 13
- channels_assnreflist,DerivedBioAssay-method
 - (zzz-DTTCompatibility-class), 13
- channels_assnreflist,Image-method
 - (zzz-DTTCompatibility-class), 13
- channels_assnreflist,MeasuredBioAssay-method

- (zzz-DTTCompatibility-class), 13
- channels_ assnreflist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- Channels_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- characteristics_ assnlist
(zzz-DTTCompatibility-class), 13
- characteristics_ assnlist,BioSample-method
(zzz-DTTCompatibility-class), 13
- characteristics_ assnlist,BioSource-method
(zzz-DTTCompatibility-class), 13
- characteristics_ assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- Characteristics_ assnlist-class
(zzz-DTTCompatibility-class), 13
- clusterBioAssayData_ assnref
(zzz-DTTCompatibility-class), 13
- clusterBioAssayData_ assnref,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- ClusterBioAssayData_ assnref-class
(zzz-DTTCompatibility-class), 13
- column (zzz-DTTCompatibility-class), 13
- column,FeatureLocation-method
(zzz-DTTCompatibility-class), 13
- column,Zone-method
(zzz-DTTCompatibility-class), 13
- componentCompounds_ assnlist
(zzz-DTTCompatibility-class), 13
- componentCompounds_ assnlist,Compound-method
(zzz-DTTCompatibility-class), 13
- ComponentCompounds_ assnlist-class
(zzz-DTTCompatibility-class), 13
- composite_ assnref
(zzz-DTTCompatibility-class), 13
- composite_ assnref,CompositePosition-method
(zzz-DTTCompatibility-class), 13
- Composite_ assnref-class
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap,CompositeCompositeMap_ assnlist
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap,DesignElementMap_ assnlist
(zzz-DTTCompatibility-class), 13
- CompositeCompositeMap-class
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap_ assnlist
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap_ assnlist,DesignElement_ package
(zzz-DTTCompatibility-class), 13
- CompositeCompositeMap_ assnlist-class
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap_ ref
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap_ ref,CompositeCompositeMaps_ assnlist
(zzz-DTTCompatibility-class), 13
- compositeCompositeMap_ ref,DesignElementMaps_ assnreflist
(zzz-DTTCompatibility-class), 13
- CompositeCompositeMap_ ref-class
(zzz-DTTCompatibility-class), 13
- compositeCompositeMaps_ assnreflist
(zzz-DTTCompatibility-class), 13
- compositeCompositeMaps_ assnreflist,CompositeSequence-met
(zzz-DTTCompatibility-class), 13
- CompositeCompositeMaps_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- compositeGroup
(zzz-DTTCompatibility-class), 13
- compositeGroup,CompositeGroup_ assnlist-method
(zzz-DTTCompatibility-class), 13
- CompositeGroup-class
(zzz-DTTCompatibility-class), 13
- compositeGroup_ assnlist
(zzz-DTTCompatibility-class), 13
- compositeGroup_ assnlist,ArrayDesign_ package-method
(zzz-DTTCompatibility-class), 13
- CompositeGroup_ assnlist-class
(zzz-DTTCompatibility-class), 13
- compositeGroup_ ref
(zzz-DTTCompatibility-class), 13
- compositeGroup_ ref,CompositeGroups_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- CompositeGroup_ ref-class
(zzz-DTTCompatibility-class), 13
- compositeGroups_ assnreflist
(zzz-DTTCompatibility-class), 13
- compositeGroups_ assnreflist,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- compositeGroups_ assnreflist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- CompositeGroups_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- compositePosition
(zzz-DTTCompatibility-class), 13
- compositePosition,CompositePositionSources_ assnlist-method
(zzz-DTTCompatibility-class), 13
- compositePosition,Coordinate_ assn-method
(zzz-DTTCompatibility-class), 13
- CompositePosition-class
(zzz-DTTCompatibility-class), 13
- compositePositionSources_ assnlist
(zzz-DTTCompatibility-class), 13
- compositePositionSources_ assnlist,CompositeCompositeMap-
(zzz-DTTCompatibility-class), 13

- CompositePositionSources_assnlist-class
(zzz-DTTCompatibility-class), 13
- compositeSequence
(zzz-DTTCompatibility-class), 13
- compositeSequence,CompositeSequence_assnlist-method
(zzz-DTTCompatibility-class), 13
- compositeSequence,DesignElement_assnlist-method
(zzz-DTTCompatibility-class), 13
- CompositeSequence-class
(zzz-DTTCompatibility-class), 13
- compositeSequence_assnlist
(zzz-DTTCompatibility-class), 13
- compositeSequence_assnlist,DesignElement_package-method
(zzz-DTTCompatibility-class), 13
- CompositeSequence_assnlist-class
(zzz-DTTCompatibility-class), 13
- compositeSequence_assnref
(zzz-DTTCompatibility-class), 13
- compositeSequence_assnref,CompositeMeasurement_method
(zzz-DTTCompatibility-class), 13
- compositeSequence_assnref,ReporterCompositeMeasurement_method
(zzz-DTTCompatibility-class), 13
- CompositeSequence_assnref-class
(zzz-DTTCompatibility-class), 13
- compositeSequence_ref
(zzz-DTTCompatibility-class), 13
- compositeSequence_ref,Composite_assnref-method
(zzz-DTTCompatibility-class), 13
- compositeSequence_ref,CompositeSequence_assnref-method
(zzz-DTTCompatibility-class), 13
- compositeSequence_ref,CompositeSequences_assnref-method
(zzz-DTTCompatibility-class), 13
- compositeSequence_ref,DesignElement_assnref-method
(zzz-DTTCompatibility-class), 13
- CompositeSequence_ref-class
(zzz-DTTCompatibility-class), 13
- compositeSequenceDimension
(zzz-DTTCompatibility-class), 13
- compositeSequenceDimension,DesignElementDimension_assnref-method
(zzz-DTTCompatibility-class), 13
- CompositeSequenceDimension-class
(zzz-DTTCompatibility-class), 13
- compositeSequenceDimension_ref
(zzz-DTTCompatibility-class), 13
- compositeSequenceDimension_ref,DesignElementDimension_assnref-method
(zzz-DTTCompatibility-class), 13
- CompositeSequenceDimension_ref-class
(zzz-DTTCompatibility-class), 13
- compositeSequences_assnreflist
(zzz-DTTCompatibility-class), 13
- compositeSequences_assnreflist,CompositeGroup_method
(zzz-DTTCompatibility-class), 13
- compositeSequences_assnreflist,CompositeSequenceDimension
(zzz-DTTCompatibility-class), 13
- CompositeSequences_assnreflist-class
(zzz-DTTCompatibility-class), 13
- compound
(zzz-DTTCompatibility-class), 13
- compound_assnlist
(zzz-DTTCompatibility-class), 13
- compound_assnlist,BioMaterial_package-method
(zzz-DTTCompatibility-class), 13
- Compound_assnlist-class
(zzz-DTTCompatibility-class), 13
- compound_assnref
(zzz-DTTCompatibility-class), 13
- CompoundAssnref_method
(zzz-DTTCompatibility-class), 13
- CompoundAssnref-class
(zzz-DTTCompatibility-class), 13
- compound_ref
(zzz-DTTCompatibility-class), 13
- compound_ref,Compound_assnref-method
(zzz-DTTCompatibility-class), 13
- compound_ref,Labels_assnreflist-method
(zzz-DTTCompatibility-class), 13
- Compound_ref-class
(zzz-DTTCompatibility-class), 13
- compoundIndices
(zzz-DTTCompatibility-class), 13
- CompoundIndices_assnlist
(zzz-DTTCompatibility-class), 13
- compoundMeasurement
(zzz-DTTCompatibility-class), 13
- compoundMeasurement,ComponentCompounds_assnlist-method
(zzz-DTTCompatibility-class), 13
- compoundMeasurement,CompoundMeasurements_assnlist-method
(zzz-DTTCompatibility-class), 13
- CompoundMeasurement-class
(zzz-DTTCompatibility-class), 13
- CompoundMeasurementAssnreflist
(zzz-DTTCompatibility-class), 13
- compoundMeasurements_assnlist,Treatment-method
(zzz-DTTCompatibility-class), 13
- CompoundMeasurements_assnlist-class
(zzz-DTTCompatibility-class), 13
- ConcentrationUnit
(zzz-DTTCompatibility-class), 13

- concentrationUnit,Unit _ assn-method
(zzz-DTTCompatibility-class), 13
- ConcentrationUnit-class
(zzz-DTTCompatibility-class), 13
- ConfidenceIndicator _ assnlist-class
(zzz-DTTCompatibility-class), 13
- ConfidenceIndicator _ ref-class
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,Error-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,Failed-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,PValue-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,Ratio-method
(zzz-DTTCompatibility-class), 13
- confidenceIndicators _ assnreflist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- ConfidenceIndicators _ assnreflist-class
(zzz-DTTCompatibility-class), 13
- contact _ assnlist
(zzz-DTTCompatibility-class), 13
- contact _ assnlist,AuditAndSecurity _ package-method
(zzz-DTTCompatibility-class), 13
- Contact _ assnlist-class
(zzz-DTTCompatibility-class), 13
- Contact _ ref-class
(zzz-DTTCompatibility-class), 13
- contacts _ assnreflist
(zzz-DTTCompatibility-class), 13
- contacts _ assnreflist,Database-method
(zzz-DTTCompatibility-class), 13
- Contacts _ assnreflist-class
(zzz-DTTCompatibility-class), 13
- containedFeatures _ assnreflist
(zzz-DTTCompatibility-class), 13
- containedFeatures _ assnreflist,FeatureDimension-method
(zzz-DTTCompatibility-class), 13
- ContainedFeatures _ assnreflist-class
(zzz-DTTCompatibility-class), 13
- control (zzz-ARRCompatibility-class), 12
- control,UserAttribute-method
(zzz-ARRCompatibility-class), 12
- Control-class
(zzz-ARRCompatibility-class), 12
- controlFeatures _ assnreflist
(zzz-DTTCompatibility-class), 13
- controlFeatures _ assnreflist,Feature-method
(zzz-DTTCompatibility-class), 13
- ControlFeatures _ assnreflist-class
(zzz-DTTCompatibility-class), 13
- controlledFeatures _ assnreflist
(zzz-DTTCompatibility-class), 13
- controlledFeatures _ assnreflist,Feature-method
(zzz-DTTCompatibility-class), 13
- ControlledFeatures _ assnreflist-class
(zzz-DTTCompatibility-class), 13
- controlType _ assn
(zzz-DTTCompatibility-class), 13
- controlType _ assn,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- controlType _ assn,Feature-method
(zzz-DTTCompatibility-class), 13
- controlType _ assn,Reporter-method
(zzz-DTTCompatibility-class), 13
- ControlType _ assn-class
(zzz-DTTCompatibility-class), 13
- coordinate _ assn
(zzz-DTTCompatibility-class), 13
- coordinate _ assn,SeqFeatureLocation-method
(zzz-DTTCompatibility-class), 13
- Coordinate _ assn-class
(zzz-DTTCompatibility-class), 13
- createdBy (zzz-ARRCompatibility-class), 12
- createdBy,ArraySetFile-method
(zzz-ARRCompatibility-class), 12
- createdBy,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- createdBy,TemplateFile-method
(zzz-ARRCompatibility-class), 12
- createdDateTime
(zzz-ARRCompatibility-class), 12
- createdDateTime,ArraySetFile-method
(zzz-ARRCompatibility-class), 12
- createdDateTime,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- createdDateTime,TemplateFile-method
(zzz-ARRCompatibility-class), 12
- createdStep (zzz-ARRCompatibility-class), 12
- createdStep,ArraySetFile-method
(zzz-ARRCompatibility-class), 12
- createdStep,PhysicalArray-method

- (zzz-ARRCompatibility-class), 12
- database (zzz-DTTCompatibility-class), 13
- database,Database_assnlist-method
(zzz-DTTCompatibility-class), 13
- Database-class
(zzz-DTTCompatibility-class), 13
- database_assnlist
(zzz-DTTCompatibility-class), 13
- database_assnlist,Description_package-method
(zzz-DTTCompatibility-class), 13
- Database_assnlist-class
(zzz-DTTCompatibility-class), 13
- database_assnref
(zzz-DTTCompatibility-class), 13
- database_assnref,DatabaseEntry-method
(zzz-DTTCompatibility-class), 13
- Database_assnref-class
(zzz-DTTCompatibility-class), 13
- database_ref
(zzz-DTTCompatibility-class), 13
- database_ref,Database_assnref-method
(zzz-DTTCompatibility-class), 13
- Database_ref-class
(zzz-DTTCompatibility-class), 13
- databaseEntry
(zzz-DTTCompatibility-class), 13
- databaseEntry,Accessions_assnlist-method
(zzz-DTTCompatibility-class), 13
- databaseEntry,DatabaseReferences_assnlist-method
(zzz-DTTCompatibility-class), 13
- databaseEntry,ExternalLIMS_assn-method
(zzz-DTTCompatibility-class), 13
- databaseEntry,IdentifierLIMS_assn-method
(zzz-DTTCompatibility-class), 13
- databaseEntry,OntologyReference_assn-method
(zzz-DTTCompatibility-class), 13
- databaseEntry,SequenceDatabases_assnlist-method
(zzz-DTTCompatibility-class), 13
- DatabaseEntry-class
(zzz-DTTCompatibility-class), 13
- databaseReferences_assnlist
(zzz-DTTCompatibility-class), 13
- databaseReferences_assnlist,Description-method
(zzz-DTTCompatibility-class), 13
- DatabaseReferences_assnlist-class
(zzz-DTTCompatibility-class), 13
- dataExternal
(zzz-DTTCompatibility-class), 13
- dataExternal,DataExternal_assn-method
(zzz-DTTCompatibility-class), 13
- DataExternal-class
(zzz-DTTCompatibility-class), 13
- dataExternal_assn
(zzz-DTTCompatibility-class), 13
- dataExternal_assn,BioDataCube-method
(zzz-DTTCompatibility-class), 13
- DataExternal_assn-class
(zzz-DTTCompatibility-class), 13
- dataFormat (zzz-DTTCompatibility-class),
13
- dataFormat,DataExternal-method
(zzz-DTTCompatibility-class), 13
- dataFormatInfoURI
(zzz-DTTCompatibility-class), 13
- dataFormatInfoURI,DataExternal-method
(zzz-DTTCompatibility-class), 13
- dataInternal (zzz-DTTCompatibility-class),
13
- dataInternal,DataInternal_assn-method
(zzz-DTTCompatibility-class), 13
- DataInternal-class
(zzz-DTTCompatibility-class), 13
- dataInternal_assn
(zzz-DTTCompatibility-class), 13
- dataInternal_assn,BioDataCube-method
(zzz-DTTCompatibility-class), 13
- DataInternal_assn-class
(zzz-DTTCompatibility-class), 13
- dataType_assn
(zzz-DTTCompatibility-class), 13
- dataType_assn,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,Error-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,Failed-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,NodeValue-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,Parameter-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,PValue-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,Ratio-method
(zzz-DTTCompatibility-class), 13
- dataType_assn,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- DataType_assn-class
(zzz-DTTCompatibility-class), 13

- datum (*zzz-DTTCompatibility-class*), 13
- datum, Datum_ assn-method (*zzz-DTTCompatibility-class*), 13
- Datum-class (*zzz-DTTCompatibility-class*), 13
- datum_ assn (*zzz-DTTCompatibility-class*), 13
- datum_ assn, QuantitationTypeTuple-method (*zzz-DTTCompatibility-class*), 13
- Datum_ assn-class (*zzz-DTTCompatibility-class*), 13
- defaultValue (*zzz-ARRCompatibility-class*), 12
- defaultValue, UserAttribute-method (*zzz-ARRCompatibility-class*), 12
- defaultValue_ assn (*zzz-DTTCompatibility-class*), 13
- defaultValue_ assn, Parameter-method (*zzz-DTTCompatibility-class*), 13
- DefaultValue_ assn-class (*zzz-DTTCompatibility-class*), 13
- defectType_ assn (*zzz-DTTCompatibility-class*), 13
- defectType_ assn, FeatureDefect-method (*zzz-DTTCompatibility-class*), 13
- defectType_ assn, ZoneDefect-method (*zzz-DTTCompatibility-class*), 13
- DefectType_ assn-class (*zzz-DTTCompatibility-class*), 13
- deltaX (*zzz-DTTCompatibility-class*), 13
- deltaX, PositionDelta-method (*zzz-DTTCompatibility-class*), 13
- deltaY (*zzz-DTTCompatibility-class*), 13
- deltaY, PositionDelta-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay, BioAssay_ assnlist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay, DerivedBioAssay_ assnlist-method (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssay-class (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssay_ assnlist-class (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay_ ref (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay_ ref, BioAssay_ assnref-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay_ ref, BioAssayMapTarget_ assnref-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay_ ref, BioAssays_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay_ ref, SourceBioAssays_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssay_ ref, TopLevelBioAssays_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssay_ ref-class (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData, BioAssayData_ assnlist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData, DerivedBioAssayData_ assnlist-method (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssayData-class (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssayData_ assnlist-class (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ assnreflist (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ assnreflist, DerivedBioAssay-method (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssayData_ assnreflist-class (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ ref (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ ref, BioAssayData_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ ref, BioAssayDataSources_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ ref, ClusterBioAssayData_ assnref-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ ref, DerivedBioAssayData_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayData_ ref, DerivedBioAssayDataTarget_ assnreflist-method (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssayData_ ref-class (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssayDataTarget_ assnref-class (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayMap_ assnreflist (*zzz-DTTCompatibility-class*), 13
- derivedBioAssayMap_ assnreflist, DerivedBioAssay-method (*zzz-DTTCompatibility-class*), 13
- DerivedBioAssayMap_ assnreflist-class (*zzz-DTTCompatibility-class*), 13
- derivedSignal (*zzz-DTTCompatibility-class*), 13
- derivedSignal, QuantitationType_ assnlist-method (*zzz-DTTCompatibility-class*), 13
- MethodSignal-class (*zzz-DTTCompatibility-class*), 13
- derivedSignal_ ref

- (zzz-DTTCompatibility-class), 13
- derivedSignal_ref,QuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- derivedSignal_ref,QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class), 13
- derivedSignal_ref,SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- derivedSignal_ref,TargetQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- DerivedSignal_ref-class (zzz-DTTCompatibility-class), 13
- description (zzz-DTTCompatibility-class), 13
- description,Descriptions_assnlist-method (zzz-DTTCompatibility-class), 13
- description,NormalizationDescription_assn-method (zzz-DTTCompatibility-class), 13
- description,OntologyEntry-method (zzz-DTTCompatibility-class), 13
- description,QualityControlDescription_assn-method (zzz-DTTCompatibility-class), 13
- description,ReplicateDescription_assn-method (zzz-DTTCompatibility-class), 13
- Description-class (zzz-DTTCompatibility-class), 13
- description_package (zzz-DTTCompatibility-class), 13
- description_package,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- Description_package-class (zzz-DTTCompatibility-class), 13
- descriptions_assnlist (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Array-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ArrayDesign-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ArrayGroup-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ArrayManufacture-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Audit-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BibliographicReference-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayCreation-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayDataCluster-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayDimension-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayMap-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioAssayTreatment-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioSample-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioSequence-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,BioSource-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Channel-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeCompositeMap-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeGroup-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeSequence-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,CompositeSequenceDimension-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Compound-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Database-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,DerivedBioAssay-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,DerivedBioAssayData-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,DerivedSignal-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Description-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Error-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ExpectedValue-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Experiment-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ExperimentalFactor-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,ExperimentDesign-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,FactorValue-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Failed-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,Feature-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,FeatureDimension-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,FeatureExtraction-method (zzz-DTTCompatibility-class), 13
- descriptions_assnlist,FeatureGroup-method (zzz-DTTCompatibility-class), 13

- (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,FeatureReporterMap-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Fiducial-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Hardware-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,HardwareApplication-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Hybridization-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Image-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ImageAcquisition-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,LabeledExtract-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ManufactureLIMS-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ManufactureLIMSBiomaterials-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,MeasuredBioAssay-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,MeasuredBioAssayData-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,MeasuredSignal-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Node-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,NodeContents-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Organization-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Parameter-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Person-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,PhysicalBioAssay-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,PresentAbsent-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Protocol-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ProtocolApplication-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,PValue-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,QuantitationTypeDimension-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,QuantitationTypeMap-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Ratio-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Reporter-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ReporterCompositeMap-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ReporterDimension-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,ReporterGroup-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Security-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,SecurityGroup-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,SeqFeature-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Software-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,SoftwareApplication-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,SpecializedQuantitationType-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Transformation-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Treatment-method (zzz-DTTCompatibility-class), 13
- descriptions__assnlist,Zone-method (zzz-DTTCompatibility-class), 13
- Descriptions__assnlist-class (zzz-DTTCompatibility-class), 13
- DesignElement__assnlist-class (zzz-DTTCompatibility-class), 13
- designElement__assnref (zzz-DTTCompatibility-class), 13
- designElement__assnref,DesignElementTuple-method (zzz-DTTCompatibility-class), 13
- DesignElement__assnref-class (zzz-DTTCompatibility-class), 13
- designElement__package (zzz-DTTCompatibility-class), 13
- designElement__package,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- DesignElement__package-class (zzz-DTTCompatibility-class), 13
- DesignElement__ref-class (zzz-DTTCompatibility-class), 13
- designElementDimension__assnlist (zzz-DTTCompatibility-class), 13
- designElementDimension__assnlist,BioAssayData__package-method (zzz-DTTCompatibility-class), 13

- (zzz-DTTCompatibility-class), 13
- DesignElementDimension_ assnlist-class
(zzz-DTTCompatibility-class), 13
- designElementDimension_ assnref
(zzz-DTTCompatibility-class), 13
- designElementDimension_ assnref,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- designElementDimension_ assnref,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- designElementDimension_ assnref,NodeContents-method
(zzz-DTTCompatibility-class), 13
- DesignElementDimension_ assnref-class
(zzz-DTTCompatibility-class), 13
- DesignElementDimension_ ref-class
(zzz-DTTCompatibility-class), 13
- DesignElementMap_ assnlist-class
(zzz-DTTCompatibility-class), 13
- DesignElementMap_ ref-class
(zzz-DTTCompatibility-class), 13
- designElementMapping
(zzz-DTTCompatibility-class), 13
- designElementMapping,DesignElementMapping_ assn-method
(zzz-DTTCompatibility-class), 13
- DesignElementMapping-class
(zzz-DTTCompatibility-class), 13
- designElementMapping_ assn
(zzz-DTTCompatibility-class), 13
- designElementMapping_ assn,Transformation-method
(zzz-DTTCompatibility-class), 13
- DesignElementMapping_ assn-class
(zzz-DTTCompatibility-class), 13
- designElementMaps_ assnreflist
(zzz-DTTCompatibility-class), 13
- designElementMaps_ assnreflist,DesignElementMapping-method
(zzz-DTTCompatibility-class), 13
- DesignElementMaps_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- designElementTuple
(zzz-DTTCompatibility-class), 13
- designElementTuple,DesignElementTuples_ assnlist-method
(zzz-DTTCompatibility-class), 13
- DesignElementTuple-class
(zzz-DTTCompatibility-class), 13
- designElementTuples_ assnlist
(zzz-DTTCompatibility-class), 13
- designElementTuples_ assnlist,BioAssayTuple-method
(zzz-DTTCompatibility-class), 13
- DesignElementTuples_ assnlist-class
(zzz-DTTCompatibility-class), 13
- designProviders_ assnreflist
(zzz-DTTCompatibility-class), 13
- designProviders_ assnreflist,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- designProviders_ assnreflist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- DesignProviders_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- DesignProviders_ assnreflist,Processors), 2
- directory,NetAffxResource-method
(NetAffxResource-class), 5
- distanceUnit
(zzz-DTTCompatibility-class), 13
- distanceUnit,DistanceUnit_ assn-method
(zzz-DTTCompatibility-class), 13
- distanceUnit,Unit_ assn-method
(zzz-DTTCompatibility-class), 13
- DistanceUnit-class
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,Fiducial-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,Position-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,PositionDelta-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,Zone-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,ZoneGroup-method
(zzz-DTTCompatibility-class), 13
- distanceUnit_ assn,ZoneLayout-method
(zzz-DTTCompatibility-class), 13
- DistanceUnit_ assn-class
(zzz-DTTCompatibility-class), 13
- DTTCompatibility, 13
- DTTCompatibility-class
(AffyCompatible-class), 3
- editor (zzz-DTTCompatibility-class), 13
- editor,BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- email (zzz-DTTCompatibility-class), 13
- email,Organization-method
(zzz-DTTCompatibility-class), 13
- email,Person-method
(zzz-DTTCompatibility-class), 13
- error (zzz-DTTCompatibility-class), 13
- error,ConfidenceIndicator_ assnlist-method
(zzz-DTTCompatibility-class), 13
- error,QuantitationType_ assnlist-method
(zzz-DTTCompatibility-class), 13

- Error-class (zzz-DTTCompatibility-class), 13
- error_ref (zzz-DTTCompatibility-class), 13
- error_ref,ConfidenceIndicators_assnreflist-method (zzz-DTTCompatibility-class), 13
- error_ref,QuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- error_ref,QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class), 13
- error_ref,SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- error_ref,TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- Error_ref-class (zzz-DTTCompatibility-class), 13
- expectedValue (zzz-DTTCompatibility-class), 13
- expectedValue,ConfidenceIndicator_assnlist-method (zzz-DTTCompatibility-class), 13
- expectedValue,QuantitationType_assnlist-method (zzz-DTTCompatibility-class), 13
- ExpectedValue-class (zzz-DTTCompatibility-class), 13
- expectedValue_ref (zzz-DTTCompatibility-class), 13
- expectedValue_ref,ConfidenceIndicators_assnreflist-method (zzz-DTTCompatibility-class), 13
- expectedValue_ref,QuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- expectedValue_ref,QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class), 13
- expectedValue_ref,SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- expectedValue_ref,TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- ExpectedValue_ref-class (zzz-DTTCompatibility-class), 13
- experiment (zzz-DTTCompatibility-class), 13
- experiment,Experiment_assnlist-method (zzz-DTTCompatibility-class), 13
- Experiment-class (zzz-DTTCompatibility-class), 13
- experiment_assnlist (zzz-DTTCompatibility-class), 13
- experiment_assnlist,Experiment_package-method (zzz-DTTCompatibility-class), 13
- Experiment_assnlist-class (zzz-DTTCompatibility-class), 13
- experiment_package (zzz-DTTCompatibility-class), 13
- experiment_package,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- Experiment_package-class (zzz-DTTCompatibility-class), 13
- Experiment_package-class (zzz-DTTCompatibility-class), 13
- experimentalFactor (zzz-DTTCompatibility-class), 13
- experimentalFactor,ExperimentalFactors_assnlist-method (zzz-DTTCompatibility-class), 13
- ExperimentalFactor-class (zzz-DTTCompatibility-class), 13
- ExperimentalFactor_assnref-class (zzz-DTTCompatibility-class), 13
- experimentalFactor_ref (zzz-DTTCompatibility-class), 13
- experimentalFactor_ref,ExperimentalFactor_assnref-method (zzz-DTTCompatibility-class), 13
- ExperimentalFactor_ref-class (zzz-DTTCompatibility-class), 13
- experimentalFactors_assnlist (zzz-DTTCompatibility-class), 13
- experimentalFactors_assnlist,ExperimentDesign-method (zzz-DTTCompatibility-class), 13
- ExperimentalFactors_assnlist-class (zzz-DTTCompatibility-class), 13
- experimentDesign (zzz-DTTCompatibility-class), 13
- experimentDesign,ExperimentDesigns_assnlist-method (zzz-DTTCompatibility-class), 13
- ExperimentDesign-class (zzz-DTTCompatibility-class), 13
- ExperimentDesigns_assnlist (zzz-DTTCompatibility-class), 13
- ExperimentDesigns_assnlist,Experiment-method (zzz-DTTCompatibility-class), 13
- ExperimentDesigns_assnlist-class (zzz-DTTCompatibility-class), 13
- exportedFromDB (zzz-DTTCompatibility-class), 13
- exportedFromDB,ExternalReference-method (zzz-DTTCompatibility-class), 13
- exportedFromServer (zzz-DTTCompatibility-class), 13
- exportedFromServer,ExternalReference-method (zzz-DTTCompatibility-class), 13
- exportID (zzz-DTTCompatibility-class), 13
- exportID,ExternalReference-method (zzz-DTTCompatibility-class), 13
- exportName (zzz-DTTCompatibility-class), 13
- exportName,ExternalReference-method (zzz-DTTCompatibility-class), 13
- externalLIMS_assn (zzz-DTTCompatibility-class), 13

- externalLIMS_assn,Compound-method
(zzz-DTTCompatibility-class), 13
- ExternalLIMS_assn-class
(zzz-DTTCompatibility-class), 13
- externalReference
(zzz-DTTCompatibility-class), 13
- externalReference,ExternalReference_assn-method
(zzz-DTTCompatibility-class), 13
- externalReference_assn
(zzz-DTTCompatibility-class), 13
- externalReference_assn,Description-method
(zzz-DTTCompatibility-class), 13
- ExternalReference_assn-class
(zzz-DTTCompatibility-class), 13

- factorValue (zzz-DTTCompatibility-class),
13
- factorValue,FactorValues_assnlist-method
(zzz-DTTCompatibility-class), 13
- FactorValue-class
(zzz-DTTCompatibility-class), 13
- factorValue_ref
(zzz-DTTCompatibility-class), 13
- factorValue_ref,BioAssayFactorValues_assnreflist-method
(zzz-DTTCompatibility-class), 13
- FactorValue_ref-class
(zzz-DTTCompatibility-class), 13
- factorValues_assnlist
(zzz-DTTCompatibility-class), 13
- factorValues_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class), 13
- FactorValues_assnlist-class
(zzz-DTTCompatibility-class), 13
- failed (zzz-DTTCompatibility-class), 13
- failed,QuantitationType_assnlist-method
(zzz-DTTCompatibility-class), 13
- Failed-class (zzz-DTTCompatibility-class),
13
- failed_ref (zzz-DTTCompatibility-class),
13
- failed_ref,QuantitationType_assnref-method
(zzz-DTTCompatibility-class), 13
- failed_ref,QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class), 13
- failed_ref,SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class), 13
- failed_ref,TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class), 13
- Failed_ref-class
(zzz-DTTCompatibility-class), 13
- failTypes_assnlist
(zzz-DTTCompatibility-class), 13
- failTypes_assnlist,Reporter-method
(zzz-DTTCompatibility-class), 13
- FailTypes_assnlist-class
(zzz-DTTCompatibility-class), 13
- fax (zzz-DTTCompatibility-class), 13
- fax,Organization-method
(zzz-DTTCompatibility-class), 13
- fax,Person-method
(zzz-DTTCompatibility-class), 13
- feature (zzz-DTTCompatibility-class), 13
- feature,DesignElement_assnlist-method
(zzz-DTTCompatibility-class), 13
- feature,Features_assnlist-method
(zzz-DTTCompatibility-class), 13
- Feature-class
(zzz-DTTCompatibility-class), 13
- feature_assnref
(zzz-DTTCompatibility-class), 13
- feature_assnref,FeatureDefect-method
(zzz-DTTCompatibility-class), 13
- feature_assnref,FeatureInformation-method
(zzz-DTTCompatibility-class), 13
- feature_assnref,ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- feature_assnref,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- Feature_assnref-class
(zzz-DTTCompatibility-class), 13
- feature_ref (zzz-DTTCompatibility-class),
13
- feature_ref,ContainedFeatures_assnreflist-method
(zzz-DTTCompatibility-class), 13
- feature_ref,ControlFeatures_assnreflist-method
(zzz-DTTCompatibility-class), 13
- feature_ref,ControlledFeatures_assnreflist-method
(zzz-DTTCompatibility-class), 13
- feature_ref,DesignElement_assnref-method
(zzz-DTTCompatibility-class), 13
- feature_ref,Feature_assnref-method
(zzz-DTTCompatibility-class), 13
- Feature_ref-class
(zzz-DTTCompatibility-class), 13
- featureDefect
(zzz-DTTCompatibility-class), 13
- featureDefect,FeatureDefects_assnlist-method
(zzz-DTTCompatibility-class), 13
- FeatureDefect-class
(zzz-DTTCompatibility-class), 13
- featureDefects_assnlist
(zzz-DTTCompatibility-class), 13
- featureDefects_assnlist,ArrayManufactureDeviation-method
(zzz-DTTCompatibility-class), 13

- FeatureDefects _ assnlist-class
(zzz-DTTCompatibility-class), 13
- featureDimension
(zzz-DTTCompatibility-class), 13
- featureDimension,DesignElementDimension _ assnlist-class
(zzz-DTTCompatibility-class), 13
- FeatureDimension-class
(zzz-DTTCompatibility-class), 13
- featureDimension _ ref
(zzz-DTTCompatibility-class), 13
- featureDimension _ ref,DesignElementDimension _ assnlist-class
(zzz-DTTCompatibility-class), 13
- FeatureDimension _ ref-class
(zzz-DTTCompatibility-class), 13
- featureExtraction
(zzz-DTTCompatibility-class), 13
- featureExtraction,FeatureExtraction _ assn-method
(zzz-DTTCompatibility-class), 13
- FeatureExtraction-class
(zzz-DTTCompatibility-class), 13
- featureExtraction _ assn
(zzz-DTTCompatibility-class), 13
- featureExtraction _ assn,MeasuredBioAssay-method
(zzz-DTTCompatibility-class), 13
- FeatureExtraction _ assn-class
(zzz-DTTCompatibility-class), 13
- featureGroup
(zzz-DTTCompatibility-class), 13
- featureGroup,FeatureGroups _ assnlist-method
(zzz-DTTCompatibility-class), 13
- FeatureGroup-class
(zzz-DTTCompatibility-class), 13
- FeatureGroup _ assnref-class
(zzz-DTTCompatibility-class), 13
- featureGroup _ ref
(zzz-DTTCompatibility-class), 13
- featureGroup _ ref,FeatureGroup _ assnref-method
(zzz-DTTCompatibility-class), 13
- FeatureGroup _ ref-class
(zzz-DTTCompatibility-class), 13
- featureGroups _ assnlist
(zzz-DTTCompatibility-class), 13
- featureGroups _ assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- featureGroups _ assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- FeatureGroups _ assnlist-class
(zzz-DTTCompatibility-class), 13
- featureHeight
(zzz-DTTCompatibility-class), 13
- featureHeight,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- featureInformation
(zzz-DTTCompatibility-class), 13
- featureInformation,FeatureInformationSources _ assnlist-method
(zzz-DTTCompatibility-class), 13
- FeatureInformation-class
(zzz-DTTCompatibility-class), 13
- featureInformationSources _ assnlist
(zzz-DTTCompatibility-class), 13
- featureInformationSources _ assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- FeatureInformationSources _ assnlist-class
(zzz-DTTCompatibility-class), 13
- featureLength
(zzz-DTTCompatibility-class), 13
- featureLength,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- featureLIMSSs _ assnlist
(zzz-DTTCompatibility-class), 13
- featureLIMSSs _ assnlist,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- FeatureLIMSSs _ assnlist-class
(zzz-DTTCompatibility-class), 13
- featureLocation
(zzz-DTTCompatibility-class), 13
- featureLocation,FeatureLocation _ assn-method
(zzz-DTTCompatibility-class), 13
- FeatureLocation-class
(zzz-DTTCompatibility-class), 13
- featureLocation _ assn
(zzz-DTTCompatibility-class), 13
- featureLocation _ assn,Feature-method
(zzz-DTTCompatibility-class), 13
- FeatureLocation _ assn-class
(zzz-DTTCompatibility-class), 13
- featureReporterMap
(zzz-DTTCompatibility-class), 13
- featureReporterMap,DesignElementMap _ assnlist-method
(zzz-DTTCompatibility-class), 13
- featureReporterMap,FeatureReporterMap _ assnlist-method
(zzz-DTTCompatibility-class), 13
- FeatureReporterMap-class
(zzz-DTTCompatibility-class), 13
- featureReporterMap _ assnlist
(zzz-DTTCompatibility-class), 13
- featureReporterMap _ assnlist,DesignElement _ package-method
(zzz-DTTCompatibility-class), 13
- FeatureReporterMap _ assnlist-class
(zzz-DTTCompatibility-class), 13
- featureReporterMap _ ref
(zzz-DTTCompatibility-class), 13
- featureReporterMap _ ref,DesignElementMaps _ assnreflist-method
(zzz-DTTCompatibility-class), 13

- featureReporterMap_ref,FeatureReporterMaps_assnreflist-method (zzz-DTTCompatibility-class), 13
- FeatureReporterMap_ref-class (zzz-DTTCompatibility-class), 13
- featureReporterMaps_assnreflist (zzz-DTTCompatibility-class), 13
- featureReporterMaps_assnreflist,Reporter-method (zzz-DTTCompatibility-class), 13
- FeatureReporterMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- features_assnlist (zzz-DTTCompatibility-class), 13
- features_assnlist,FeatureGroup-method (zzz-DTTCompatibility-class), 13
- Features_assnlist-class (zzz-DTTCompatibility-class), 13
- featureShape_assn (zzz-DTTCompatibility-class), 13
- featureShape_assn,FeatureGroup-method (zzz-DTTCompatibility-class), 13
- FeatureShape_assn-class (zzz-DTTCompatibility-class), 13
- featureWidth (zzz-DTTCompatibility-class), 13
- featureWidth,FeatureGroup-method (zzz-DTTCompatibility-class), 13
- fiducial (zzz-DTTCompatibility-class), 13
- fiducial,Fiducials_assnlist-method (zzz-DTTCompatibility-class), 13
- Fiducial-class (zzz-DTTCompatibility-class), 13
- fiducials_assnlist (zzz-DTTCompatibility-class), 13
- fiducials_assnlist,ArrayGroup-method (zzz-DTTCompatibility-class), 13
- Fiducials_assnlist-class (zzz-DTTCompatibility-class), 13
- fiducialType_assn (zzz-DTTCompatibility-class), 13
- fiducialType_assn,Fiducial-method (zzz-DTTCompatibility-class), 13
- FiducialType_assn-class (zzz-DTTCompatibility-class), 13
- filenameURI (zzz-DTTCompatibility-class), 13
- filenameURI,DataExternal-method (zzz-DTTCompatibility-class), 13
- firstName (zzz-DTTCompatibility-class), 13
- firstName,Person-method (zzz-DTTCompatibility-class), 13
- format_assn (zzz-DTTCompatibility-class), 13
- format_assn,Image-method (zzz-DTTCompatibility-class), 13
- Format_assn-class (zzz-DTTCompatibility-class), 13
- guid (zzz-ARRCompatibility-class), 12
- guid,ArraySetFile-method (zzz-ARRCompatibility-class), 12
- guid,PhysicalArray-method (zzz-ARRCompatibility-class), 12
- guid,TemplateFile-method (zzz-ARRCompatibility-class), 12
- hardware (zzz-DTTCompatibility-class), 13
- hardware,Hardware_assnlist-method (zzz-DTTCompatibility-class), 13
- Hardware-class (zzz-DTTCompatibility-class), 13
- hardware_assnlist (zzz-DTTCompatibility-class), 13
- hardware_assnlist,Protocol_package-method (zzz-DTTCompatibility-class), 13
- Hardware_assnlist-class (zzz-DTTCompatibility-class), 13
- hardware_assnref (zzz-DTTCompatibility-class), 13
- hardware_assnref,HardwareApplication-method (zzz-DTTCompatibility-class), 13
- hardware_assnref,Software-method (zzz-DTTCompatibility-class), 13
- Hardware_assnref-class (zzz-DTTCompatibility-class), 13
- hardware_ref (zzz-DTTCompatibility-class), 13
- hardware_ref,Hardware_assnref-method (zzz-DTTCompatibility-class), 13
- hardware_ref,Hardwares_assnreflist-method (zzz-DTTCompatibility-class), 13
- Hardware_ref-class (zzz-DTTCompatibility-class), 13
- hardwareApplication (zzz-DTTCompatibility-class), 13
- hardwareApplication,HardwareApplications_assnlist-method (zzz-DTTCompatibility-class), 13
- HardwareApplication-class (zzz-DTTCompatibility-class), 13
- hardwareApplications_assnlist (zzz-DTTCompatibility-class), 13
- hardwareApplications_assnlist,ProtocolApplication-method (zzz-DTTCompatibility-class), 13
- HardwareApplications_assnlist-class (zzz-DTTCompatibility-class), 13

- hardwareManufacturers_assnreflist
(zzz-DTTCompatibility-class), 13
- hardwareManufacturers_assnreflist,Hardware-method
(zzz-DTTCompatibility-class), 13
- HardwareManufacturers_assnreflist-class
(zzz-DTTCompatibility-class), 13
- hardwares_assnreflist
(zzz-DTTCompatibility-class), 13
- hardwares_assnreflist,Protocol-method
(zzz-DTTCompatibility-class), 13
- Hardwares_assnreflist-class
(zzz-DTTCompatibility-class), 13
- higherLevelAnalysis_package
(zzz-DTTCompatibility-class), 13
- higherLevelAnalysis_package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- HigherLevelAnalysis_package-class
(zzz-DTTCompatibility-class), 13
- hybridization
(zzz-DTTCompatibility-class), 13
- hybridization,BioAssayCreation_assn-method
(zzz-DTTCompatibility-class), 13
- Hybridization-class
(zzz-DTTCompatibility-class), 13

- identifier (zzz-DTTCompatibility-class), 13
- identifier,Array-method
(zzz-DTTCompatibility-class), 13
- identifier,Array_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- identifier,ArrayDesign_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- identifier,ArrayGroup_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- identifier,ArrayManufacture_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayData_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayDataCluster_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayDimension-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayDimension_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayMap_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13
- identifier,BioMaterial_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioSample-method
(zzz-DTTCompatibility-class), 13
- identifier,BioSample_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioSequence-method
(zzz-DTTCompatibility-class), 13
- identifier,BioSequence_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,BioSource-method
(zzz-DTTCompatibility-class), 13
- identifier,BioSource_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Channel-method
(zzz-DTTCompatibility-class), 13
- identifier,Channel_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeCompositeMap-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeCompositeMap_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeGroup-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeGroup_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeSequence_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class), 13
- identifier,CompositeSequenceDimension_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Compound-method
(zzz-DTTCompatibility-class), 13
- identifier,Compound_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ConfidenceIndicator_ref-method
(zzz-DTTCompatibility-class), 13
- identifier>Contact_ref-method
(zzz-DTTCompatibility-class), 13

- identifier,Database-method
(zzz-DTTCompatibility-class), 13
- identifier,Database_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,DerivedBioAssay-method
(zzz-DTTCompatibility-class), 13
- identifier,DerivedBioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- identifier,DerivedBioAssayData_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- identifier,DerivedSignal_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,DesignElement_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,DesignElementDimension_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,DesignElementMap_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Error-method
(zzz-DTTCompatibility-class), 13
- identifier,Error_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- identifier,ExpectedValue_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Experiment-method
(zzz-DTTCompatibility-class), 13
- identifier,ExperimentalFactor-method
(zzz-DTTCompatibility-class), 13
- identifier,ExperimentalFactor_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,FactorValue-method
(zzz-DTTCompatibility-class), 13
- identifier,FactorValue_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Failed-method
(zzz-DTTCompatibility-class), 13
- identifier,Failed_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Feature-method
(zzz-DTTCompatibility-class), 13
- identifier,Feature_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureDimension-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureDimension_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureExtraction-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureGroup_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- identifier,FeatureReporterMap_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Hardware-method
(zzz-DTTCompatibility-class), 13
- identifier,Hardware_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Hybridization-method
(zzz-DTTCompatibility-class), 13
- identifier,Image-method
(zzz-DTTCompatibility-class), 13
- identifier,Image_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- identifier,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- identifier,LabeledExtract_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- identifier,MeasuredBioAssay-method
(zzz-DTTCompatibility-class), 13
- identifier,MeasuredBioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- identifier,MeasuredBioAssayData_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- identifier,MeasuredSignal_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Organization-method
(zzz-DTTCompatibility-class), 13
- identifier,Organization_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Parameter-method
(zzz-DTTCompatibility-class), 13
- identifier,Parameter_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Person-method
(zzz-DTTCompatibility-class), 13
- identifier,Person_ref-method
(zzz-DTTCompatibility-class), 13

- identifier,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- identifier,PhysicalArrayDesign_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- identifier,PhysicalBioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- identifier,PresentAbsent_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Protocol-method
(zzz-DTTCompatibility-class), 13
- identifier,Protocol_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,PValue-method
(zzz-DTTCompatibility-class), 13
- identifier,PValue_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,QuantitationType_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class), 13
- identifier,QuantitationTypeDimension_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,QuantitationTypeMap-method
(zzz-DTTCompatibility-class), 13
- identifier,QuantitationTypeMap_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Ratio-method
(zzz-DTTCompatibility-class), 13
- identifier,Ratio_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Reporter-method
(zzz-DTTCompatibility-class), 13
- identifier,Reporter_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ReporterCompositeMap-method
(zzz-DTTCompatibility-class), 13
- identifier,ReporterCompositeMap_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ReporterDimension-method
(zzz-DTTCompatibility-class), 13
- identifier,ReporterDimension_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,ReporterGroup-method
(zzz-DTTCompatibility-class), 13
- identifier,ReporterGroup_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Security-method
(zzz-DTTCompatibility-class), 13
- identifier,Security_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,SecurityGroup-method
(zzz-DTTCompatibility-class), 13
- identifier,SecurityGroup_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Software-method
(zzz-DTTCompatibility-class), 13
- identifier,Software_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- identifier,SpecializedQuantitationType_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,StandardQuantitationType_ref-method
(zzz-DTTCompatibility-class), 13
- identifier,Transformation-method
(zzz-DTTCompatibility-class), 13
- identifier,Treatment-method
(zzz-DTTCompatibility-class), 13
- identifier,Zone-method
(zzz-DTTCompatibility-class), 13
- identifier,Zone_ref-method
(zzz-DTTCompatibility-class), 13
- identifierLIMS_assn
(zzz-DTTCompatibility-class), 13
- identifierLIMS_assn,ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- identifierLIMS_assn,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- IdentifierLIMS_assn-class
(zzz-DTTCompatibility-class), 13
- Image-class (zzz-DTTCompatibility-class),
13
- image_ref (zzz-DTTCompatibility-class),
13
- image_ref,Images_assnreflist-method
(zzz-DTTCompatibility-class), 13
- Image_ref-class
(zzz-DTTCompatibility-class), 13
- imageAcquisition
(zzz-DTTCompatibility-class), 13
- imageAcquisition,BioAssayTreatments_assnlist-method
(zzz-DTTCompatibility-class), 13
- ImageAcquisition-class
(zzz-DTTCompatibility-class), 13
- images_assnreflist
(zzz-DTTCompatibility-class), 13
- images_assnreflist,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- Images_assnreflist-class
(zzz-DTTCompatibility-class), 13

- immobilizedCharacteristics_assnreflist
(zzz-DTTCompatibility-class), 13
- immobilizedCharacteristics_assnreflist,Reporter-method
(zzz-DTTCompatibility-class), 13
- ImmobilizedCharacteristics_assnreflist-class
(zzz-DTTCompatibility-class), 13
- information_assnref
(zzz-DTTCompatibility-class), 13
- information_assnref,Array-method
(zzz-DTTCompatibility-class), 13
- Information_assnref-class
(zzz-DTTCompatibility-class), 13
- Internal-class, 3
- isApproximateLength
(zzz-DTTCompatibility-class), 13
- isApproximateLength,BioSequence-method
(zzz-DTTCompatibility-class), 13
- isBackground
(zzz-DTTCompatibility-class), 13
- isBackground,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- isBackground,Error-method
(zzz-DTTCompatibility-class), 13
- isBackground,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- isBackground,Failed-method
(zzz-DTTCompatibility-class), 13
- isBackground,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- isBackground,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- isBackground,PValue-method
(zzz-DTTCompatibility-class), 13
- isBackground,Ratio-method
(zzz-DTTCompatibility-class), 13
- isBackground,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- isCircular (zzz-DTTCompatibility-class), 13
- isCircular,BioSequence-method
(zzz-DTTCompatibility-class), 13
- isSolvent (zzz-DTTCompatibility-class), 13
- isSolvent,Compound-method
(zzz-DTTCompatibility-class), 13
- issue (zzz-DTTCompatibility-class), 13
- issue,BibliographicReference-method
(zzz-DTTCompatibility-class), 13

- kindCV (zzz-DTTCompatibility-class), 13
- kindCV,Measurement-method
(zzz-DTTCompatibility-class), 13

- labeledExtract
(zzz-DTTCompatibility-class), 13
- labeledExtract,BioMaterial_assnlist-method
(zzz-DTTCompatibility-class), 13
- labeledExtract-class
(zzz-DTTCompatibility-class), 13
- labeledExtract_ref
(zzz-DTTCompatibility-class), 13
- labeledExtract_ref,BioMaterial_assnref-method
(zzz-DTTCompatibility-class), 13
- LabeledExtract_ref-class
(zzz-DTTCompatibility-class), 13
- labels_assnreflist
(zzz-DTTCompatibility-class), 13
- labels_assnreflist,Channel-method
(zzz-DTTCompatibility-class), 13
- labels_assnreflist,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- Labels_assnreflist-class
(zzz-DTTCompatibility-class), 13
- lastName (zzz-DTTCompatibility-class), 13
- lastName,Person-method
(zzz-DTTCompatibility-class), 13
- length,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- length,BioSequence-method
(zzz-DTTCompatibility-class), 13
- libraryPackageName
(zzz-ARRCompatibility-class), 12
- libraryPackageName,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- list, 4
- lowerRightX (zzz-DTTCompatibility-class),
13
- lowerRightX,Zone-method
(zzz-DTTCompatibility-class), 13
- lowerRightY (zzz-DTTCompatibility-class),
13
- lowerRightY,Zone-method
(zzz-DTTCompatibility-class), 13

- MAGE-ML-class
(zzz-DTTCompatibility-class), 13
- make (zzz-DTTCompatibility-class), 13
- make,Hardware-method
(zzz-DTTCompatibility-class), 13
- manufactureLIMS
(zzz-DTTCompatibility-class), 13
- manufactureLIMS,FeatureLIMSs_assnlist-method
(zzz-DTTCompatibility-class), 13
- ManufactureLIMS-class
(zzz-DTTCompatibility-class), 13
- manufactureLIMSBiomaterial
(zzz-DTTCompatibility-class), 13

- (zzz-DTTCompatibility-class), 13
- measurement,ActionMeasurement_assn-method
(zzz-DTTCompatibility-class), 13
- measurement,DefaultValue_assn-method
(zzz-DTTCompatibility-class), 13
- measurement,Measurement_assn-method
(zzz-DTTCompatibility-class), 13
- Measurement-class
(zzz-DTTCompatibility-class), 13
- measurement_assn
(zzz-DTTCompatibility-class), 13
- measurement_assn,BioMaterialMeasurement-method
(zzz-DTTCompatibility-class), 13
- measurement_assn,CompoundMeasurement-method
(zzz-DTTCompatibility-class), 13
- measurement_assn,FactorValue-method
(zzz-DTTCompatibility-class), 13
- Measurement_assn-class
(zzz-DTTCompatibility-class), 13
- measurement_package
(zzz-DTTCompatibility-class), 13
- measurement_package,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- Measurement_package-class
(zzz-DTTCompatibility-class), 13
- mediaCol (zzz-ARRCompatibility-class), 12
- mediaCol,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- mediaFileGUID
(zzz-ARRCompatibility-class), 12
- mediaFileGUID,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- mediaFileName
(zzz-ARRCompatibility-class), 12
- mediaFileName,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- mediaRow (zzz-ARRCompatibility-class),
12
- mediaRow,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- mediaType (zzz-ARRCompatibility-class),
12
- mediaType,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- members_assnreflist
(zzz-DTTCompatibility-class), 13
- members_assnreflist,SecurityGroup-method
(zzz-DTTCompatibility-class), 13
- Members_assnreflist-class
(zzz-DTTCompatibility-class), 13
- midInitials (zzz-DTTCompatibility-class),
13
- midInitials,Person-method
(zzz-DTTCompatibility-class), 13
- mismatchInformation
(zzz-DTTCompatibility-class), 13
- mismatchInformation,MismatchInformation_assnlist-method
(zzz-DTTCompatibility-class), 13
- MismatchInformation-class
(zzz-DTTCompatibility-class), 13
- mismatchInformation_assnlist
(zzz-DTTCompatibility-class), 13
- mismatchInformation_assnlist,CompositePosition-method
(zzz-DTTCompatibility-class), 13
- mismatchInformation_assnlist,FeatureInformation-method
(zzz-DTTCompatibility-class), 13
- mismatchInformation_assnlist,ReporterPosition-method
(zzz-DTTCompatibility-class), 13
- MismatchInformation_assnlist-class
(zzz-DTTCompatibility-class), 13
- model (zzz-DTTCompatibility-class), 13
- model,Hardware-method
(zzz-DTTCompatibility-class), 13
- name (zzz-DTTCompatibility-class), 13
- name,Array-method
(zzz-DTTCompatibility-class), 13
- name,Array_ref-method
(zzz-DTTCompatibility-class), 13
- name,ArrayAttribute-method
(zzz-ARRCompatibility-class), 12
- name,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- name,ArrayDesign_ref-method
(zzz-DTTCompatibility-class), 13
- name,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- name,ArrayGroup_ref-method
(zzz-DTTCompatibility-class), 13
- name,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- name,ArrayManufacture_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayData_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayDataCluster_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayDimension-method
(zzz-DTTCompatibility-class), 13

- name,BioAssayDimension_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayMap_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13
- name,BioMaterial_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioSample-method
(zzz-DTTCompatibility-class), 13
- name,BioSample_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioSequence-method
(zzz-DTTCompatibility-class), 13
- name,BioSequence_ref-method
(zzz-DTTCompatibility-class), 13
- name,BioSource-method
(zzz-DTTCompatibility-class), 13
- name,BioSource_ref-method
(zzz-DTTCompatibility-class), 13
- name,Channel-method
(zzz-DTTCompatibility-class), 13
- name,Channel_ref-method
(zzz-DTTCompatibility-class), 13
- name,CompositeCompositeMap-method
(zzz-DTTCompatibility-class), 13
- name,CompositeCompositeMap_ref-method
(zzz-DTTCompatibility-class), 13
- name,CompositeGroup-method
(zzz-DTTCompatibility-class), 13
- name,CompositeGroup_ref-method
(zzz-DTTCompatibility-class), 13
- name,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- name,CompositeSequence_ref-method
(zzz-DTTCompatibility-class), 13
- name,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class), 13
- name,CompositeSequenceDimension_ref-method
(zzz-DTTCompatibility-class), 13
- name,Compound-method
(zzz-DTTCompatibility-class), 13
- name,Compound_ref-method
(zzz-DTTCompatibility-class), 13
- name,ConfidenceIndicator_ref-method
(zzz-DTTCompatibility-class), 13
- name>Contact_ref-method
(zzz-DTTCompatibility-class), 13
- name,Database-method
(zzz-DTTCompatibility-class), 13
- name,Database_ref-method
(zzz-DTTCompatibility-class), 13
- name,DerivedBioAssay-method
(zzz-DTTCompatibility-class), 13
- name,DerivedBioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- name,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- name,DerivedBioAssayData_ref-method
(zzz-DTTCompatibility-class), 13
- name,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- name,DerivedSignal_ref-method
(zzz-DTTCompatibility-class), 13
- name,DesignElement_ref-method
(zzz-DTTCompatibility-class), 13
- name,DesignElementDimension_ref-method
(zzz-DTTCompatibility-class), 13
- name,DesignElementMap_ref-method
(zzz-DTTCompatibility-class), 13
- name,Error-method
(zzz-DTTCompatibility-class), 13
- name,Error_ref-method
(zzz-DTTCompatibility-class), 13
- name,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- name,ExpectedValue_ref-method
(zzz-DTTCompatibility-class), 13
- name,Experiment-method
(zzz-DTTCompatibility-class), 13
- name,ExperimentalFactor-method
(zzz-DTTCompatibility-class), 13
- name,ExperimentalFactor_ref-method
(zzz-DTTCompatibility-class), 13
- name,FactorValue-method
(zzz-DTTCompatibility-class), 13
- name,FactorValue_ref-method
(zzz-DTTCompatibility-class), 13
- name,Failed-method
(zzz-DTTCompatibility-class), 13
- name,Failed_ref-method
(zzz-DTTCompatibility-class), 13
- name,Feature-method
(zzz-DTTCompatibility-class), 13
- name,Feature_ref-method
(zzz-DTTCompatibility-class), 13
- name,FeatureDimension-method
(zzz-DTTCompatibility-class), 13
- name,FeatureDimension_ref-method
(zzz-DTTCompatibility-class), 13
- name,FeatureExtraction-method
(zzz-DTTCompatibility-class), 13

- name,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- name,FeatureGroup_ref-method
(zzz-DTTCompatibility-class), 13
- name,FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- name,FeatureReporterMap_ref-method
(zzz-DTTCompatibility-class), 13
- name,Hardware-method
(zzz-DTTCompatibility-class), 13
- name,Hardware_ref-method
(zzz-DTTCompatibility-class), 13
- name,Hybridization-method
(zzz-DTTCompatibility-class), 13
- name,Image-method
(zzz-DTTCompatibility-class), 13
- name,Image_ref-method
(zzz-DTTCompatibility-class), 13
- name,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- name,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- name,LabeledExtract_ref-method
(zzz-DTTCompatibility-class), 13
- name,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- name,MeasuredBioAssay-method
(zzz-DTTCompatibility-class), 13
- name,MeasuredBioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- name,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- name,MeasuredBioAssayData_ref-method
(zzz-DTTCompatibility-class), 13
- name,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- name,MeasuredSignal_ref-method
(zzz-DTTCompatibility-class), 13
- name,NameValueType-method
(zzz-DTTCompatibility-class), 13
- name,NodeValue-method
(zzz-DTTCompatibility-class), 13
- name,Organization-method
(zzz-DTTCompatibility-class), 13
- name,Organization_ref-method
(zzz-DTTCompatibility-class), 13
- name,Parameter-method
(zzz-DTTCompatibility-class), 13
- name,Parameter_ref-method
(zzz-DTTCompatibility-class), 13
- name,Person-method
(zzz-DTTCompatibility-class), 13
- name,Person_ref-method
(zzz-DTTCompatibility-class), 13
- name,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- name,PhysicalArrayDesign_ref-method
(zzz-DTTCompatibility-class), 13
- name,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- name,PhysicalBioAssay_ref-method
(zzz-DTTCompatibility-class), 13
- name,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- name,PresentAbsent_ref-method
(zzz-DTTCompatibility-class), 13
- name,Protocol-method
(zzz-DTTCompatibility-class), 13
- name,Protocol_ref-method
(zzz-DTTCompatibility-class), 13
- name,PValue-method
(zzz-DTTCompatibility-class), 13
- name,PValue_ref-method
(zzz-DTTCompatibility-class), 13
- name,QuantitationType_ref-method
(zzz-DTTCompatibility-class), 13
- name,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class), 13
- name,QuantitationTypeDimension_ref-method
(zzz-DTTCompatibility-class), 13
- name,QuantitationTypeMap-method
(zzz-DTTCompatibility-class), 13
- name,QuantitationTypeMap_ref-method
(zzz-DTTCompatibility-class), 13
- name,Ratio-method
(zzz-DTTCompatibility-class), 13
- name,Ratio_ref-method
(zzz-DTTCompatibility-class), 13
- name,Reporter-method
(zzz-DTTCompatibility-class), 13
- name,Reporter_ref-method
(zzz-DTTCompatibility-class), 13
- name,ReporterCompositeMap-method
(zzz-DTTCompatibility-class), 13
- name,ReporterCompositeMap_ref-method
(zzz-DTTCompatibility-class), 13
- name,ReporterDimension-method
(zzz-DTTCompatibility-class), 13
- name,ReporterDimension_ref-method
(zzz-DTTCompatibility-class), 13
- name,ReporterGroup-method
(zzz-DTTCompatibility-class), 13
- name,ReporterGroup_ref-method
(zzz-DTTCompatibility-class), 13

- name,Security-method
(zzz-DTTCompatibility-class), 13
- name,Security_ref-method
(zzz-DTTCompatibility-class), 13
- name,SecurityGroup-method
(zzz-DTTCompatibility-class), 13
- name,SecurityGroup_ref-method
(zzz-DTTCompatibility-class), 13
- name,Software-method
(zzz-DTTCompatibility-class), 13
- name,Software_ref-method
(zzz-DTTCompatibility-class), 13
- name,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- name,SpecializedQuantitationType_ref-method
(zzz-DTTCompatibility-class), 13
- name,StandardQuantitationType_ref-method
(zzz-DTTCompatibility-class), 13
- name,Transformation-method
(zzz-DTTCompatibility-class), 13
- name,Treatment-method
(zzz-DTTCompatibility-class), 13
- name,UserAttribute-method
(zzz-ARRCompatibility-class), 12
- name,Zone-method
(zzz-DTTCompatibility-class), 13
- name,Zone_ref-method
(zzz-DTTCompatibility-class), 13
- names,AffxArray-method
(zzz-NetAffxCompatibility-class), 14
- names,NetAffxResource-method
(NetAffxResource-class), 5
- namespace (zzz-ARRCompatibility-class), 12
- namespace,UserAttribute-method
(zzz-ARRCompatibility-class), 12
- nameValueType
(zzz-DTTCompatibility-class), 13
- nameValueType,PropertySets_assnlist-method
(zzz-DTTCompatibility-class), 13
- nameValueType,QualityControlStatistics_assnlist-method
(zzz-DTTCompatibility-class), 13
- nameValueType,SummaryStatistics_assnlist-method
(zzz-DTTCompatibility-class), 13
- NameValueType-class
(zzz-DTTCompatibility-class), 13
- NetAffxCompatibility, 5, 14
- NetAffxCompatibility-class
(AffyCompatible-class), 3
- NetAffxResource, 3, 4, 4, 5, 9, 14
- NetAffxResource-accessors (accessors), 2
- NetAffxResource-class, 5
- newSequence
(zzz-DTTCompatibility-class), 13
- newSequence,MismatchInformation-method
(zzz-DTTCompatibility-class), 13
- node (zzz-DTTCompatibility-class), 13
- node,Nodes_assnlist-method
(zzz-DTTCompatibility-class), 13
- Node-class (zzz-DTTCompatibility-class), 13
- nodeContents
(zzz-DTTCompatibility-class), 13
- nodeContents,NodeContents_assnlist-method
(zzz-DTTCompatibility-class), 13
- NodeContents-class
(zzz-DTTCompatibility-class), 13
- nodeContents_assnlist
(zzz-DTTCompatibility-class), 13
- nodeContents_assnlist,Node-method
(zzz-DTTCompatibility-class), 13
- NodeContents_assnlist-class
(zzz-DTTCompatibility-class), 13
- nodes_assnlist
(zzz-DTTCompatibility-class), 13
- nodes_assnlist,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- nodes_assnlist,Node-method
(zzz-DTTCompatibility-class), 13
- Nodes_assnlist-class
(zzz-DTTCompatibility-class), 13
- nodeValue (zzz-DTTCompatibility-class), 13
- nodeValue,NodeValue_assnlist-method
(zzz-DTTCompatibility-class), 13
- NodeValue-class
(zzz-DTTCompatibility-class), 13
- nodeValue_assnlist
(zzz-DTTCompatibility-class), 13
- nodeValue_assnlist,Node-method
(zzz-DTTCompatibility-class), 13
- NodeValue_assnlist-class
(zzz-DTTCompatibility-class), 13
- normalizationDescription_assn
(zzz-DTTCompatibility-class), 13
- normalizationDescription_assn,ExperimentDesign-method
(zzz-DTTCompatibility-class), 13
- NormalizationDescription_assn-class
(zzz-DTTCompatibility-class), 13
- numArrays (zzz-DTTCompatibility-class), 13
- numArrays,ArrayGroup-method
(zzz-DTTCompatibility-class), 13

- numberOfFeatures
 - (zzz-DTTCompatibility-class), 13
- numberOfFeatures,ArrayDesign-method
 - (zzz-DTTCompatibility-class), 13
- numberOfFeatures,PhysicalArrayDesign-method
 - (zzz-DTTCompatibility-class), 13
- numFeaturesPerCol
 - (zzz-DTTCompatibility-class), 13
- numFeaturesPerCol,ZoneLayout-method
 - (zzz-DTTCompatibility-class), 13
- numFeaturesPerRow
 - (zzz-DTTCompatibility-class), 13
- numFeaturesPerRow,ZoneLayout-method
 - (zzz-DTTCompatibility-class), 13

- ontologyEntries_assnlist
 - (zzz-DTTCompatibility-class), 13
- ontologyEntries_assnlist,BioSequence-method
 - (zzz-DTTCompatibility-class), 13
- OntologyEntries_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Action_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Annotations_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Associations_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Category_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Characteristics_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,CompoundIndices_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,ControlType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,DataType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,DefectType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,FailTypes_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,FeatureShape_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,FiducialType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Format_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,MaterialType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,OntologyEntries_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Parameters_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,PolymerType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Roles_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Scale_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Species_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,SubstrateType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,SurfaceType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,TechnologyType_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Type_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Types_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,Value_assn-method
 - (zzz-DTTCompatibility-class), 13
- ontologyEntry,WarningType_assn-method
 - (zzz-DTTCompatibility-class), 13
- OntologyEntry-class
 - (zzz-DTTCompatibility-class), 13
- ontologyReference_assn
 - (zzz-DTTCompatibility-class), 13
- ontologyReference_assn,OntologyEntry-method
 - (zzz-DTTCompatibility-class), 13
- OntologyReference_assn-class
 - (zzz-DTTCompatibility-class), 13
- organization (zzz-DTTCompatibility-class), 13
- organization,Contact_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- organization,Organization_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- Organization-class
 - (zzz-DTTCompatibility-class), 13
- Organization_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- organization_ref
 - (zzz-DTTCompatibility-class), 13
- organization_ref,Affiliation_assnref-method
 - (zzz-DTTCompatibility-class), 13
- organization_ref,ArrayManufacturers_assnreflist-method
 - (zzz-DTTCompatibility-class), 13
- organization_ref,Contacts_assnreflist-method
 - (zzz-DTTCompatibility-class), 13
- organization_ref,DesignProviders_assnreflist-method
 - (zzz-DTTCompatibility-class), 13

- organization_ref,HardwareManufacturers_assnreflist-method (zzz-DTTCompatibility-class), 13
- organization_ref,Members_assnreflist-method (zzz-DTTCompatibility-class), 13
- organization_ref,Owner_assnreflist-method (zzz-DTTCompatibility-class), 13
- organization_ref,Parent_assnref-method (zzz-DTTCompatibility-class), 13
- organization_ref,Performer_assnref-method (zzz-DTTCompatibility-class), 13
- organization_ref,Providers_assnreflist-method (zzz-DTTCompatibility-class), 13
- organization_ref,SoftwareManufacturers_assnreflist-method (zzz-DTTCompatibility-class), 13
- organization_ref,SourceContact_assnreflist-method (zzz-DTTCompatibility-class), 13
- Organization_ref-class (zzz-DTTCompatibility-class), 13
- orientationMark (zzz-DTTCompatibility-class), 13
- orientationMark,ArrayGroup-method (zzz-DTTCompatibility-class), 13
- orientationMarkPosition (zzz-DTTCompatibility-class), 13
- orientationMarkPosition,ArrayGroup-method (zzz-DTTCompatibility-class), 13
- originalProjectName (zzz-ARRCompatibility-class), 12
- originalProjectName,ArraySetFile-method (zzz-ARRCompatibility-class), 12
- originRelativeTo (zzz-DTTCompatibility-class), 13
- originRelativeTo,Array-method (zzz-DTTCompatibility-class), 13
- otherKind (zzz-DTTCompatibility-class), 13
- otherKind,Measurement-method (zzz-DTTCompatibility-class), 13
- owner_assnreflist (zzz-DTTCompatibility-class), 13
- owner_assnreflist,Security-method (zzz-DTTCompatibility-class), 13
- Owner_assnreflist-class (zzz-DTTCompatibility-class), 13
- pages (zzz-DTTCompatibility-class), 13
- pages,BibliographicReference-method (zzz-DTTCompatibility-class), 13
- parameter (zzz-DTTCompatibility-class), 13
- parameter,ParameterTypes_assnlist-method (zzz-DTTCompatibility-class), 13
- Parameter-class (zzz-DTTCompatibility-class), 13
- parameter_ref (zzz-DTTCompatibility-class), 13
- parameter_ref,ParameterType_assnref-method (zzz-DTTCompatibility-class), 13
- Parameter_ref-class (zzz-DTTCompatibility-class), 13
- parameters_assnlist (zzz-DTTCompatibility-class), 13
- parameters_assnlist,BibliographicReference-method (zzz-DTTCompatibility-class), 13
- Parameters_assnlist-class (zzz-DTTCompatibility-class), 13
- parameterType_assnref (zzz-DTTCompatibility-class), 13
- parameterType_assnref,ParameterValue-method (zzz-DTTCompatibility-class), 13
- ParameterType_assnref-class (zzz-DTTCompatibility-class), 13
- parameterTypes_assnlist (zzz-DTTCompatibility-class), 13
- parameterTypes_assnlist,Hardware-method (zzz-DTTCompatibility-class), 13
- parameterTypes_assnlist,Protocol-method (zzz-DTTCompatibility-class), 13
- parameterTypes_assnlist,Software-method (zzz-DTTCompatibility-class), 13
- ParameterTypes_assnlist-class (zzz-DTTCompatibility-class), 13
- parameterValue (zzz-DTTCompatibility-class), 13
- parameterValue,ParameterValues_assnlist-method (zzz-DTTCompatibility-class), 13
- ParameterValue-class (zzz-DTTCompatibility-class), 13
- parameterValues_assnlist (zzz-DTTCompatibility-class), 13
- parameterValues_assnlist,HardwareApplication-method (zzz-DTTCompatibility-class), 13
- parameterValues_assnlist,ProtocolApplication-method (zzz-DTTCompatibility-class), 13
- parameterValues_assnlist,SoftwareApplication-method (zzz-DTTCompatibility-class), 13
- ParameterValues_assnlist-class (zzz-DTTCompatibility-class), 13
- parent_assnref (zzz-DTTCompatibility-class), 13
- parent_assnref,Organization-method (zzz-DTTCompatibility-class), 13
- Parent_assnref-class (zzz-DTTCompatibility-class), 13

- password (accessors), 2
- password,NetAffxResource-method
(NetAffxResource-class), 5
- patassignmentMethod
(zzz-ARRCompatibility-class), 12
- patassignmentMethod,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- performer_ assnref
(zzz-DTTCompatibility-class), 13
- performer_ assnref,Audit-method
(zzz-DTTCompatibility-class), 13
- Performer_ assnref-class
(zzz-DTTCompatibility-class), 13
- performers_ assnreflist
(zzz-DTTCompatibility-class), 13
- performers_ assnreflist,ProtocolApplication-method
(zzz-DTTCompatibility-class), 13
- Performers_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- Person-class (zzz-DTTCompatibility-class), 13
- Person_ assnlist-class
(zzz-DTTCompatibility-class), 13
- person_ ref (zzz-DTTCompatibility-class), 13
- person_ ref,ArrayManufacturers_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,Contacts_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,DesignProviders_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,HardwareManufacturers_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,Members_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,Owner_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,Performer_ assnref-method
(zzz-DTTCompatibility-class), 13
- person_ ref,Performers_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,Providers_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,SoftwareManufacturers_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- person_ ref,SourceContact_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- Person_ ref-class
(zzz-DTTCompatibility-class), 13
- phone (zzz-DTTCompatibility-class), 13
- phone,Organization-method
(zzz-DTTCompatibility-class), 13
- phone,Person-method
(zzz-DTTCompatibility-class), 13
- physicalArray
(zzz-ARRCompatibility-class), 12
- physicalArray,PhysicalArrays-method
(zzz-ARRCompatibility-class), 12
- PhysicalArray-class
(zzz-ARRCompatibility-class), 12
- physicalArrayDesign
(zzz-DTTCompatibility-class), 13
- physicalArrayDesign,ArrayDesign_ assnlist-method
(zzz-DTTCompatibility-class), 13
- PhysicalArrayDesign-class
(zzz-DTTCompatibility-class), 13
- physicalArrayDesign_ ref
(zzz-DTTCompatibility-class), 13
- physicalArrayDesign_ ref,ArrayDesign_ assnref-method
(zzz-DTTCompatibility-class), 13
- PhysicalArrayDesign_ ref-class
(zzz-DTTCompatibility-class), 13
- physicalArrays
(zzz-ARRCompatibility-class), 12
- physicalArrays,ArraySetFile-method
(zzz-ARRCompatibility-class), 12
- PhysicalArrays-class
(zzz-ARRCompatibility-class), 12
- physicalBioAssay
(zzz-DTTCompatibility-class), 13
- physicalBioAssay,BioAssay_ assnlist-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay,PhysicalBioAssay_ assnlist-method
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssay-class
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssay_ assnlist-class
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssay_ assnref-class
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref,BioAssay_ assnref-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref,BioAssays_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref,PhysicalBioAssay_ assnref-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref,PhysicalBioAssaySource_ assnref-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref,PhysicalBioAssayTarget_ assnref-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ ref,SourceBioAssays_ assnreflist-method
(zzz-DTTCompatibility-class), 13

- physicalBioAssay_ref,Target_assnref-method
(zzz-DTTCompatibility-class), 13
- physicalBioAssay_ref,TopLevelBioAssays_assnreflist-method
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssay_ref-class
(zzz-DTTCompatibility-class), 13
- physicalBioAssayData_assnlist
(zzz-DTTCompatibility-class), 13
- physicalBioAssayData_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssayData_assnlist-class
(zzz-DTTCompatibility-class), 13
- physicalBioAssaySource_assnref
(zzz-DTTCompatibility-class), 13
- physicalBioAssaySource_assnref,FeatureExtraction-method
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssaySource_assnref-class
(zzz-DTTCompatibility-class), 13
- PhysicalBioAssayTarget_assnref-class
(zzz-DTTCompatibility-class), 13
- polymerType_assn
(zzz-DTTCompatibility-class), 13
- polymerType_assn,BioSequence-method
(zzz-DTTCompatibility-class), 13
- PolymerType_assn-class
(zzz-DTTCompatibility-class), 13
- position (zzz-DTTCompatibility-class), 13
- position,Position_assn-method
(zzz-DTTCompatibility-class), 13
- Position-class
(zzz-DTTCompatibility-class), 13
- position_assn
(zzz-DTTCompatibility-class), 13
- position_assn,Feature-method
(zzz-DTTCompatibility-class), 13
- position_assn,Fiducial-method
(zzz-DTTCompatibility-class), 13
- Position_assn-class
(zzz-DTTCompatibility-class), 13
- positionDelta
(zzz-DTTCompatibility-class), 13
- positionDelta,PositionDelta_assn-method
(zzz-DTTCompatibility-class), 13
- PositionDelta-class
(zzz-DTTCompatibility-class), 13
- positionDelta_assn
(zzz-DTTCompatibility-class), 13
- positionDelta_assn,FeatureDefect-method
(zzz-DTTCompatibility-class), 13
- positionDelta_assn,ZoneDefect-method
(zzz-DTTCompatibility-class), 13
- PositionDelta_assn-class
(zzz-DTTCompatibility-class), 13
- presentAbsent
(zzz-DTTCompatibility-class), 13
- presentAbsent,QuantitationType_assnlist-method
(zzz-DTTCompatibility-class), 13
- PresentAbsent-class
(zzz-DTTCompatibility-class), 13
- presentAbsent_ref
(zzz-DTTCompatibility-class), 13
- presentAbsent_ref,QuantitationType_assnref-method
(zzz-DTTCompatibility-class), 13
- presentAbsent_ref,QuantitationTypes_assnreflist-method
(zzz-DTTCompatibility-class), 13
- presentAbsent_ref,SourcesQuantitationType_assnreflist-method
(zzz-DTTCompatibility-class), 13
- presentAbsent_ref,TargetQuantitationType_assnref-method
(zzz-DTTCompatibility-class), 13
- PresentAbsent_ref-class
(zzz-DTTCompatibility-class), 13
- producerTransformation_assn
(zzz-DTTCompatibility-class), 13
- producerTransformation_assn,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- ProducerTransformation_assn-class
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Array-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ArrayManufactureDeviation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Audit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioAssayDimension-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioAssayMapping-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13

- (zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioDataCube-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioDataTuples-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioMaterialMeasurement-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioSample-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioSequence-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,BioSource-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Channel-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,CompositeCompositeMap-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,CompositeGroup-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,CompositePosition-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,CompositeSequenceDimension-property
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Compound-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,CompoundMeasurement-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ConcentrationUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Database-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,DatabaseEntry-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,DerivedBioAssay-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Description-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,DesignElementMapping-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,DistanceUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Error-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Experiment-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ExperimentalFactor-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ExperimentDesign-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ExternalReference-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FactorValue-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Failed-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Feature-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureDefect-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureDimension-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureExtraction-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureInformation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureLocation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Fiducial-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Hardware-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,HardwareApplication-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Hybridization-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Image-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,MassUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,MeasuredBioAssay-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13

- (zzz-DTTCompatibility-class), 13
- propertySets_assnlist,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Measurement-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,MismatchInformation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,NameValueType-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Node-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,NodeContents-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,NodeValue-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,OntologyEntry-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Organization-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Parameter-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ParameterValue-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Person-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Position-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,PositionDelta-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Protocol-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ProtocolApplication-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,PValue-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,QuantitationTypeDimension-property
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,QuantitationTypeMap-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,QuantitationTypeMapping-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,QuantityUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Ratio-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Reporter-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ReporterCompositeMap-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ReporterDimension-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ReporterPosition-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Security-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,SecurityGroup-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,SeqFeature-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,SeqFeatureLocation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,SequencePosition-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Software-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,SoftwareApplication-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,TemperatureUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,TimeUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Transformation-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Treatment-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,VolumeUnit-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,Zone-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ZoneDefect-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ZoneGroup-method
(zzz-DTTCompatibility-class), 13
- propertySets_assnlist,ZoneLayout-method
(zzz-DTTCompatibility-class), 13
- PropertySets_assnlist-class
(zzz-DTTCompatibility-class), 13
- protocol,Protocol_assnlist-method
(zzz-DTTCompatibility-class), 13
- Protocol-class
(zzz-DTTCompatibility-class), 13
- protocol_assnlist
(zzz-DTTCompatibility-class), 13

- protocol_ assnlist, Protocol_ package-method
(zzz-DTTCompatibility-class), 13
- Protocol_ assnlist-class
(zzz-DTTCompatibility-class), 13
- protocol_ assnref
(zzz-DTTCompatibility-class), 13
- protocol_ assnref, ProtocolApplication-method
(zzz-DTTCompatibility-class), 13
- Protocol_ assnref-class
(zzz-DTTCompatibility-class), 13
- protocol_ package
(zzz-DTTCompatibility-class), 13
- protocol_ package, MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- Protocol_ package-class
(zzz-DTTCompatibility-class), 13
- protocol_ ref
(zzz-DTTCompatibility-class), 13
- protocol_ ref, Protocol_ assnref-method
(zzz-DTTCompatibility-class), 13
- Protocol_ ref-class
(zzz-DTTCompatibility-class), 13
- protocolApplication
(zzz-DTTCompatibility-class), 13
- protocolApplication, ProtocolApplications_ assnlist-method
(zzz-DTTCompatibility-class), 13
- ProtocolApplication-class
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, CompositeCompositeMap-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, FeatureExtraction-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, Hybridization-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, QuantitationTypeMap-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, ReporterCompositeMap-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, Transformation-method
(zzz-DTTCompatibility-class), 13
- protocolApplications_ assnlist, Treatment-method
(zzz-DTTCompatibility-class), 13
- ProtocolApplications_ assnlist-class
(zzz-DTTCompatibility-class), 13
- providers_ assnreflist
(zzz-DTTCompatibility-class), 13
- providers_ assnreflist, Experiment-method
(zzz-DTTCompatibility-class), 13
- Providers_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- publication (zzz-DTTCompatibility-class),
13
- publication, BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- publisher (zzz-DTTCompatibility-class), 13
- publisher, BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- pvalue (zzz-DTTCompatibility-class), 13
- pvalue, ConfidenceIndicator_ assnlist-method
(zzz-DTTCompatibility-class), 13
- pvalue, QuantitationType_ assnlist-method
(zzz-DTTCompatibility-class), 13
- PValue-class (zzz-DTTCompatibility-class),
13
- pvalue_ ref (zzz-DTTCompatibility-class),
13
- pvalue_ ref, ConfidenceIndicators_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- pvalue_ ref, QuantitationType_ assnref-method
(zzz-DTTCompatibility-class), 13
- pvalue_ ref, QuantitationTypes_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- pvalue_ ref, SourcesQuantitationType_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- pvalue_ ref, TargetQuantitationType_ assnref-method
(zzz-DTTCompatibility-class), 13
- PValue_ ref-class
(zzz-DTTCompatibility-class), 13
- quality (zzz-DTTCompatibility-class), 13
- quality, ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- quality, ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- qualityControlDescription_ assn
(zzz-DTTCompatibility-class), 13

- qualityControlDescription_assn,ExperimentDesign-QuantitationTypeDimension_assnlist-class (zzz-DTTCompatibility-class), 13
- QualityControlDescription_assn-class (zzz-DTTCompatibility-class), 13
- qualityControlStatistics_assnlist (zzz-DTTCompatibility-class), 13
- qualityControlStatistics_assnlist,ArrayManufactureQuantitationTypeDimension_assnref,MeasuredBioAssayData-ref-class (zzz-DTTCompatibility-class), 13
- qualityControlStatistics_assnlist,BioSample-methodQuantitationTypeDimension_assnref-class (zzz-DTTCompatibility-class), 13
- qualityControlStatistics_assnlist,BioSource-methodQuantitationTypeDimension_ref (zzz-DTTCompatibility-class), 13
- qualityControlStatistics_assnlist,LabeledExtract-methodQuantitationTypeDimension_ref,QuantitationDimension_assnlist (zzz-DTTCompatibility-class), 13
- QualityControlStatistics_assnlist-class (zzz-DTTCompatibility-class), 13
- quantitationDimension_assnref (zzz-DTTCompatibility-class), 13
- quantitationDimension_assnref,NodeContents-method (zzz-DTTCompatibility-class), 13
- QuantitationDimension_assnref-class (zzz-DTTCompatibility-class), 13
- quantitationType_assnlist (zzz-DTTCompatibility-class), 13
- quantitationType_assnlist,QuantitationType_packageQuantitationTypeMap_assnlist (zzz-DTTCompatibility-class), 13
- QuantitationType_assnlist-class (zzz-DTTCompatibility-class), 13
- quantitationType_assnref (zzz-DTTCompatibility-class), 13
- quantitationType_assnref,QuantitationTypeTuple-method (zzz-DTTCompatibility-class), 13
- QuantitationType_assnref-class (zzz-DTTCompatibility-class), 13
- quantitationType_package (zzz-DTTCompatibility-class), 13
- quantitationType_package,MAGE-ML-method (zzz-DTTCompatibility-class), 13
- QuantitationType_package-class (zzz-DTTCompatibility-class), 13
- QuantitationType_ref-class (zzz-DTTCompatibility-class), 13
- quantitationTypeDimension (zzz-DTTCompatibility-class), 13
- quantitationTypeDimension,QuantitationTypeDimension_assnlistTypeMapping_assn,Transformation-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension-class (zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnlist (zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_assnlist,BioAssayData_packageTypeMaps_assnreflist,DerivedSignal-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension_ref,QuantitationDimension_assnlist (zzz-DTTCompatibility-class), 13
- quantitationTypeDimension_ref,QuantitationTypeDimension_ref-class (zzz-DTTCompatibility-class), 13
- QuantitationTypeDimension_ref-class (zzz-DTTCompatibility-class), 13
- quantitationTypeMap (zzz-DTTCompatibility-class), 13
- quantitationTypeMap,QuantitationTypeMap_assnlist-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeMap-class (zzz-DTTCompatibility-class), 13
- quantitationTypeMap_assnlist (zzz-DTTCompatibility-class), 13
- quantitationTypeMap_assnlist,BioAssayData_package-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeMap_assnlist-class (zzz-DTTCompatibility-class), 13
- quantitationTypeMap_ref (zzz-DTTCompatibility-class), 13
- quantitationTypeMap_ref,QuantitationTypeMaps_assnreflist (zzz-DTTCompatibility-class), 13
- QuantitationTypeMap_ref-class (zzz-DTTCompatibility-class), 13
- quantitationTypeMapping (zzz-DTTCompatibility-class), 13
- quantitationTypeMapping,QuantitationTypeMapping_assn-m (zzz-DTTCompatibility-class), 13
- QuantitationTypeMapping-class (zzz-DTTCompatibility-class), 13
- quantitationTypeMapping_assn (zzz-DTTCompatibility-class), 13
- quantitationTypeMapping_assn,Transformation-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeMapping_assn-class (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,DerivedSignal-method (zzz-DTTCompatibility-class), 13

- quantitationTypeMaps_assnreflist,Error-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,ExpectedValue-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,Failed-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,MeasuredSignal-Ratio-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,PresentAbsent-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,PValue-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,QuantitationTypeMapping-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,Ratio-method (zzz-DTTCompatibility-class), 13
- quantitationTypeMaps_assnreflist,SpecializedQuantitationType-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- quantitationTypes_assnreflist (zzz-DTTCompatibility-class), 13
- quantitationTypes_assnreflist,QuantitationTypeDimension-method (zzz-DTTCompatibility-class), 13
- QuantitationTypes_assnreflist-class (zzz-DTTCompatibility-class), 13
- quantitationTypeTuple (zzz-DTTCompatibility-class), 13
- quantitationTypeTuple,QuantitationTypeTuples_assnlist-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeTuple-class (zzz-DTTCompatibility-class), 13
- quantitationTypeTuples_assnlist (zzz-DTTCompatibility-class), 13
- quantitationTypeTuples_assnlist,DesignElementTuple-method (zzz-DTTCompatibility-class), 13
- QuantitationTypeTuples_assnlist-class (zzz-DTTCompatibility-class), 13
- quantityUnit (zzz-DTTCompatibility-class), 13
- quantityUnit,Unit_assn-method (zzz-DTTCompatibility-class), 13
- QuantityUnit-class (zzz-DTTCompatibility-class), 13
- ratio (zzz-DTTCompatibility-class), 13
- ratio,QuantitationType_assnlist-method (zzz-DTTCompatibility-class), 13
- Ratio-class (zzz-DTTCompatibility-class), 13
- ratio_ref (zzz-DTTCompatibility-class), 13
- ratio_ref,QuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- ratio_ref,QuantitationTypes_assnreflist-method (zzz-DTTCompatibility-class), 13
- ratio_ref,SourcesQuantitationType_assnreflist-method (zzz-DTTCompatibility-class), 13
- ratio_ref,TargetQuantitationType_assnref-method (zzz-DTTCompatibility-class), 13
- Ratio-class (zzz-DTTCompatibility-class), 13
- readAnnotation (NetAffxResource-class), 5
- readAnnotation,NetAffxResource,character,character-method (NetAffxResource-class), 5
- readAnnotation,NetAffxResource,character,numeric-method (NetAffxResource-class), 5
- readAnnotation,NetAffxResource,missing,AffxAnnotation-method (NetAffxResource-class), 5
- readAnnotation,NetAffxResource,numeric,character-method (NetAffxResource-class), 5
- readAnnotation,NetAffxResource,numeric,numeric-method (NetAffxResource-class), 5
- readArr, 7
- readMage, 8
- readNetAffx, 9
- readXml, 7, 8, 9, 10, 11
- regions_assnlist (zzz-DTTCompatibility-class), 13
- regions_assnlist,SeqFeature-method (zzz-DTTCompatibility-class), 13
- Regions_assnlist-class (zzz-DTTCompatibility-class), 13
- releaseDate (zzz-DTTCompatibility-class), 13
- releaseDate,SoftwareApplication-method (zzz-DTTCompatibility-class), 13
- replacedLength (zzz-DTTCompatibility-class), 13
- replacedLength,MismatchInformation-method (zzz-DTTCompatibility-class), 13
- replicateDescription_assn (zzz-DTTCompatibility-class), 13
- replicateDescription_assn,ExperimentDesign-method (zzz-DTTCompatibility-class), 13
- ReplicateDescription_assn-class (zzz-DTTCompatibility-class), 13
- reporter (zzz-DTTCompatibility-class), 13
- reporter,DesignElement_assnlist-method (zzz-DTTCompatibility-class), 13
- reporter,Reporter_assnlist-method (zzz-DTTCompatibility-class), 13
- Reporter-class (zzz-DTTCompatibility-class), 13
- reporter_assnlist (zzz-DTTCompatibility-class), 13

- reporter_assnlist,DesignElement_package-method (zzz-DTTCompatibility-class), 13
- Reporter_assnlist-class (zzz-DTTCompatibility-class), 13
- reporter_assnref (zzz-DTTCompatibility-class), 13
- reporter_assnref,FeatureReporterMap-method (zzz-DTTCompatibility-class), 13
- reporter_assnref,ReporterPosition-method (zzz-DTTCompatibility-class), 13
- Reporter_assnref-class (zzz-DTTCompatibility-class), 13
- reporter_ref (zzz-DTTCompatibility-class), 13
- reporter_ref,DesignElement_assnref-method (zzz-DTTCompatibility-class), 13
- reporter_ref,Reporter_assnref-method (zzz-DTTCompatibility-class), 13
- reporter_ref,Reporters_assnreflist-method (zzz-DTTCompatibility-class), 13
- Reporter_ref-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMap (zzz-DTTCompatibility-class), 13
- reporterCompositeMap,DesignElementMap_assnlist-method (zzz-DTTCompatibility-class), 13
- reporterCompositeMap,ReporterCompositeMap_assnlist-method (zzz-DTTCompatibility-class), 13
- ReporterCompositeMap-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_assnlist (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_assnlist,DesignElement_package-method (zzz-DTTCompatibility-class), 13
- ReporterCompositeMap_assnlist-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_ref (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_ref,DesignElementMaps_assnreflist-method (zzz-DTTCompatibility-class), 13
- reporterCompositeMap_ref,ReporterCompositeMaps_assnreflist-method (zzz-DTTCompatibility-class), 13
- ReporterCompositeMap_ref-class (zzz-DTTCompatibility-class), 13
- reporterCompositeMaps_assnreflist (zzz-DTTCompatibility-class), 13
- reporterCompositeMaps_assnreflist,CompositeSequence_package-method (zzz-DTTCompatibility-class), 13
- ReporterCompositeMaps_assnreflist-class (zzz-DTTCompatibility-class), 13
- reporterDimension (zzz-DTTCompatibility-class), 13
- reporterDimension,DesignElementDimension_assnlist-method (zzz-DTTCompatibility-class), 13
- ReporterDimension-class (zzz-DTTCompatibility-class), 13
- reporterDimension_ref (zzz-DTTCompatibility-class), 13
- reporterDimension_ref,DesignElementDimension_assnref-method (zzz-DTTCompatibility-class), 13
- ReporterDimension_ref-class (zzz-DTTCompatibility-class), 13
- reporterGroup (zzz-DTTCompatibility-class), 13
- reporterGroup,ReporterGroup_assnlist-method (zzz-DTTCompatibility-class), 13
- ReporterGroup-class (zzz-DTTCompatibility-class), 13
- reporterGroup_assnlist (zzz-DTTCompatibility-class), 13
- reporterGroup_assnlist,ArrayDesign_package-method (zzz-DTTCompatibility-class), 13
- ReporterGroup_assnlist-class (zzz-DTTCompatibility-class), 13
- reporterGroup_ref (zzz-DTTCompatibility-class), 13
- reporterGroup_ref,ReporterGroups_assnreflist-method (zzz-DTTCompatibility-class), 13
- ReporterGroup_ref-class (zzz-DTTCompatibility-class), 13
- reporterGroups_assnreflist (zzz-DTTCompatibility-class), 13
- reporterGroups_assnreflist,ArrayDesign-method (zzz-DTTCompatibility-class), 13
- reporterGroups_assnreflist,PhysicalArrayDesign-method (zzz-DTTCompatibility-class), 13
- ReporterGroups_assnreflist-class (zzz-DTTCompatibility-class), 13
- reporterPosition (zzz-DTTCompatibility-class), 13
- reporterPosition,Coordinate_assn-method (zzz-DTTCompatibility-class), 13
- reporterPosition,ReporterPositionSources_assnlist-method (zzz-DTTCompatibility-class), 13
- ReporterPosition-class (zzz-DTTCompatibility-class), 13
- reporterPositionSources_assnlist (zzz-DTTCompatibility-class), 13
- reporterPositionSources_assnlist,ReporterCompositeMap-method (zzz-DTTCompatibility-class), 13
- ReporterPositionSources_assnlist-class (zzz-DTTCompatibility-class), 13
- reporters_assnreflist (zzz-DTTCompatibility-class), 13

- reporters__assnreflist,ReporterDimension-method
(zzz-DTTCompatibility-class), 13
- reporters__assnreflist,ReporterGroup-method
(zzz-DTTCompatibility-class), 13
- Reporters__assnreflist-class
(zzz-DTTCompatibility-class), 13
- required (zzz-ARRCompatibility-class), 12
- required,UserAttribute-method
(zzz-ARRCompatibility-class), 12
- roles__assnlist
(zzz-DTTCompatibility-class), 13
- roles__assnlist,Organization-method
(zzz-DTTCompatibility-class), 13
- roles__assnlist,Person-method
(zzz-DTTCompatibility-class), 13
- Roles__assnlist-class
(zzz-DTTCompatibility-class), 13

- scale__assn (zzz-DTTCompatibility-class),
13
- scale__assn,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- scale__assn,Error-method
(zzz-DTTCompatibility-class), 13
- scale__assn,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- scale__assn,Failed-method
(zzz-DTTCompatibility-class), 13
- scale__assn,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- scale__assn,NodeValue-method
(zzz-DTTCompatibility-class), 13
- scale__assn,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- scale__assn,PValue-method
(zzz-DTTCompatibility-class), 13
- scale__assn,Ratio-method
(zzz-DTTCompatibility-class), 13
- scale__assn,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- Scale__assn-class
(zzz-DTTCompatibility-class), 13
- security (zzz-DTTCompatibility-class), 13
- security,Security__assnlist-method
(zzz-DTTCompatibility-class), 13
- Security-class
(zzz-DTTCompatibility-class), 13
- security__assnlist
(zzz-DTTCompatibility-class), 13
- security__assnlist,AuditAndSecurity__package-method
(zzz-DTTCompatibility-class), 13
- Security__assnlist-class
(zzz-DTTCompatibility-class), 13

- security__assnref
(zzz-DTTCompatibility-class), 13
- security__assnref,Array-method
(zzz-DTTCompatibility-class), 13
- security__assnref,ArrayDesign-method
(zzz-DTTCompatibility-class), 13
- security__assnref,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- security__assnref,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- security__assnref,Audit-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioAssayCreation-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioAssayDataCluster-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioAssayDimension-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioAssayMap-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioSample-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioSequence-method
(zzz-DTTCompatibility-class), 13
- security__assnref,BioSource-method
(zzz-DTTCompatibility-class), 13
- security__assnref,Channel-method
(zzz-DTTCompatibility-class), 13
- security__assnref,CompositeCompositeMap-method
(zzz-DTTCompatibility-class), 13
- security__assnref,CompositeGroup-method
(zzz-DTTCompatibility-class), 13
- security__assnref,CompositeSequence-method
(zzz-DTTCompatibility-class), 13
- security__assnref,CompositeSequenceDimension-method
(zzz-DTTCompatibility-class), 13
- security__assnref,Compound-method
(zzz-DTTCompatibility-class), 13
- security__assnref,Database-method
(zzz-DTTCompatibility-class), 13
- security__assnref,DerivedBioAssay-method
(zzz-DTTCompatibility-class), 13
- security__assnref,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- security__assnref,DerivedSignal-method
(zzz-DTTCompatibility-class), 13
- security__assnref,Description-method
(zzz-DTTCompatibility-class), 13

- security_ assnref,Error-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Experiment-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ExperimentalFactor-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ExperimentDesign-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,FactorValue-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Failed-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Feature-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,FeatureDimension-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,FeatureExtraction-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,FeatureReporterMap-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Fiducial-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Hardware-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,HardwareApplication-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Hybridization-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Image-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,LabeledExtract-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,MAGE-ML-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ManufactureLIMS-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ManufactureLIMSBiomaterial-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,MeasuredBioAssay-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,MeasuredSignal-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Node-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,NodeContents-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Organization-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Parameter-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Person-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,PhysicalBioAssay-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,PresentAbsent-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Protocol-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ProtocolApplication-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,PValue-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,QuantitationTypeDimension-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,QuantitationTypeMap-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Ratio-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Reporter-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ReporterCompositeMap-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ReporterDimension-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,ReporterGroup-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Security-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,SecurityGroup-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,SeqFeature-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Software-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,SoftwareApplication-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,SpecializedQuantitationType-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Transformation-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Treatment-method
(zzz-DTTCompatibility-class), 13
- security_ assnref,Zone-method
(zzz-DTTCompatibility-class), 13

- Security_ assnref-class
(zzz-DTTCompatibility-class), 13
- security_ref (zzz-DTTCompatibility-class), 13
- security_ref,Security_ assnref-method
(zzz-DTTCompatibility-class), 13
- Security_ref-class
(zzz-DTTCompatibility-class), 13
- securityGroup
(zzz-DTTCompatibility-class), 13
- securityGroup,SecurityGroup_ assnlist-method
(zzz-DTTCompatibility-class), 13
- SecurityGroup-class
(zzz-DTTCompatibility-class), 13
- securityGroup_ assnlist
(zzz-DTTCompatibility-class), 13
- securityGroup_ assnlist,AuditAndSecurity_ packageshowMethod
(zzz-DTTCompatibility-class), 13
- SecurityGroup_ assnlist-class
(zzz-DTTCompatibility-class), 13
- securityGroup_ref
(zzz-DTTCompatibility-class), 13
- securityGroup_ref,SecurityGroups_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- SecurityGroup_ref-class
(zzz-DTTCompatibility-class), 13
- securityGroups_ assnreflist
(zzz-DTTCompatibility-class), 13
- securityGroups_ assnreflist,Security-method
(zzz-DTTCompatibility-class), 13
- SecurityGroups_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- seqFeature (zzz-DTTCompatibility-class), 13
- seqFeature,SeqFeatures_ assnlist-method
(zzz-DTTCompatibility-class), 13
- SeqFeature-class
(zzz-DTTCompatibility-class), 13
- seqFeatureLocation
(zzz-DTTCompatibility-class), 13
- seqFeatureLocation,Regions_ assnlist-method
(zzz-DTTCompatibility-class), 13
- seqFeatureLocation,Subregions_ assnlist-method
(zzz-DTTCompatibility-class), 13
- SeqFeatureLocation-class
(zzz-DTTCompatibility-class), 13
- seqFeatures_ assnlist
(zzz-DTTCompatibility-class), 13
- seqFeatures_ assnlist,BioSequence-method
(zzz-DTTCompatibility-class), 13
- SeqFeatures_ assnlist-class
(zzz-DTTCompatibility-class), 13
- sequenceDatabases_ assnlist
(zzz-DTTCompatibility-class), 13
- sequenceDatabases_ assnlist,BioSequence-method
(zzz-DTTCompatibility-class), 13
- SequenceDatabases_ assnlist-class
(zzz-DTTCompatibility-class), 13
- sequencePosition
(zzz-DTTCompatibility-class), 13
- sequencePosition,Coordinate_ assn-method
(zzz-DTTCompatibility-class), 13
- SequencePosition-class
(zzz-DTTCompatibility-class), 13
- serialNumber
(zzz-DTTCompatibility-class), 13
- serialNumber,HardwareApplication-method
(zzz-DTTCompatibility-class), 13
- show,TypedAssociation-method
(Internal-class), 3
- show,AffyCompatible-method
(AffyCompatible-class), 3
- show,NetAffxResource-method
(NetAffxResource-class), 5
- software (zzz-DTTCompatibility-class), 13
- software,Software_ assnlist-method
(zzz-DTTCompatibility-class), 13
- Software-class
(zzz-DTTCompatibility-class), 13
- software_ assnlist
(zzz-DTTCompatibility-class), 13
- software_ assnlist,Protocol_ package-method
(zzz-DTTCompatibility-class), 13
- Software_ assnlist-class
(zzz-DTTCompatibility-class), 13
- software_ assnref
(zzz-DTTCompatibility-class), 13
- software_ assnref,SoftwareApplication-method
(zzz-DTTCompatibility-class), 13
- Software_ assnref-class
(zzz-DTTCompatibility-class), 13
- software_ref
(zzz-DTTCompatibility-class), 13
- software_ref,Software_ assnref-method
(zzz-DTTCompatibility-class), 13
- software_ref,Softwares_ assnreflist-method
(zzz-DTTCompatibility-class), 13
- Software_ref-class
(zzz-DTTCompatibility-class), 13
- softwareApplication
(zzz-DTTCompatibility-class), 13
- softwareApplication,SoftwareApplications_ assnlist-method
(zzz-DTTCompatibility-class), 13
- SoftwareApplication-class

- (zzz-DTTCompatibility-class), 13
- softwareApplications_assnlist
 - (zzz-DTTCompatibility-class), 13
- softwareApplications_assnlist,ProtocolApplication-spacing
 - (zzz-DTTCompatibility-class), 13
- SoftwareApplications_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- softwareManufacturers_assnreflist
 - (zzz-DTTCompatibility-class), 13
- softwareManufacturers_assnreflist,Software-method
 - (zzz-DTTCompatibility-class), 13
- SoftwareManufacturers_assnreflist-class
 - (zzz-DTTCompatibility-class), 13
- softwares_assnreflist
 - (zzz-DTTCompatibility-class), 13
- softwares_assnreflist,Hardware-method
 - (zzz-DTTCompatibility-class), 13
- softwares_assnreflist,Protocol-method
 - (zzz-DTTCompatibility-class), 13
- softwares_assnreflist,Software-method
 - (zzz-DTTCompatibility-class), 13
- Softwares_assnreflist-class
 - (zzz-DTTCompatibility-class), 13
- sourceBioAssays_assnreflist
 - (zzz-DTTCompatibility-class), 13
- sourceBioAssays_assnreflist,BioAssayMap-method
 - (zzz-DTTCompatibility-class), 13
- SourceBioAssays_assnreflist-class
 - (zzz-DTTCompatibility-class), 13
- sourceBioMaterialMeasurements_assnlist
 - (zzz-DTTCompatibility-class), 13
- sourceBioMaterialMeasurements_assnlist,BioAssayMap-method
 - (zzz-DTTCompatibility-class), 13
- sourceBioMaterialMeasurements_assnlist,HybridizationAssay-class
 - (zzz-DTTCompatibility-class), 13
- sourceBioMaterialMeasurements_assnlist,TreatmentAssn,BioSequence-method
 - (zzz-DTTCompatibility-class), 13
- SourceBioMaterialMeasurements_assnlist-class
 - (zzz-DTTCompatibility-class), 13
- sourceContact_assnreflist
 - (zzz-DTTCompatibility-class), 13
- sourceContact_assnreflist,BioSource-method
 - (zzz-DTTCompatibility-class), 13
- SourceContact_assnreflist-class
 - (zzz-DTTCompatibility-class), 13
- sourcesQuantitationType_assnreflist
 - (zzz-DTTCompatibility-class), 13
- sourcesQuantitationType_assnreflist,QuantitationType_ref-class
 - (zzz-DTTCompatibility-class), 13
- SourcesQuantitationType_assnreflist-class
 - (zzz-DTTCompatibility-class), 13
- spacingBetweenCols
 - (zzz-DTTCompatibility-class), 13
- spacingBetweenCols,ZoneLayout-method
 - (zzz-DTTCompatibility-class), 13
- spacingBetweenRows
 - (zzz-DTTCompatibility-class), 13
- spacingBetweenRows,ZoneLayout-method
 - (zzz-DTTCompatibility-class), 13
- spacingsBetweenZonesX
 - (zzz-DTTCompatibility-class), 13
- spacingsBetweenZonesX,ZoneGroup-method
 - (zzz-DTTCompatibility-class), 13
- spacingsBetweenZonesY
 - (zzz-DTTCompatibility-class), 13
- spacingsBetweenZonesY,ZoneGroup-method
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType,QuantitationType_assnlist-method
 - (zzz-DTTCompatibility-class), 13
- SpecializedQuantitationType-class
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType_ref
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType_ref,QuantitationType_assnreflist
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType_ref,QuantitationTypes_assnreflist
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType_ref,SourcesQuantitationType_assnreflist
 - (zzz-DTTCompatibility-class), 13
- specializedQuantitationType_ref,TargetQuantitationType_assnreflist
 - (zzz-DTTCompatibility-class), 13
- SpecializedQuantitationType_ref-class
 - (zzz-DTTCompatibility-class), 13
- species_assn,BioSequence-method
 - (zzz-DTTCompatibility-class), 13
- species_assn,CompositeGroup-method
 - (zzz-DTTCompatibility-class), 13
- species_assn,FeatureGroup-method
 - (zzz-DTTCompatibility-class), 13
- species_assn,ReporterGroup-method
 - (zzz-DTTCompatibility-class), 13
- Species_assn-class
 - (zzz-DTTCompatibility-class), 13
- StandardQuantitationType_ref-class
 - (zzz-DTTCompatibility-class), 13
- StandardQuantitationType_ref-class
 - (zzz-DTTCompatibility-class), 13
- startCoord,MismatchInformation-method
 - (zzz-DTTCompatibility-class), 13
- strandType (zzz-DTTCompatibility-class),

- 13
- strandType,SeqFeatureLocation-method
(zzz-DTTCompatibility-class), 13
- subregions_ assnlist
(zzz-DTTCompatibility-class), 13
- subregions_ assnlist,SeqFeatureLocation-method
(zzz-DTTCompatibility-class), 13
- Subregions_ assnlist-class
(zzz-DTTCompatibility-class), 13
- substrateType_ assn
(zzz-DTTCompatibility-class), 13
- substrateType_ assn,ArrayGroup-method
(zzz-DTTCompatibility-class), 13
- SubstrateType_ assn-class
(zzz-DTTCompatibility-class), 13
- summaryStatistics_ assnlist
(zzz-DTTCompatibility-class), 13
- summaryStatistics_ assnlist,DerivedBioAssayData-method
(zzz-DTTCompatibility-class), 13
- summaryStatistics_ assnlist,MeasuredBioAssayData-method
(zzz-DTTCompatibility-class), 13
- SummaryStatistics_ assnlist-class
(zzz-DTTCompatibility-class), 13
- surfaceType_ assn
(zzz-DTTCompatibility-class), 13
- surfaceType_ assn,PhysicalArrayDesign-method
(zzz-DTTCompatibility-class), 13
- SurfaceType_ assn-class
(zzz-DTTCompatibility-class), 13
- target_ assnref
(zzz-DTTCompatibility-class), 13
- target_ assnref,BioAssayTreatment-method
(zzz-DTTCompatibility-class), 13
- target_ assnref,ImageAcquisition-method
(zzz-DTTCompatibility-class), 13
- Target_ assnref-class
(zzz-DTTCompatibility-class), 13
- targetQuantitationType_ assnref
(zzz-DTTCompatibility-class), 13
- targetQuantitationType_ assnref,Error-method
(zzz-DTTCompatibility-class), 13
- targetQuantitationType_ assnref,ExpectedValue-method
(zzz-DTTCompatibility-class), 13
- targetQuantitationType_ assnref,PValue-method
(zzz-DTTCompatibility-class), 13
- targetQuantitationType_ assnref,QuantitationType-method
(zzz-DTTCompatibility-class), 13
- TargetQuantitationType_ assnref-class
(zzz-DTTCompatibility-class), 13
- technologyType_ assn
(zzz-DTTCompatibility-class), 13
- technologyType_ assn,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- TechnologyType_ assn-class
(zzz-DTTCompatibility-class), 13
- temperatureUnit
(zzz-DTTCompatibility-class), 13
- temperatureUnit,Unit_ assn-method
(zzz-DTTCompatibility-class), 13
- TemperatureUnit-class
(zzz-DTTCompatibility-class), 13
- TemplateFile-class
(zzz-ARRCompatibility-class), 12
- timeUnit (zzz-DTTCompatibility-class), 13
- timeUnit,Unit_ assn-method
(zzz-DTTCompatibility-class), 13
- TimeUnit-class
(zzz-DTTCompatibility-class), 13
- tolerance,ArrayManufacture-method
(zzz-DTTCompatibility-class), 13
- tollFreePhone
(zzz-DTTCompatibility-class), 13
- tollFreePhone,Organization-method
(zzz-DTTCompatibility-class), 13
- tollFreePhone,Person-method
(zzz-DTTCompatibility-class), 13
- topLevelBioAssays_ assnreflist
(zzz-DTTCompatibility-class), 13
- topLevelBioAssays_ assnreflist,ExperimentDesign-method
(zzz-DTTCompatibility-class), 13
- TopLevelBioAssays_ assnreflist-class
(zzz-DTTCompatibility-class), 13
- transformation
(zzz-DTTCompatibility-class), 13
- transformation,ProducerTransformation_ assn-method
(zzz-DTTCompatibility-class), 13
- Transformation-class
(zzz-DTTCompatibility-class), 13
- treatment (zzz-DTTCompatibility-class), 13
- treatment,Treatments_ assnlist-method
(zzz-DTTCompatibility-class), 13
- Treatment-class
(zzz-DTTCompatibility-class), 13
- treatments_ assnlist
(zzz-DTTCompatibility-class), 13
- treatments_ assnlist,BioSample-method
(zzz-DTTCompatibility-class), 13
- treatments_ assnlist,BioSource-method
(zzz-DTTCompatibility-class), 13
- treatments_ assnlist,LabeledExtract-method
(zzz-DTTCompatibility-class), 13

- Treatments_ assnlist-class
(zzz-DTTCompatibility-class), 13
- type (zzz-DTTCompatibility-class), 13
- type,ArraySetFile-method
(zzz-ARRCompatibility-class), 12
- type,Measurement-method
(zzz-DTTCompatibility-class), 13
- type,NameValueType-method
(zzz-DTTCompatibility-class), 13
- type,PhysicalArray-method
(zzz-ARRCompatibility-class), 12
- type,TemplateFile-method
(zzz-ARRCompatibility-class), 12
- type,UserAttribute-method
(zzz-ARRCompatibility-class), 12
- type_ assn (zzz-DTTCompatibility-class),
13
- type_ assn,BioSample-method
(zzz-DTTCompatibility-class), 13
- type_ assn,BioSequence-method
(zzz-DTTCompatibility-class), 13
- type_ assn,DatabaseEntry-method
(zzz-DTTCompatibility-class), 13
- type_ assn,DerivedBioAssay-method
(zzz-DTTCompatibility-class), 13
- type_ assn,Hardware-method
(zzz-DTTCompatibility-class), 13
- type_ assn,NodeValue-method
(zzz-DTTCompatibility-class), 13
- type_ assn,Protocol-method
(zzz-DTTCompatibility-class), 13
- type_ assn,Software-method
(zzz-DTTCompatibility-class), 13
- Type_ assn-class
(zzz-DTTCompatibility-class), 13
- types_ assnlist
(zzz-DTTCompatibility-class), 13
- types_ assnlist,CompositeGroup-method
(zzz-DTTCompatibility-class), 13
- types_ assnlist,ExperimentDesign-method
(zzz-DTTCompatibility-class), 13
- types_ assnlist,FeatureGroup-method
(zzz-DTTCompatibility-class), 13
- types_ assnlist,ReporterGroup-method
(zzz-DTTCompatibility-class), 13
- Types_ assnlist-class
(zzz-DTTCompatibility-class), 13
- unit_ assn (zzz-DTTCompatibility-class),
13
- unit_ assn,Measurement-method
(zzz-DTTCompatibility-class), 13
- Unit_ assn-class
(zzz-DTTCompatibility-class), 13
- unitName (zzz-DTTCompatibility-class),
13
- unitName,ConcentrationUnit-method
(zzz-DTTCompatibility-class), 13
- unitName,DistanceUnit-method
(zzz-DTTCompatibility-class), 13
- unitName,MassUnit-method
(zzz-DTTCompatibility-class), 13
- unitName,QuantityUnit-method
(zzz-DTTCompatibility-class), 13
- unitName,TemperatureUnit-method
(zzz-DTTCompatibility-class), 13
- unitName,TimeUnit-method
(zzz-DTTCompatibility-class), 13
- unitName,VolumeUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV
(zzz-DTTCompatibility-class), 13
- unitNameCV,ConcentrationUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV,DistanceUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV,MassUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV,QuantityUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV,TemperatureUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV,TimeUnit-method
(zzz-DTTCompatibility-class), 13
- unitNameCV,VolumeUnit-method
(zzz-DTTCompatibility-class), 13
- upperLeftX (zzz-DTTCompatibility-class),
13
- upperLeftX,Zone-method
(zzz-DTTCompatibility-class), 13
- upperLeftY (zzz-DTTCompatibility-class),
13
- upperLeftY,Zone-method
(zzz-DTTCompatibility-class), 13
- uri (zzz-DTTCompatibility-class), 13
- uri,BibliographicReference-method
(zzz-DTTCompatibility-class), 13
- uri,Database-method
(zzz-DTTCompatibility-class), 13
- uri,DatabaseEntry-method
(zzz-DTTCompatibility-class), 13
- uri,Description-method
(zzz-DTTCompatibility-class), 13
- uri,Hardware-method

- (zzz-DTTCompatibility-class), 13
- uri,Image-method
 - (zzz-DTTCompatibility-class), 13
- uri,Organization-method
 - (zzz-DTTCompatibility-class), 13
- uri,Person-method
 - (zzz-DTTCompatibility-class), 13
- uri,Protocol-method
 - (zzz-DTTCompatibility-class), 13
- uri,Software-method
 - (zzz-DTTCompatibility-class), 13
- user (accessors), 2
- user,NetAffxResource-method
 - (NetAffxResource-class), 5
- userAttribute
 - (zzz-ARRCompatibility-class), 12
- userAttribute,UserAttributes-method
 - (zzz-ARRCompatibility-class), 12
- UserAttribute-class
 - (zzz-ARRCompatibility-class), 12
- userAttributes
 - (zzz-ARRCompatibility-class), 12
- userAttributes,ArraySetFile-method
 - (zzz-ARRCompatibility-class), 12
- userAttributes,TemplateFile-method
 - (zzz-ARRCompatibility-class), 12
- UserAttributes-class
 - (zzz-ARRCompatibility-class), 12
- userAttributeValue
 - (zzz-ARRCompatibility-class), 12
- userAttributeValue,UserAttribute-method
 - (zzz-ARRCompatibility-class), 12
- UserAttributeValue-class
 - (zzz-ARRCompatibility-class), 12
- value (zzz-ARRCompatibility-class), 12
- value,Control-method
 - (zzz-ARRCompatibility-class), 12
- value,Datum-method
 - (zzz-DTTCompatibility-class), 13
- value,Measurement-method
 - (zzz-DTTCompatibility-class), 13
- value,NameValueType-method
 - (zzz-DTTCompatibility-class), 13
- value,NodeValue-method
 - (zzz-DTTCompatibility-class), 13
- value,OntologyEntry-method
 - (zzz-DTTCompatibility-class), 13
- value,ParameterValue-method
 - (zzz-DTTCompatibility-class), 13
- value_ assn (zzz-DTTCompatibility-class), 13
- value_ assn,FactorValue-method
 - (zzz-DTTCompatibility-class), 13
- Value_ assn-class
 - (zzz-DTTCompatibility-class), 13
- vector, 4
- version (zzz-ARRCompatibility-class), 12
- version,ArrayDesign-method
 - (zzz-DTTCompatibility-class), 13
- version,ArraySetFile-method
 - (zzz-ARRCompatibility-class), 12
- version,Database-method
 - (zzz-DTTCompatibility-class), 13
- version,PhysicalArrayDesign-method
 - (zzz-DTTCompatibility-class), 13
- version,SoftwareApplication-method
 - (zzz-DTTCompatibility-class), 13
- version,TemplateFile-method
 - (zzz-ARRCompatibility-class), 12
- volume (zzz-DTTCompatibility-class), 13
- volume,BibliographicReference-method
 - (zzz-DTTCompatibility-class), 13
- volumeUnit (zzz-DTTCompatibility-class), 13
- volumeUnit,Unit_ assn-method
 - (zzz-DTTCompatibility-class), 13
- VolumeUnit-class
 - (zzz-DTTCompatibility-class), 13
- warningType_ assn
 - (zzz-DTTCompatibility-class), 13
- warningType_ assn,Reporter-method
 - (zzz-DTTCompatibility-class), 13
- WarningType_ assn-class
 - (zzz-DTTCompatibility-class), 13
- width (zzz-DTTCompatibility-class), 13
- width,ArrayGroup-method
 - (zzz-DTTCompatibility-class), 13
- x (zzz-DTTCompatibility-class), 13
- x,Position-method
 - (zzz-DTTCompatibility-class), 13
- xclass, 7–10, 10
- xmlValue.XMLAttributeValue, 11
- y (zzz-DTTCompatibility-class), 13
- y,Position-method
 - (zzz-DTTCompatibility-class), 13
- year (zzz-DTTCompatibility-class), 13
- year,BibliographicReference-method
 - (zzz-DTTCompatibility-class), 13
- zone (zzz-DTTCompatibility-class), 13
- zone,ZoneLocations_ assnlist-method
 - (zzz-DTTCompatibility-class), 13

- Zone-class (*zzz-DTTCompatibility-class*),
13
- zone_ assnref
(*zzz-DTTCompatibility-class*), 13
- zone_ assnref,Feature-method
(*zzz-DTTCompatibility-class*), 13
- zone_ assnref,ZoneDefect-method
(*zzz-DTTCompatibility-class*), 13
- Zone_ assnref-class
(*zzz-DTTCompatibility-class*), 13
- zone_ ref (*zzz-DTTCompatibility-class*), 13
- zone_ ref,Zone_ assnref-method
(*zzz-DTTCompatibility-class*), 13
- Zone_ ref-class
(*zzz-DTTCompatibility-class*), 13
- zoneDefect (*zzz-DTTCompatibility-class*),
13
- zoneDefect,Adjustments_ assnlist-method
(*zzz-DTTCompatibility-class*), 13
- ZoneDefect-class
(*zzz-DTTCompatibility-class*), 13
- zoneGroup (*zzz-DTTCompatibility-class*),
13
- zoneGroup,ZoneGroups_ assnlist-method
(*zzz-DTTCompatibility-class*), 13
- ZoneGroup-class
(*zzz-DTTCompatibility-class*), 13
- zoneGroups_ assnlist
(*zzz-DTTCompatibility-class*), 13
- zoneGroups_ assnlist,PhysicalArrayDesign-method
(*zzz-DTTCompatibility-class*), 13
- ZoneGroups_ assnlist-class
(*zzz-DTTCompatibility-class*), 13
- zoneLayout (*zzz-DTTCompatibility-class*),
13
- zoneLayout,ZoneLayout_ assn-method
(*zzz-DTTCompatibility-class*), 13
- ZoneLayout-class
(*zzz-DTTCompatibility-class*), 13
- zoneLayout_ assn
(*zzz-DTTCompatibility-class*), 13
- zoneLayout_ assn,ZoneGroup-method
(*zzz-DTTCompatibility-class*), 13
- ZoneLayout_ assn-class
(*zzz-DTTCompatibility-class*), 13
- zoneLocations_ assnlist
(*zzz-DTTCompatibility-class*), 13
- zoneLocations_ assnlist,ZoneGroup-method
(*zzz-DTTCompatibility-class*), 13
- ZoneLocations_ assnlist-class
(*zzz-DTTCompatibility-class*), 13
- zonesPerX (*zzz-DTTCompatibility-class*),
13
- zonesPerX,ZoneGroup-method
(*zzz-DTTCompatibility-class*), 13
- zonesPerY (*zzz-DTTCompatibility-class*),
13
- zonesPerY,ZoneGroup-method
(*zzz-DTTCompatibility-class*), 13
- zzz-.TypedAssociation-class
(Internal-class), 3
- zzz-AffyCompatible-class
(AffyCompatible-class), 3
- zzz-ARRCompatibility-class, 12
- zzz-DTTCompatibility-class, 13
- zzz-NetAffxCompatibility-class, 14