

Package ‘dexus’

October 7, 2014

Type Package

Title DEXUS - Identifying Differential Expression in RNA-Seq Studies with Unknown Conditions or without Replicates

Description DEXUS identifies differentially expressed genes in RNA-Seq data under all possible study designs such as studies without replicates, without sample groups, and with unknown conditions. DEXUS works also for known conditions, for example for RNA-Seq data with two or multiple conditions. RNA-Seq read count data can be provided both by the S4 class Count Data Set and by read count matrices. Differentially expressed transcripts can be visualized by heatmaps, in which unknown conditions, replicates, and samples groups are also indicated. This software is fast since the core algorithm is written in C. For very large data sets, a parallel version of DEXUS is provided in this package. DEXUS is a statistical model that is selected in a Bayesian framework by an EM algorithm. DEXUS does not need replicates to detect differentially expressed transcripts, since the replicates (or conditions) are estimated by the EM method for each transcript. The method provides an informative/non-informative value to extract differentially expressed transcripts at a desired significance level or power.

Version 1.4.0

Date 2013-11-05

Maintainer Guenter Klambauer <klambauer@bioinf.jku.at>

Author Guenter Klambauer

License LGPL (>= 2.0)

Depends R (>= 2.15), methods, BiocGenerics

Suggests parallel, statmod, stats, DESeq, RColorBrewer

Collate 'AllClasses.R' 'AllGenerics.R' 'binomTest.R' 'normalization.R'
'dexus.R' 'getSizeNB.R' 'functions.R' 'plot-methods.R'
'show-methods.R' 'methodsAccess.R' 'dexss.R'

biocViews

Sequencing, Mus_musculus, Homo_sapiens, Zea_Mays, Macaca_mulatta, Pan_troglodytes, RNASeq, GeneExpression, DifferentialExpression, CellBiology, HapMap, RNAseqData, Classification, QualityControl

R topics documented:

accessors	2
countsBottomly	3
countsGilad	4
countsLi	4
countsMontgomery	5
countsPickrell	5
dexss	6
dexus	8
dexus.parallel	11
DEXUSResult-class	12
getSizeNB	14
INI	15
INIThreshold<-	16
normalizeData	16
plot	17
sort	18
['	18

Index	20
--------------	-----------

accessors	<i>Accessors for a "DEXUSResult".</i>
-----------	---------------------------------------

Description

These generic functions return the slots of an RNA-Seq analysis performed by DEXUS. The results of DEXUS are stored as an instance of `DEXUSResult-class`.

Arguments

object An instance of "DEXUSResult".

Value

The accessor functions return a the matrices or vectors contained in the corresponding slot of the "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
result <- dexus(countsBottomly[1:20,1:10])
transcriptNames(result)
sampleNames(result)
inputData(result)
normalizedData(result)
sizeFactors(result)
INIValues(result)
INIThreshold(result)
INICalls(result)
pvals(result)
responsibilities(result)
posteriorProbs(result)
logFC(result)
conditionSizes(result)
sizeParameters(result)
means(result)
dispersions(result)
params(result)
```

countsBottomly

RNA-Seq data of two mice strains.

Description

The two common mice strains C57BL/6J (B6) and DBA/2J (D2) were used for comparing gene expression measures of RNA-Seq and microarrays.

Usage

```
countsBottomly
```

Format

A data matrix of 36229 rows (genes) and 21 columns (samples).

Source

<http://bowtie-bio.sourceforge.net/recount/>

References

Bottomly, D., Walter, N. A. R., Hunter, J. E., Darakjian, P., Kawane, S., Buck, K. J., Searles, R. P., Mooney, M., McWeeney, S. K., and Hitzemann, R. (2011). *Evaluating gene expression in C57BL/6J and DBA/2J mouse striatum using RNA-Seq and microarrays*. Plos One, 6(3), e17820.

countsGilad

RNA-Seq data of humans, chimpanzees and rhesus macaques.

Description

Liver RNA samples of three males and three females from each of the species human, chimpanzee and rhesus macaques were sequenced.

Usage

countsGilad

Format

A data matrix of 20689 rows (genes) and 18 columns (samples).

Source

ftp://ftp.ncbi.nlm.nih.gov/pub/geo/DATA/supplementary/series/GSE17274/GSE17274_ReadCountPerLane.txt.gz

References

Blekhman, R., Marioni, J. C., Zumbo, P., Stephens, M., and Gilad, Y. (2010). *Sex-specific and lineage-specific alternative splicing in primates*. *Genome Res*, 20(2), 180-189.

countsLi

RNA-Seq data of the developmental zones of maize leaves.

Description

RNA-Sequencing was performed on different locations of the maize plant leaf.

Usage

countsLi

Format

A data matrix of 110185 rows (genes) and 12 columns (samples).

Source

<http://www.ncbi.nlm.nih.gov/sra/> accession number:SRP002265

References

Li, P., Ponnala, L., Gandotra, N., Wang, L., Si, Y., Tausta, S. L., Kebrom, T. H., Provard, N., Patel, R., Myers, C. R., Reidel, E. J., Turgeon, R., Liu, P., Sun, Q., Nelson, T., and Brutnell, T. P. (2010). *The developmental dynamics of the maize leaf transcriptome*. Nat Genet, 42(12), 1060-1067.

countsMontgomery *RNA-Seq data of 60 European HapMap individuals.*

Description

The RNA of lymphoblastoid cell lines of 60 HapMap individuals was sequenced in order to study eQTLs.

Usage

countsMontgomery

Format

A data matrix of 12984 rows (genes) and 60 columns (samples).

Source

<http://bowtie-bio.sourceforge.net/recount/>

References

Montgomery, S. B., Sammeth, M., Gutierrez-Arcelus, M., Lach, R. P., Ingle, C., Nisbett, J., Guigo, R., and Dermitzakis, E. T. (2010). *Transcriptome genetics using second generation sequencing in a caucasian population*. Nature, 464(7289), 773-777.

countsPickrell *RNA-Seq data of 69 Nigerian HapMap individuals.*

Description

The RNA of lymphoblastoid cell lines of 69 HapMap individuals was sequenced in order to study eQTLs.

Usage

countsPickrell

Format

A data matrix of 12984 rows (genes) and 69 columns (samples).

Source

<http://bowtie-bio.sourceforge.net/recount/>

References

Pickrell, J. K., Marioni, J. C., Pai, A. A., Degner, J. F., Engelhardt, B. E., Nkadori, E., Veyrieras, J.-B., Stephens, M., Gilad, Y., and Pritchard, J. K. (2010). *Understanding mechanisms underlying human gene expression variation with RNA sequencing*. *Nature*, 464(7289), 768-772.

dexss

Detection of Differential Expression in a semi-supervised Setting

Description

Performs the DEXSS algorithm for detection of differentially expressed genes in RNA-seq data for a semi-supervised setting, i.e. that the condition of some samples is known, and for some samples the condition is unknown.

Usage

```
dexss(X, nclasses = 2, G = 1, alphaInit, cyc = 20,
      labels, normalization = "RLE", kmeansIter = 10,
      ignoreIfAllCountsSmaller = 1, theta = 2.5, minMu = 0.5,
      rmax = 13, initialization = "kmeans",
      multiclassPhiPoolingFunction = NULL, quiet = FALSE,
      resultObject = "S4")
```

Arguments

X	either a vector of counts or a raw data matrix, where columns are interpreted as samples and rows as genomic regions. An instance of "countDataSet" is also accepted.
nclasses	The number of conditions, i.e. mixture components. (Default = 2)
G	The weight of the prior distribution of the mixture weights. Not used in the supervised case. (Default = 1).
cyc	Positive integer that sets the number of cycles of the EM algorithm. (Default = 20).
alphaInit	The initial estimates of the condition sizes, i.e., mixture weights. Not used in the supervised case. (Default = c(0.5,0.5)) .
labels	The labels for the classes, will be coerced into an integer. For this semi-supervised version the known labels/conditions must be coded as integers starting with 1. The samples with the label 1 will be considered as being in the "major condition". For the samples with unknown labels/conditions an "NA" must be set.
normalization	method used for normalizing the reads. "RLE" is the method used by (Anders and Huber, 2010), "upperquartile" is the Upper-Quartile method by (Bullard et al., 2010), and none deactivates normalization. (Default = "RLE").

kmeansIter	number of times the K-Means algorithm is run. (Default = 10).
ignoreIfAllCountsSmaller	Ignores transcript for which all read counts are smaller than this value. These transcripts are considered as "not expressed" (Default = 1).
theta	The weight of the prior on the size parameter or inverse dispersion parameter. Theta is adjusted to each transcript by dividing by the mean read count of the transcript. The higher theta, the lower r and the higher the overdispersion will be. (Default = 2.5).
minMu	Minimal mean for all negative binomial distributions. (Default = 0.5).
rmax	Maximal value for the size parameter. The inverse of this parameter is the lower bound on the dispersion. In analogy to (Anders and Huber, 2010) we use 13 as default. (Default = 13).
initialization	Method used to find the initial clusters. Dexss can either use the quantiles of the readcounts of each gene or run k-means on the counts. (Default = "kmeans").
multiclassPhiPoolingFunction	In "multiClass" mode the dispersion is either estimated across all classes at once (NULL), or separately for each condition, i.e., class. The size parameters or dispersion per class are then joined to one estimate by the mean ("mean"), minimum ("min") or maximum ("max"). In our investigations estimation across all classes at once performed best. (Default = NULL).
quiet	Logical that indicates whether dexss should report the steps of the algorithm. Suppresses messages from the program if set to TRUE. (Default = FALSE).
resultObject	Type of the result object; can either be a list ("list") or an instance of "DEXUS-Result" ("S4"). (Default="S4").

Details

The read count x is explained by a finite mixture of negative binomials:

$$p(x) = \sum_{i=1}^n \alpha_i \text{NB}(x; \mu_i, r_i),$$

where α_i is the weight of the mixture component, NB is the negative binomial with mean parameter μ_i and size parameter r_i . The parameters are selected by an EM algorithm in a Bayesian framework.

Each component in the mixture model corresponds to one condition.

- If the groups, conditions, replicate status, or labels are unknown, DEXUS tries to estimate these conditions. For each transcript DEXUS tries to explain the read counts by one negative binomial distribution. If this is possible, the transcript is explained by one condition and therefore it is not differentially expressed. If more than one negative binomial distribution is needed to explain the read counts of a transcript, this transcript indicates that it is differentially expressed. Evidence for differential expression is strong if a large amount of samples participate in each condition and the mean expression values are well separated. Both of these criteria are measured by the informative/non-informative (I/NI) call.
- If there are more than two groups given by the vector labels, DEXUS uses a generalized linear model to explain the data in analogy to (McCarthy, 2012).

- If there are two groups given by the vector labels, DEXUS uses the exact test for count data to test between the sample groups, as implemented by (Anders and Huber, 2010) in the package "DESeq".

Value

"list" or "DEXUSResult". A list containing the results and the parameters of the algorithm or an instance of "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Untertiner <unterthiner@bioinf.jku.at>

References

- Anders, S. and Huber, W. (2010). *Differential expression analysis for sequence count data*. Genome Biol, 11(10), R106.
- Bullard, J. H., Purdom, E., Hansen, K. D., and Dudoit, S. (2010). *Evaluation of statistical methods for normalization and differential expression in mRNA-seq experiments*. BMC Bioinformatics, 11, 94.
- McCarthy, D. J., Chen, Y., and Smyth, G. K. (2012). *Differential expression analysis of multifactor RNA-Seq experiments with respect to biological variation*. Nucleic Acids Res, 40(10), 4288-4297.

Examples

```
data(dexus)
labels1 <- substr(colnames(countsBottomly),1,2)
labels2 <- c()
labels2[which(labels1=="D2")] <- 1
labels2[which(labels1=="B6")] <- 2
labels2[c(3,7,8,10,12,15)] <- NA
res <- dexss(countsBottomly[1:100, ],labels=labels2,nclasses=2,G=0)
```

dexus

Detection of Differential Expression in an Unsupervised Setting

Description

Performs the DEXUS algorithm for detection of differentially expressed genes in RNA-seq data for a) unknown conditions, b) multiple known conditions, and c) two known conditions.

Usage

```
dexus(X, nclasses = 2, alphaInit, G = 1, cyc = 20,
      labels = NULL, normalization = "RLE", kmeansIter = 10,
      ignoreIfAllCountsSmaller = 1, theta = 2.5, minMu = 0.5,
      rmax = 13, initialization = "kmeans",
      multiclassPhiPoolingFunction = NULL, quiet = FALSE,
      resultObject = "S4")
```


Arguments

X	either a vector of counts or a raw data matrix, where columns are interpreted as samples and rows as genomic regions. An instance of "countDataSet" is also accepted.
nclasses	The number of conditions, i.e. mixture components. (Default = 2)
alphaInit	The initial estimates of the condition sizes, i.e., mixture weights. Not used in the supervised case. (Default = c(0.5,0.5)) .
G	The weight of the prior distribution of the mixture weights. Not used in the supervised case. (Default = 1).
cyc	Positive integer that sets the number of cycles of the EM algorithm. (Default = 20).
labels	labels for the classes, will be coerced into a factor by as.factor. Can either be a factor, character or integer. If this vector is given, supervised detection is used. If this vector is set to NULL the unsupervised detection is performed. (Default=NULL).
normalization	method used for normalizing the reads. "RLE" is the method used by (Anders and Huber, 2010), "upperquartile" is the Upper-Quartile method by (Bullard et al., 2010), and none deactivates normalization. (Default = "RLE").
kmeansIter	number of times the K-Means algorithm is run. (Default = 10).
ignoreIfAllCountsSmaller	Ignores transcript for which all read counts are smaller than this value. These transcripts are considered as "not expressed" (Default = 1).
theta	The weight of the prior on the size parameter or inverse dispersion parameter. Theta is adjusted to each transcript by dividing by the mean read count of the transcript. The higher theta, the lower r and the higher the overdispersion will be. (Default = 2.5).
minMu	Minimal mean for all negative binomial distributions. (Default = 0.5).
rmax	Maximal value for the size parameter. The inverse of this parameter is the lower bound on the dispersion. In analogy to (Anders and Huber, 2010) we use 13 as default. (Default = 13).
initialization	Method used to find the initial clusters. Dexus can either use the quantiles of the readcounts of each gene or run k-means on the counts. (Default = "kmeans").
multiclassPhiPoolingFunction	In "multiClass" mode the dispersion is either estimated across all classes at once (NULL), or separately for each condition, i.e., class. The size parameters or dispersion per class are then joined to one estimate by the mean ("mean"), minimum ("min") or maximum ("max"). In our investigations estimation across all classes at once performed best. (Default = NULL).
quiet	Logical that indicates whether dexus should report the steps of the algorithm. Suppresses messages from the program if set to TRUE. (Default = FALSE).
resultObject	Type of the result object; can either be a list ("list") or an instance of "DEXUS-Result" ("S4"). (Default="S4").

Details

The read count x is explained by a finite mixture of negative binomials:

$$p(x) = \sum_{i=1}^n \alpha_i \text{NB}(x; \mu_i, r_i),$$

where α_i is the weight of the mixture component, NB is the negative binomial with mean parameter μ_i and size parameter r_i . The parameters are selected by an EM algorithm in a Bayesian framework. Each component in the mixture model corresponds to one condition.

- If the groups, conditions, replicate status, or labels are unknown, DEXUS tries to estimate these conditions. For each transcript DEXUS tries to explain the read counts by one negative binomial distribution. If this is possible, the transcript is explained by one condition and therefore it is not differentially expressed. If more than one negative binomial distribution is needed to explain the read counts of a transcript, this transcript indicates that it is differentially expressed. Evidence for differential expression is strong if a large amount of samples participate in each condition and the mean expression values are well separated. Both of these criteria are measured by the informative/non-informative (I/NI) call.
- If there are more than two groups given by the vector labels, DEXUS uses a generalized linear model to explain the data in analogy to (McCarthy, 2012).
- If there are two groups given by the vector labels, DEXUS uses the exact test for count data to test between the sample groups, as implemented by (Anders and Huber, 2010) in the package "DESeq".

Value

"list" or "DEXUSResult". A list containing the results and the parameters of the algorithm or an instance of "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Untertiner <unterthiner@bioinf.jku.at>

References

- Anders, S. and Huber, W. (2010). *Differential expression analysis for sequence count data*. Genome Biol, 11(10), R106.
- Bullard, J. H., Purdom, E., Hansen, K. D., and Dudoit, S. (2010). *Evaluation of statistical methods for normalization and differential expression in mRNA-seq experiments*. BMC Bioinformatics, 11, 94.
- McCarthy, D. J., Chen, Y., and Smyth, G. K. (2012). *Differential expression analysis of multifactor RNA-Seq experiments with respect to biological variation*. Nucleic Acids Res, 40(10), 4288-4297.

Examples

```
data(dexus)
result <- dexus(countsMontgomery[1:10, ])
```

dexus.parallel *A parallel version of DEXUS.*

Description

Speeds up DEXUS by using multiple processors. Uses the parallel package to parallelize a DEXUS call.

Usage

```
dexus.parallel(X, ncores = 2, normalization = "RLE",
               ignoreIfAllCountsSmaller = 1, resultObject = "S4", ...)
```

Arguments

X	Either a vector of counts or a raw data matrix, where columns are interpreted as samples and rows as genomic regions.
ncores	The number of cores (CPUs) that will be used by the parallelization.
normalization	Normalization method to be used. (Default="RLE")
ignoreIfAllCountsSmaller	A transcript is considered as not expressed if all counts are smaller than the given value. (Default=1)
resultObject	Type of the result object; can either be a list ("list") or an instance of "DEXUS-Result" ("S4"). (Default="S4").
...	Other options to be passed to dexus().

Value

"list"

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
result <- dexus.parallel(countsPickrell[1:10, ], ncores=1)
```

DEXUSResult-class Class "DEXUSResult"

Description

This class contains the result of an RNA-Seq data analysis. The class contains the transcript names together with the parameters per condition, i.e., overdispersion and mean. Further it contains informative/non-informative values or p-values.

Objects from the Class

Objects can be created by calls of the form `new("DEXUSResult", ...)`.

Slots

`transcriptNames` The names of the transcripts, genes, exons, or regions of interest

`sampleNames` The sample names as they were given in the input matrix.

`inputData` The original read count matrix.

`normalizedData` The normalized read count matrix.

`sizeFactors` The size factors that were calculated for the normalization. This is that factor that scales each column or sample.

`INIValues` An informative/non-informative value for each sample that measures the evidence for differential expression.

`INIThreshold` The threshold for the I/NI values. Transcript with I/NI values above the threshold will be considered as differentially expressed.

`INICalls` A binary value for each transcript indicating whether it is differentially expressed.

`pvals` In case of two known conditions or multiple known conditions it is possible to calculate a *p*-value for each transcript. This value is given in this slot.

`responsibilities` A matrix of the size of the input matrix. It indicates the condition for each sample and transcript. The condition named "1" is the major condition. All other conditions are minor conditions. In case of supervised (two known conditions or multiple known conditions) analyses this clustering matrix will be the same for all transcripts.

`posteriorProbs` An array of the dimension of transcripts times samples times conditions. It gives the probability that a certain read count *x* was generated under a condition.

`logFC` The log foldchanges between the conditions. The reference is always condition "1".

`conditionSizes` The ratio of samples belonging to that condition. These are the α_i values of the model.

`sizeParameters` The size parameter estimates for each condition. These are the r_i values of the model.

`means` The mean of each condition. The μ_i values of the model.

`dispersions` The dispersion estimates for each condition. The inverse size parameters.

`params` The input parameters of the DEXUS algorithm.

Methods

[Subsetting of a DEXUSResult.

as.data.frame Converts the result object into a data frame.

conditionSizes Returns the condition sizes or α_i parameters of the model.

dispersions Returns the dispersion, i.e. the inverse size parameters, of the model.

INI I/NI filtering of the result object.

INICalls Returns a logical value indication whether this transcript is differentially expressed or not.

INIThreshold Returns the thresholds for the I/NI values.

INIThreshold<- Sets the I/NI threshold. I/NI calls will be changed accordingly.

INIValues Returns the I/NI values.

inputData Returns the input read counts.

logFC Returns the log foldchange with respect to the first condition.

means Returns the mean per condition.

normalizedData Returns the normalized data.

params Returns a list of input parameters of DEXUS.

plot Plots a heatmap of the read counts of the top genes.

posteriorProbs Returns an array of posterior probabilities.

pvals Returns the p -values per transcript in supervised mode.

responsibilities Returns the clustering vector.

sampleNames Returns the sample names.

show Displays a data frame of results.

sizeFactors Returns the size factors used for normalization.

sizeParameters Returns the size parameters, i.e. the r_i values of the model.

sort Sorts the result object by I/NI values or p -values.

transcriptNames Returns the transcript names.

Author(s)

Guenter Klambauer

Examples

```
showClass("DEXUSResult")
```

getSizeNB	<i>Maximum-likelihood and maximum-a-posteriori estimators for the negative binomial distribution.</i>
-----------	---

Description

Estimates the size parameter of a a negative binomial distribution from given data.

Usage

```
getSizeNB(x, maxCyc = 1000, eta = 0, rmax = Inf,
method = "bisection")
```

Arguments

x	The input data. Must be a numeric vector.
maxCyc	The maximum number of cycles of the numeric procedure to find the estimator. (Default = 1000).
eta	The weight of the exponential prior. The higher eta, the lower the estimate for the size parameter. Setting eta = 0 means that the prior is not used and, therefore, the maximum-likelihood estimator is calculated. (Default = 0).
rmax	Upper bound on the size parameter. This corresponds to a truncated exponential prior. If not used there is a non-zero probability that the estimator for the size parameter is ∞ . (Default = Inf).
method	The procedure used to solve the equation

$$\sum_{k=1}^N \psi(x_i + r) - N\psi(r) + N \log \left(\frac{r}{r + 1/N \sum_{i=1}^N x_i} \right) - \eta = 0$$

for r .

This can either be "bisection" or "regula falsi". (Default="bisection").

Details

Depending on the parameters you can either obtain the *Maximum-likelihood estimator* or the *maximum-a-posteriori estimator* using an exponential prior.

maximum-likelihood estimator	eta = 0
maximum-a-posteriori estimator	eta > 0

By setting the variable rmax to a positive value one can enforce an upper bound on the parameter.

The inverse of the size parameter is the overdispersion parameter.

Value

"numeric" An estimate of the size parameter of the negative binomial distribution. The overdispersion parameter is the inverse of the size parameter of a negative binomial distribution

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
x <- rnbinom(mu=50, size=5, n=10)
getSizeNB(x)
```

INI

INI filtering of a DEXUS result.

Description

This function filters the result object for informative transcripts. Transcripts with an I/NI value below the given threshold are filtered out.

Arguments

object	An instance of "DEXUSResult".
threshold	A numeric determining the threshold for the I/NI values.

Value

An instance of "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
res <- dexus(countsBottomly[1:100, ])
INI(res)
```

INIThreshold<- *Set the I/NI threshold.*

Description

This generic function sets the threshold of the I/NI value. Transcripts with I/NI values above the I/NI threshold are considered as differentially expressed. The results of DEXUS are stored as an instance of `DEXUSResult-class`.

Arguments

object An instance of "DEXUSResult".
value A numeric to be used for thresholding the I/NI values.

Value

INIThreshold<- returns an instance of "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
result <- dexus(countsBottomly[1:20,1:10])
INIThreshold(result) <- 0.1
```

normalizeData *Normalization of RNA-Seq count data.*

Description

Normalizes RNA-seq count data using previously published approaches. Each samples' read counts are corrected by a normalizing factor. The options are "RLE" by (Anders and Huber, 2010), and "upperquartile" by (Bullard et al., 2010).

Usage

```
normalizeData(X, normalization)
```

Arguments

X data a raw data matrix, where' columns are interpreted as samples and rows as genomic regions.
normalization method used for normalizing the reads. RLE is the method used by (Anders and Huber, 2010), upperquartile is the Upper-Quartile method from (Bullard et al., 2010), and none deactivates normalization. (Default = "RLE").

Value

"list" A list containing the normalized data (in its "X" component) as well as the size-factors used for the normalization ("sizeFactors").

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

References

Anders, S. and Huber, W. (2010). *Differential expression analysis for sequence count data*. *Genome Biol*, 11(10), R106.

Bullard, J. H., Purdom, E., Hansen, K. D., and Dudoit, S. (2010). *Evaluation of statistical methods for normalization and differential expression in mRNA-seq experiments*. *BMC Bioinformatics*, 11, 94.

Examples

```
data(dexus)
norm <- normalizeData(countsBottomly, "RLE")
```

 plot

Visualization of a result of the DEXUS algorithm.

Description

Plots a heatmap of the log read counts of the top ranked genes or of selected genes.

Arguments

x	An instance of "CNVDetectionResult"
idx	The indices or the transcript names of the transcripts that should be visualized as heatmap.
cexSamples	Size of the column labels, i.e. the samples.
cexGenes	Size of the row labels, i.e. the transcripts.
newColNames	renames the samples.
type	Mark the samples, that do not belong to the major class by crosses ("crosses"), or boxes ("boxes").

Value

Generates a heatmap of the expression values of the top-ranked transcripts.

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
r <- dexus(countsBottomly[1:100, ])
plot(r)
```

sort	<i>Sorting a DEXUS result.</i>
------	--------------------------------

Description

This function sorts the result object by I/NI values or p-values such that the transcripts with the highest I/NI value or the lowest p-value are ranked first.

Arguments

object An instance of "DEXUSResult".

Value

An instance of "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
res <- dexus(countsBottomly[1:100, ])
sort(res)
```

‘[‘	<i>Subsetting a "DEXUSResult".</i>
-----	------------------------------------

Description

Information about specific transcripts can be accessed in the "DEXUSResult" object by using the standard brackets "[idx]" for subsetting. Either transcript names or transcript indices can be used.

Arguments

x "DEXUSResult"
i Either a numeric vector of indices or a character vector containing the transcript names.

Value

An instance of "DEXUSResult".

Author(s)

Guenter Klambauer <klambauer@bioinf.jku.at> and Thomas Unterthiner <unterthiner@bioinf.jku.at>

Examples

```
data(dexus)
res <- dexus(countsBottomly[1:100, ])
res["ENSMUSG00000000486"]
res[50:55]
```

Index

- *Topic **classes**
 - DEXUSResult-class, 12
- *Topic **datasets**
 - countsBottomly, 3
 - countsGilad, 4
 - countsLi, 4
 - countsMontgomery, 5
 - countsPickrell, 5
- [,DEXUSResult,character-method (DEXUSResult-class), 12
- [,DEXUSResult,logical-method (DEXUSResult-class), 12
- [,DEXUSResult,numeric-method (DEXUSResult-class), 12
- ['[, 18
- ['[,DEXUSResult,character-method ('['[, 18
- ['[,DEXUSResult,logical-method ('['[, 18
- ['[,DEXUSResult,numeric-method ('['[, 18
- accessors, 2
- as.data.frame,DEXUSResult-method (DEXUSResult-class), 12
- conditionSizes (accessors), 2
- conditionSizes,DEXUSResult-method (DEXUSResult-class), 12
- countsBottomly, 3
- countsGilad, 4
- countsLi, 4
- countsMontgomery, 5
- countsPickrell, 5
- dexss, 6
- DEXSS, (dexss), 6
- dexus, 8
- DEXUS, (dexus), 8
- dexus.parallel, 11
- DEXUSResult-class, 12
- dispersions (accessors), 2
- dispersions,DEXUSResult-method (DEXUSResult-class), 12
- getSizeNB, 14
- INI, 15
- INI,DEXUSResult-method (DEXUSResult-class), 12
- INICalls (accessors), 2
- INICalls,DEXUSResult-method (DEXUSResult-class), 12
- INIThreshold (accessors), 2
- INIThreshold,DEXUSResult-method (DEXUSResult-class), 12
- INIThreshold-set (INIThreshold<-), 16
- INIThreshold<-, 16
- INIThreshold<-,DEXUSResult-method (DEXUSResult-class), 12
- INIValues (accessors), 2
- INIValues,DEXUSResult-method (DEXUSResult-class), 12
- inputData (accessors), 2
- inputData,DEXUSResult-method (DEXUSResult-class), 12
- logFC (accessors), 2
- logFC,DEXUSResult-method (DEXUSResult-class), 12
- means (accessors), 2
- means,DEXUSResult-method (DEXUSResult-class), 12
- normalizeData, 16
- normalizedData (accessors), 2
- normalizedData,DEXUSResult-method (DEXUSResult-class), 12
- params (accessors), 2
- params,DEXUSResult-method (DEXUSResult-class), 12

plot, [17](#)
plot, DEXUSResult, missing-method
 (DEXUSResult-class), [12](#)
posteriorProbs (accessors), [2](#)
posteriorProbs, DEXUSResult-method
 (DEXUSResult-class), [12](#)
pvals (accessors), [2](#)
pvals, DEXUSResult-method
 (DEXUSResult-class), [12](#)

responsibilities (accessors), [2](#)
responsibilities, DEXUSResult-method
 (DEXUSResult-class), [12](#)

sampleNames (accessors), [2](#)
sampleNames, DEXUSResult-method
 (DEXUSResult-class), [12](#)
show, DEXUSResult-method
 (DEXUSResult-class), [12](#)
sizeFactors (accessors), [2](#)
sizeFactors, DEXUSResult-method
 (DEXUSResult-class), [12](#)
sizeParameters (accessors), [2](#)
sizeParameters, DEXUSResult-method
 (DEXUSResult-class), [12](#)
sort, [18](#)
sort, DEXUSResult-method
 (DEXUSResult-class), [12](#)

transcriptNames (accessors), [2](#)
transcriptNames, DEXUSResult-method
 (DEXUSResult-class), [12](#)